



## Company Confidential

**Report Number:** 34-95-111

**Title:** Enforcement Residue Analytical Method for GOAL Herbicide (oxyfluorfen) in Crop Commodities with GC/MS Confirmation

**Author(s):** Martin, D. A., Zhang, Q.

**Date:** 04/29/1996

**Common Name:** oxyfluorfen

**Product Name(s):** Goal, Technical

**Performing Lab:** Centre Analytical La, Spring House Labs., McKenzie

**EPA Guideline No:** 860.1340

**Old Guideline:** 171-4D

**EPA Guideline Code:** 171-4(c)

**Description:** Residue Analytical Methods, Plant

**MRID Number:** 44400202

**Accession Number:** -

**Submitted Date:** 10/14/1997

**TIMS Number:** 7434

**Temporary Location:** 078

**STUDY TITLE**

Enforcement Residue Analytical Method  
for GOAL® Herbicide (oxyfluorfen) in Crop Commodities  
with GC/MS Confirmation

**DATA REQUIREMENT**

Guideline 171-4(c)

**STUDY DIRECTOR**

D. A. Martin

**AUTHORS**

D. A. Martin, Q. Zhang

**STUDY COMPLETED ON**

4/29/96

**SPONSOR**

Rohm and Haas Company  
727 Norristown Road, Spring House, PA 19477-0904

**PERFORMING LABORATORIES**

Rohm and Haas Company  
727 Norristown Road, Spring House, PA 19477-0904

Centre Analytical Laboratories, Inc.  
3048 Research Drive, State College, PA 16801

McKenzie Laboratories  
3725 E. Atlanta Ave., Suite One, Phoenix, AZ 85040

**STUDY PROTOCOL**

34P-95-92

**LABORATORY PROJECT ID**

Rohm and Haas Technical Report No. 34-95-111

This page intentionally left blank for statement of data confidentiality claims.

### GLP Compliance Statement

The information contained in this report pertains to test method development and, as such, is not required to adhere to 40 CFR 160. The data, however, were compiled from the following:

- GLP compliant residue studies submitted to the EPA,
- additional fortification results for 3 of the above studies received after the submission due date, and,
- a new study conducted to validate the Oxyfluorfen (Goal®) Crop Residue Method, TR 34-94-150 using GC/MS as the confirmatory detection.

This report and work should be considered in compliance with Good Laboratory Practice standards as defined by EPA (40 CFR 160) with the following exceptions:

1. This report is a compilation of multiple studies by several study directors rather than a final report for a single study.
2. Additional fortification data for figs, olives, and pomegranates were obtained after submission of the respective final reports and were authorized by letter by the study director rather than via protocol amendment.

This report is a true and accurate representation of the work performed as attested to in previously submitted reports and by the undersigned.

Dennis A. Martin 4/29/96  
Dennis A. Martin Date

Study Director, Author  
Rohm and Haas Company

Stanley S. Stavinski April 29/96  
Stanley S. Stavinski Date

Sponsor  
Rohm and Haas Company

Qipan Zhang 4/29/96  
Qipan Zhang Date

Author  
Rohm and Haas Company

Stanley S. Stavinski April 29/96  
Stanley S. Stavinski Date

Submitter/Applicant  
Rohm and Haas Company

### QUALITY ASSURANCE STATEMENT

This report and the original raw data have been reviewed by the Quality Assurance Unit of the Rohm and Haas Company Agricultural Products - Development Research Department and has been found to be a true and accurate reflection of the study performed.

Activity	Date Conducted	Date Reported <sup>1</sup>	OAU <sup>2</sup>
<b>Additional work for submitted reports</b>			
<b>Figs</b>			
Protocol Review	7/10/95	8/18/95 (7/20/95)	CAL
Extraction	7/11/95	8/18/95 (8/16/95)	CAL
Raw Data Review	8/7-8/95	8/18/95 (8/16/95)	CAL
Raw Data Review	8/18/95	12/18/95 (9/25/95)	CAL
<b>Olives</b>			
Protocol Review	7/7/95	8/18/95 (7/20/95)	CAL
Extraction	7/7/95	8/18/95 (8/16/95)	CAL
Raw Data Review	8/8/95	8/18/95 (8/16/95)	CAL
Raw Data Review	8/18/95	12/18/95 (9/25/95)	CAL
<b>Pomegranates</b>			
Protocol Review	7/10/95	8/18/95 (7/20/95)	CAL
Extraction	7/7/95	8/18/95 (8/16/95)	CAL
Raw Data Review	8/8/95	8/18/95 (8/16/95)	CAL
Raw Data Review	8/18/95	12/18/95 (9/25/95)	CAL
<b>Peach report to be submitted</b>			
<b>PA Trial</b>			
Inspection, harvest	8/17/94	8/18/94	R&H
Field report audit	7/21/95	7/25/95	R&H
Analytical subprotocol review	6/7/95	6/7/95	R&H
Inspections, analysis	7/6, 7/95	7/17/95 (7/12/95)	McK
Data audit, analysis	8/4/95	8/4/95 (8/4/95)	McK
<b>GA Trial</b>			
Field report audit	1/30/96	2/1/96	R&H
Inspections, analysis	11/21, 22, 28/95	12/11/95 (12/7/95)	McK
Inspections, logs	11/29/95	12/11/95 (12/7/95)	McK
Final report audit	4/1/96	4/1/96	R&H
<b>Current study (34P-95-92)</b>			
protocol review	10/17/95	10/17/95	R&H
protocol review	11/01/95	12/18/95 (11/30/95)	CAL
extraction	11/02/95	12/18/95 (11/30/95)	CAL
extraction	11/09/95	12/18/95 (11/30/95)	CAL
raw data review	12/06/95	2/02/96 (1/10/96)	CAL
<b>Final report audits</b>			
	3/20/96	3/20/96	R&H
	4/4/96	4/4/96	R&H
	4/26/96	4/26/96	R&H

1 Dates in parentheses pertain to notification of contract laboratory management.

2 CAL= Centre Analytical Laboratories  
R&H= Rohm and Haas Company  
McK= McKenzie Laboratories

*H. E. Aller*      *April 29, 1996*  
\_\_\_\_\_  
H. E. Aller      Date  
QA Officer

Rohm and Haas Company

**TABLE OF CONTENTS**

	<b>Page</b>
1. Objective	9
2. Summary	9
3. Introduction	9
4. Experimental Compound	9
5. Chemicals, Supplies, and Solutions	9
6. Equipment	11
7. Method	12
7.1 Flow Diagram	12
7.2 Sample Processing	13
7.3 Sample Extraction	13
7.4 Liquid/Liquid Partition	14
7.5 Column Chromatography	14
7.6 Basic Alumina Solid Phase Extraction (SPE)	15
8. Primary Detection Method: GC/ECD Quantitation	16
8.1 Preparation of Standard Curves	16
8.2 Instrumentation	16
8.3 Quantitation	18
8.4 Fortification	18
9. Confirmatory Detection Method: GC/MS Quantitation	19
10. Results	19
11. Discussion	21
12. Project Information	22

**Tables**

1.	Apple Oxyfluorfen Fortification Summary	25
2.	Artichoke Oxyfluorfen Fortification Summary	26
3.	Avocado Oxyfluorfen Fortification Summary	27
4.	Cherry Oxyfluorfen Fortification Summary	28
5.	Fig Oxyfluorfen Fortification Summary	29
6.	Grape Oxyfluorfen Fortification Summary	30
7.	Kiwi Oxyfluorfen Fortification Summary	32
8.	Olive Oxyfluorfen Fortification Summary	33
9.	Peach Oxyfluorfen Fortification Summary	34
10.	Peanut Nutmeat Oxyfluorfen Fortification Summary	35
11.	Pomegranate Oxyfluorfen Fortification Summary	36
12.	Soybean Oxyfluorfen Fortification Summary	37
13.	Peanut Hay Oxyfluorfen Fortification Summary	38
14.	Peanut Shell Oxyfluorfen Fortification Summary	39
15.	Peanut Vine Oxyfluorfen Fortification Summary	40
16.	Grape Oxyfluorfen Fortification Summary Rtx-50 Confirmatory Column	41
17.	Soybean Oxyfluorfen Fortification Summary Rtx-50 Confirmatory Column	42
18.	Oxyfluorfen GC/MS Confirmatory Data Summary	43
19.	Oxyfluorfen Data Summary by Crop	44

**Figures**

1.	Apple Control 95-1408-003	45
2.	Apple Fortification 0.01 PPM 95-1408-003	46
3.	Apple Fortification 0.10 PPM 95-1408-003	47
4.	Apple Fortification 0.50 PPM 95-1408-003	48
5.	Artichoke Control 94-0060-001	49
6.	Artichoke Fortification 0.01 PPM 94-0060-001	50
7.	Artichoke Fortification 0.10 PPM 94-0060-001	51
8.	Artichoke Fortification 0.50 PPM 94-0060-001	52
9.	Avocado Control 94-0141-001	53
10.	Avocado Fortification 0.01 PPM 94-0141-001	54
11.	Avocado Fortification 0.10 PPM 94-0141-001	55
12.	Avocado Fortification 0.50 PPM 94-0141-001	56
13.	Cherry Control 92-0041-002	57
14.	Cherry Fortification 0.01 PPM 92-0041-002	58
15.	Cherry Fortification 0.10 PPM 92-0041-002	59
16.	Cherry Fortification 0.50 PPM 92-0041-002	60
17.	Fig Control 94-0142-001	61
18.	Fig Fortification 0.01 PPM 94-0142-001	62

19.	Fig Fortification 0.10 PPM 94-0142-001	63
20.	Fig Fortification 0.50 PPM 94-0142-001	64
21.	Grape Control 93-0012-001	65
22.	Grape Fortification 0.01 PPM 93-0012-001	66
23.	Grape Fortification 0.10 PPM 93-0012-001	67
24.	Grape Fortification 0.50 PPM 93-0012-001	68
25.	Kiwi Control 94-0146-001	69
26.	Kiwi Fortification 0.01 PPM 94-0146-001	70
27.	Kiwi Fortification 0.10 PPM 94-0146-001	71
28.	Kiwi Fortification 0.50 PPM 94-0146-001	72
29.	Olive Control 94-0172-001	73
30.	Olive Fortification 0.01 PPM 94-0172-001	74
31.	Olive Fortification 0.10 PPM 94-0172-001	75
32.	Olive Fortification 0.50 PPM 94-0172-001	76
33.	Peach Control 94-0117-001	77
34.	Peach Fortification 0.01 PPM 94-0117-001	78
35.	Peach Fortification 0.10 PPM 94-0117-001	79
36.	Peach Fortification 0.50 PPM 94-0117-001	80
37.	Peanut Nutmeat Control 94-0176-010	81
38.	Peanut Nutmeat Fortification 0.01 PPM 94-0176-010	82
39.	Peanut Nutmeat Fortification 0.10 PPM 94-0176-010	83
40.	Peanut Nutmeat Fortification 0.50 PPM 94-0176-010	84
41.	Pomegranate Control 94-0143-001	85
42.	Pomegranate Fortification 0.01 PPM 94-0143-001	86
43.	Pomegranate Fortification 0.10 PPM 94-0143-001	87
44.	Pomegranate Fortification 0.50 PPM 94-0143-001	88
45.	Soybean Control 92-0107-014	89
46.	Soybean Fortification 0.01 PPM 92-0107-014	90
47.	Soybean Fortification 0.10 PPM 92-0107-014	91
48.	Soybean Fortification 0.50 PPM 92-0107-014	92
49.	Peanut Hay Control 93-0142-003	93
50.	Peanut Hay Fortification 0.01 PPM 93-0142-003	94
51.	Peanut Hay Fortification 0.10 PPM 93-0142-003	95
52.	Peanut Hay Fortification 0.50 PPM 93-0142-003	96
53.	Peanut Shell Control 93-0157-004	97
54.	Peanut Shell Fortification 0.01 PPM 93-0157-004	98
55.	Peanut Shell Fortification 0.10 PPM 93-0157-004	99
56.	Peanut Shell Fortification 0.50 PPM 93-0157-004	100
57.	Peanut Vine Control 93-0142-001	101
58.	Peanut Vine Fortification 0.01 PPM 93-0142-001	102
59.	Peanut Vine Fortification 0.10 PPM 93-0142-001	103
60.	Peanut Vine Fortification 0.50 PPM 93-0142-001	104
61.	Grape Control (Rtx-50 confirmatory)	105
62.	Grape Fortification 0.01 PPM (Rtx-50 confirmatory)	106
63.	Grape Fortification 0.10 PPM (Rtx-50 confirmatory)	107



64.	Grape Fortification 0.50 PPM (Rtx-50 confirmatory)	108
65.	Soybean Control (Rtx-50 confirmatory)	109
66.	Soybean Fortification 0.01 PPM (Rtx-50 confirmatory)	110
67.	Soybean Fortification 0.10 PPM (Rtx-50 confirmatory)	111
68.	Soybean Fortification 0.50 PPM (Rtx-50 confirmatory)	112
69.	Standard 0.005 µg/ml (GC/ECD)	113
70.	Standard 0.01 µg/ml (GC/ECD)	114
71.	Standard 0.02 µg/ml (GC/ECD)	115
72.	Standard 0.10 µg/ml (GC/ECD)	116
73.	Standard 0.15 µg/ml (GC/ECD)	117
74.	Oxyfluorfen Calibration Curve (GC/ECD)	118
75.	Grape Control (GC/MS)	119
76.	Grape Fortification 0.025 ppm (GC/MS)	120
77.	Grape Fortification 0.05 ppm (GC/MS)	121
78.	Grape Fortification 0.125 ppm (GC/MS)	122
79.	Peach Control (GC/MS)	123
80.	Peach Fortification 0.025 ppm (GC/MS)	124
81.	Peach Fortification 0.05 ppm (GC/MS)	125
82.	Peach Fortification 0.125 ppm (GC/MS)	126
83.	Peanut Nutmeat Control (GC/MS)	127
84.	Peanut Nutmeat Fortification 0.025 ppm (GC/MS)	128
85.	Peanut Nutmeat Fortification 0.05 ppm (GC/MS)	129
86.	Peanut Nutmeat Fortification 0.125 ppm (GC/MS)	130
87.	Peanut Shell Control (GC/MS)	131
88.	Peanut Shell Fortification 0.025 ppm (GC/MS)	132
89.	Peanut Shell Fortification 0.05 ppm (GC/MS)	133
90.	Peanut Shell Fortification 0.125 ppm (GC/MS)	134
91.	Standard 0.025 µg/ml (GC/MS)	135
92.	Standard 0.05 µg/ml (GC/MS)	136
93.	Standard 0.10 µg/ml (GC/MS)	137
94.	Standard 0.50 µg/ml (GC/MS)	138
95.	Oxyfluorfen Calibration Curve (GC/MS)	139

## Appendices

I.	GC/MS Method Validation Protocol	141
	Certificate of Analysis (RH-32915)	152
	GLP Statements/Quality Assurance Statements	153
II.	Representative Analytical Run (GC/ECD)	168
	Representative Analytical Run (GC/MS)	191

## 1. Objective

The original analytical method for crops, TR 34-94-150, was issued December 21, 1994 and later submitted to the EPA as Appendix 8 in MRID 42793300. This enforcement method is a reissue of the original with minor modifications, including the addition of GC/MS confirmatory detection.

## 2. Summary


This report outlines the oxyfluorfen (parent only) residue analytical method for crop commodities. Oxyfluorfen is extracted from crops with acetonitrile. The residues are then further purified by liquid-liquid partition, silica gel or Florisil column chromatography and basic alumina solid phase extraction. Quantitation is performed by gas-liquid chromatography on a Restex Rtx-200 or Rtx-50 column with electron capture detection. Confirmatory detection is accomplished through GC/MS.

## 3. Introduction

Oxyfluorfen (GOAL® Herbicide) is a registered herbicide used on a wide variety of crops. In order to maintain the commercial registration, an analytical residue method for crops is required to generate residue data for the setting and enforcement of tolerances.

## 4. Experimental Compound

This analytical method is developed for oxyfluorfen (GOAL Herbicide), RH-32915, which has a CAS registry number 42874-03-3. The Chemical Abstracts name is 2-chloro-1-(3-ethoxy-4-nitrophenoxy)-4-(trifluoromethyl) benzene. Its structure is shown below.



## 5. Chemicals, Supplies, Solutions

<u>Chemicals and Supplies</u>	<u>Supplier *</u>
Acetonitrile, HPLC Reagent	J. T. Baker
Cotton, sterile and absorbent	Johnson & Johnson
Ethyl Ether, HPLC Reagent	J. T. Baker

Ethyl Acetate, HPLC Reagent	J. T. Baker
Florisil, 60/100 mesh	J. T. Baker
Hexane, HPLC Reagent	J. T. Baker
Methanol, HPLC Reagent	J. T. Baker
Petroleum Ether, Omnisolv	EM Science**
RH-32915 Analytical Standard	Rohm and Haas
Silica Gel, 40-140 mesh	J. T. Baker Catalog # 3404-01
Sodium Chloride	J. T. Baker
Sodium Sulfate, anhydrous granular	Fisher
Supelclean LC-Alumina-B SPE tubes	Supelco Catalog # 5-7085
Toluene, HPLC Reagent	J. T. Baker
Water, Milli-Q	Millipore System

### Solutions

Saturated sodium chloride aqueous solution

5/95 (v/v) ethyl ether/petroleum ether

30/70 (v/v) ethyl ether/petroleum ether

50/50 (v/v) ethyl ether/petroleum ether

40/60 (v/v) ethyl ether/petroleum ether

50/50 (v/v) ethyl acetate/hexane

5/95 (v/v) ethyl acetate/hexane

\* Equivalent chemicals, equipment and instrumentation can be substituted, when shown to be appropriate.

\*\* EM Science brand is recommended, for details see Section 10.

## 6. Equipment

<u>Item</u>	<u>Supplier*</u>
Buchner Funnel, porcelain	Coors
Chromatography Columns: Glass 250 mm length, 16 mm I.D.	Ace Glass
Centrifuge Bottles, 250 ml	Fisher
Centrifuge	Beckman
Filter Paper, Whatman #41	Whatman
G10 Gyrotory Shaker	New Brunswick Scientific
Pipets, glass disposable, 10 ml, 5 ml, 1 ml	Fisher
Rotary Evaporator, with dry ice trap	Brinkman
Round Bottom Flask (100, 250, 500 ml, 24/40 ST)	Kimax, Pyrex
Test tubes, 40 ml 12.0 mm I.D.	Kontes
Separatory Funnels, 500 ml	Kimax, Pyrex
Standard laboratory equipment, balance, beakers, etc.	
Visiprep Solid Phase Extraction Vacuum Manifold - 60 ml reservoirs - adapter for sample reservoirs	Supelco Catalog # 5-7044
Volumetric Flasks (100ml)	Pyrex

\* Equivalent chemicals, equipment and instrumentation can be substituted, when shown to be appropriate.

## 7. Method


The crop commodities analyzed with this method can be separated into two groups based on the path of analysis. They are designated Group I and Group II for ease in detailing the differences in the methodology, specifically during the Column Chromatography. The groups are:

Group I: apple, artichoke, avocado, cherry, fig, grape, kiwi, olive, peach, peanut nutmeat, pomegranate, soybean

Group II: peanut hay, shell, vine

### 7.1 Flow Diagram

Below is a flow diagram giving a general description of the method. The individual steps are detailed following the diagram.



## 7.2 Sample Processing

All crop samples are finely and homogenously chopped in preparation for sample extraction. *Apple, artichoke, fig, grape, kiwi, peanut hay, peanut nutmeat, pomegranate, and peanut vine* samples are chopped with dry ice in a Hobart food chopper. For *avocado, cherry, peach, and olive* samples, the pits are removed prior to processing with dry ice in a Hobart food chopper. Any suitable commercial food chopper can be substituted if a Hobart is not available. All samples with dry ice are left open and allowed to sublime overnight in the freezer. *Peanut shell* and *soybean* samples are processed dry in a Wiley Mill. All samples are maintained frozen prior to analysis.

## 7.3 Sample Extraction

*All the glassware and centrifuge bottles must be pre-washed with methanol twice before any experiment.*

For all crops, weigh a homogenized 5.0 g sample of the crop into a 250 ml centrifuge bottle; and, fortify as needed for fresh spikes. Add 100 ml of acetonitrile and then homogenize with a Gyrotory Shaker at 600-800 rpm for 20 min. Centrifuge each extraction mixture at 5000 rpm for 20 minutes at 4°C\*. Decant the supernatant into a 500 ml separatory funnel. Add another 100 ml of acetonitrile to the centrifuge bottle containing the crop pellet, then homogenize with the Gyrotory Shaker at 600-800 rpm for another 20 minutes and centrifuge again as described above. Decant the supernatant into the 500 ml separatory funnel with the first acetonitrile extract. Add 200 ml petroleum ether to the 500 ml separatory funnel containing the 200 ml acetonitrile extract. Shake for 20 seconds and vent into a fume hood. Allow layers to separate and drain the acetonitrile (lower) phase into a 500 ml 24/40 ST round bottom flask. Discard the petroleum ether layer to waste.

Evaporate the acetonitrile gently to approximately 50 ml by rotary evaporator at 60°C under diminished pressure. Transfer the acetonitrile extract to a new 500 ml separatory funnel for the next step, liquid/liquid partition.

\* Option: Filter the samples under suction with a vacuum pump (about 50 mm Hg) through a Buchner funnel with Whatman #41 filter paper which has been prewashed with acetonitrile and dried.

#### 7.4 Liquid/Liquid Partition

For all crops, add 100 ml petroleum ether to the 500 ml round bottom flask which previously contained the acetonitrile extract. Transfer the petroleum ether to the 500 ml separatory funnel containing the acetonitrile extract. Add 10 ml of saturated NaCl solution and 200 ml water to the separatory funnel. Shake vigorously for 2 minutes and vent into a fume hood. Allow the layers to separate and drain the acetonitrile/water (lower) phase and reserve for further partitioning. Then collect the petroleum ether phase in a new 500 ml 24/40 ST round bottom flask. Add another 100 ml petroleum ether to the 500 ml flask that contained the acetonitrile extract, transfer to the separatory funnel and partition the reserved acetonitrile/water phase again, as described above.

Rinse the separatory funnels with methanol, discarding to waste, and allow the funnels to dry. Transfer the petroleum ether fractions to the separatory funnel and then partition with 200 ml water. Shake for one minute, allow layers to separate and drain the water (lower) phase to waste. Collect the petroleum ether layer into a 500 ml 24/40 ST round bottom flask. Evaporate gently to approximately 5 ml of petroleum ether by rotary evaporator at 40°C under diminished pressure. The sample extract is then ready for the next step, column chromatography.

#### 7.5 Column Chromatography

*Group I follows the Silica Gel Open Column step. Group II follows the Florisil Open Column step.*

Silica Gel Open Column (Group I: apple, artichoke, avocado, cherry, fig, grape, kiwi, olive, peach, peanut nutmeat, pomegranate, soybean)

Activate the silica gel, mesh 40-140, by heating for 24 hours at 200°C. Remove from the oven and store in tightly capped jars in a desiccator. Pack a 250 mm X 16 mm I.D. glass column plugged with cotton with 20 ml of the activated silica gel. Top the column with 10-15 g (approx. 2-3 cm column high) of anhydrous granular sodium sulfate. Precondition the column with 30 ml of petroleum ether.

Add the 5 ml of the remaining petroleum ether extract of the sample, section 6.4, to the column and elute to the top of the silica gel bed. Add 50 ml of petroleum ether to the 500 ml round bottom flask, transfer to the column, and elute to the top of the silica gel bed. Add 20 ml of 95/5 (v/v) petroleum ether/ethyl ether to the 500 ml round bottom flask, transfer to the column, and elute to the top of the silica gel bed. Discard all washes to this point.

Elute the oxyfluorfen by adding 75 ml of 60/40 (v/v) petroleum ether/ethyl ether to the 500 ml round bottom flask, transferring to the column, and eluting to the top of the silica gel bed. Collect the eluate in a 250 ml 24/40 ST round bottom flask. Evaporate to dryness at 45 °C under diminished pressure. Redissolve the residue in 5 ml of hexane. The sample is now ready for the next step, solid phase extraction.

#### Florisil Open Column (Group II: peanut hay, shell, vine)

Activate the Florisil, mesh 60-100, by heating for 24 hours at 150°C. Remove from the oven and store in tightly capped jars in a desiccator. Pack a 250 mm X 16 mm I.D. glass column plugged with cotton with 20 ml of the activated Florisil. Top the column with 10-15 g (approx. 2-3 cm column high) of anhydrous sodium sulfate. Precondition the column with 30 ml of petroleum ether.

Add the 5 ml of the remaining petroleum ether extract of the sample, section 6.4, to the column and elute to the top of the Florisil bed. Add 50 ml of petroleum ether to the 500 ml round bottom flask, transfer to the column, and elute to the top of the Florisil bed. Add 20 ml of 70/30 (v/v) petroleum ether/ethyl ether to the 500 ml round bottom flask, transfer to the column, and elute to the top of the Florisil bed. Discard all washes to this point.

Elute the oxyfluorfen residue by adding 50 ml of 50/50 (v/v) petroleum ether/ethyl ether to the 500 ml round bottom flask, transferring to the column, and eluting to the top of the Florisil bed. Collect the eluate in a 250 ml 24/40 ST round bottom flask. Evaporate to dryness at 45 °C under diminished pressure. Redissolve the residue in 5 ml of hexane. The sample is now ready for the next step, solid phase extraction.

#### 7.6 Basic Alumina Solid Phase Extraction (SPE)

Connect a 6 ml Basic Alumina SPE tube with a 60 ml reservoir using an adapter to a Supelco Visiprep Solid Phase Extraction Vacuum manifold or similar apparatus. Precondition the 6 ml Basic Alumina SPE tube with 10 ml ethyl acetate followed by 10 ml hexane. Flow rate through the column should be 2-3 mls per minute. For all crops, add the sample in 5 ml hexane to the column, allow it to elute until the sample is 2 mm above the top of the absorbent bed. Rinse and wash the sample flask 3-4 times with a total 40 ml of hexane, add to the column, and elute as described above. Discard all the hexane eluent at this point.





Gas Flows: Nitrogen (makeup) 40 ml/min  
Helium (column) 3.5 ml/min

Injection Volume: 1  $\mu$ l

Temperatures: Injector 265°C  
Detector 270°C  
Column 215°C  
Initial Hold 1.0 min

<u>Final Temp</u>	<u>Rate (°C/min)</u>	<u>Hold (min)</u>	<u>Total Time(min)</u>
250°C	10.0	17.0	21.5

Electron Capture Detection:	<u>Attenuation</u>	<u>Range</u>
	32	10

Under the conditions mentioned above, the typical retention time is 16.1 min. Conditions may need to be adjusted to achieve the optimum response of the instrument.

Column 2, Restex Rtx-50: 0.32 mm I.D., 60 m, 1.0  $\mu$  film  
Catalog #10557

Guard Column: Restex deactivated, uncoated fused silica gel, 5 m,  
0.53 mm I.D., Catalog # 10045

Column Connection: Supelco Glass Seal Connector. Catalog #2-0479

Gas Flows: Nitrogen (makeup) 40 ml/min  
Helium (column) 3.5 ml/min

Injection Volume: 1  $\mu$ l

Temperatures: Injector 265°C  
Detector 270°C  
  
Column 215°C  
Initial Hold 1.0 min

<u>Final Temp</u>	<u>Rate (°C/min)</u>	<u>Hold (min)</u>	<u>Total Time (min)</u>
250°C	5.0	25.0	33.0



## 9. Confirmatory Detection Method: GC/MS Quantitation

In order to confirm that residues found using the procedure described previously in this report are in fact oxyfluorfen, RH-32915, residues, a confirmatory detection method was established using GC/MS with selected ion monitoring. This method has a Limit of Quantitation of 0.025 ppm. Analytical results are shown in Table 18. Representative chromatograms are shown in Figures 76-96.

The detailed instrument information and conditions is given below:

Instrumentation: Hewlett-Packard model 5890 Series II GC; model 5972 mass selective detector

MS/Selected Ion Monitoring:

252 molecular weight ion

Dwell Time: 100 msec

Ionization Potential: 70 eV

Electron Multiplier Voltage: 1900-2318 V

Column: Rtx-50, 0.25 mm i.d. X 30 m, 0.5  $\mu$ m film thickness

Oven Temps.: Hold at 100°C for 2 min.  
100°-270°C at 10°C/min.  
Hold for 15 min.

Retention Time: 19.6 minutes

Injector Temp.: 280°C

Detector Temp.: 300°C

Carrier Gas/

Flow Rate: Helium/2.18 ml/min.

Head Pressure: 20 psi

Injection Mode: splitless

Injection Liner: silanized single taper

Injection Purge

Delay: 1.0 min.

Septum Purge: 2 ml/min.

Injection Volume: 2  $\mu$ l

## 10. Results

The results detailed in this method are covered in various reports and protocols. All reports used are listed in Section 12. Project Information. Many of these reports have been submitted and accepted by the EPA previously. Distinctions in individual results are made on the respective recovery tables, Tables 1-18.

For figs, olives, and pomegranates, there are a few fortifications reported here that were not reported in their respective technical reports, which were already submitted and accepted. Due to the timeline established for submission of those reports, the results were received too late to be included in them (refer to Tables 5, 8, and 11 for identification of data).

Six recoveries have been designated as outliers and are not included in the final calculation of both the average and standard deviation (see Tables 1, 7, 10, 14, 17). These points did not fall within two standard deviations when used in the calculation. It should also be noted that throughout all crop commodities the control samples did not have any detectable residues, that is above 0.003 ppm.

In GC/ECD, the limit of quantitation (LOQ) for all crop commodities except for peanut vine, shell, and hay is 0.01 ppm established in TR 34-94-150, Appendix 8 in MRID 42793300, and also by fortifications at this level, with a limit of detection of 0.003 ppm. For peanut vine, shell, and hay the limit of quantitation is 0.02 ppm, established by fortifications at this level, with a limit of detection of 0.007 ppm. Table 19 summarizes the LOQ and LOD values, along with average recoveries for each crop.

Tables 1-12 summarize the fortification data for the crop commodities detailed in this method as Group I. This data was obtained with GC/ECD detection using a Rtx-200 column. Average recoveries ranged from 85.4% in apple to 104.9% in soybean.

Tables 13-15 summarize the fortification data for the crop commodities detailed in this method as Group II. This data was also obtained with GC/ECD detection using a Rtx-200 column. Average recoveries were 92.2% in peanut hay, 95.3% in peanut shell, and 81.8% in peanut vine.

Tables 16 and 17 detail the results obtained from grape and soybean, Group I, analyses on the Rtx-50 confirmatory column. Recoveries averaged 92.6% in grapes and 104% in soybean.

In GC/MS, the limit of quantitation for all crop commodities is 0.025 ppm, established by fortifications at this level. The limit of detection is 0.008 ppm.

Table 18 displays the results of analyses of crops using GC/MS as the confirmatory detection method compared to analysis on GC/ECD. Individual recoveries ranged from 81% in peanut nutmeat to 117% in peanut shell.

Table 19 is a summary of the results for all crops detailed in this method, including average recovery, LOQ, and LOD values.

Representative chromatograms are shown in Figures 1-95, including standards, controls and three fortifications (alphabetical within each group). Figures 1-48 show the crop commodities identified as Group I; those from Group II are shown in Figures 49-60. Figures 61-68 are representative of the grape and soybean analyzed on the Rtx-50

confirmatory column. Figures 69-74 are representative chromatograms of standards on the GC/ECD. Figures 75-95 are representative of the GC/MS analysis of designated crop commodities.

NOTE: Representative complete analytical runs, including a run sheet and standard curve for both GC/ECD and GC/MS, can be found in Appendix 2.

## 11. Discussion

This method is intended to replace the current enforcement method for oxyfluorfen in crop commodities. As detailed in this report, a significant amount of analytical data was obtained in support of this method. This method effectively and consistently extracts and quantifies oxyfluorfen residues. The results show both methods of detection, GC/ECD and GC/MS, to be comparable in the recovery of oxyfluorfen residues. A comparison to the current enforcement method follows.

The extraction is similar to the current enforcement method, as published in PAM vol. II, method II. The techniques are comparable including similar solvents. Although blending is not used in this proposed enforcement method, the procedure using a shaker and centrifuge have proven to be effective in the other GOAL methods described below. Since the procedures are very similar we again conclude that the proposed extraction method provides an efficient mechanism in analyzing GOAL residues in crop commodities.

The extraction methodology found in this method is similar to that in both the "Oxyfluorfen (GOAL®) Meat and Fat Analytical Method" (TR 34-93-72) and the "Oxyfluorfen (GOAL®) Soil Analytical Method" (TR 34-94-47). It was shown to be quite efficient (~90%) at extracting GOAL residues via <sup>14</sup>C-labeled matrices in the previously mentioned methods. We conclude that the same outcome is expected of the crop residue samples.

Below are some points to be taken into consideration when performing this method:

- The Rtx-200 and Rtx-50 capillary columns should be conditioned according to the manufacturer's instructions. It was also necessary to prime the column before analysis by making several injections of standards.
- The purity of the petroleum ether can be critical. A self inspection of the solvent being used can be made by the following procedure: rotovap 200 ml of petroleum ether to dryness, then add 5 ml hexane. Inject the hexane extract on the GC with the same conditions as a sample. If the chromatograms give interfering peaks, this batch of petroleum ether cannot be used in the study. Therefore, a substitute for this product from another manufacturer or another batch from the same manufacturer must be used.

- In certain grape matrices, the pellet after the first centrifuge step may adhere very tightly to the side of the centrifuge bottle. If this does occur, it may be necessary to break it apart for the second extraction step to result in optimum performance.

In conclusion, this method is both robust and effective in the extraction of oxyfluorfen residues from crop commodities as written; and, it should be transferrable to other crops as needed.

## 12. Project Information

Title: Enforcement Residue Analytical Method  
for GOAL® Herbicide (oxyfluorfen) in Crop Commodities  
with GC/MS Confirmatory Detection

Sponsor: Rohm and Haas Company  
100 Independence Mall West, Philadelphia, PA, 19106-2399

### Testing Laboratories:

Centre Analytical Laboratories, Inc.  
3048 Research Drive, State College, PA 16801

McKenzie Laboratories  
3725 E. Atlanta Ave., Suite One, Phoenix, AZ 85040

Rohm and Haas Company  
727 Norristown Road, Spring House, PA 19477

Study Director: Dennis A. Martin, B.S.

Study Protocol: 34P-95-92, "Method Validation- GC/MS Confirmatory  
Detection in Oxyfluorfen (Goal®) Crop Residue Analytical  
Method, TR 34-94-150"

### Principal Investigators:

Dennis A. Martin, B.S.  
Kevin R. Kleinert, B.S.  
Rohm and Haas Company  
727 Norristown Road, Spring House, PA 19477

Shaozhi Zheng, M.S.  
Centre Analytical Laboratories, Inc.  
3048 Research Drive, State College, PA 16801

Melinda Lalko  
 McKenzie Laboratories  
 3725 E. Atlanta Ave., Suite One, Phoenix, AZ 85040

Archives:

The Rohm and Haas Company, Agricultural Products  
 Development Department Archives at Spring House

Data Handling:

Data system used was the Nelson Analytical, Inc. system with Version 1.2 of the 4400 Chromatography Software as its basis. This software provides the following functions:

- automatic collection and processing of chromatographic data
- detection and integration of peaks for chromatographic standards, samples and fortifications
- calculation and display of summary information for each standard, sample and fortification run
- display results from reports on the Laserjet printer
- editing and insertion of the summary data into IBM System 2000 data base

Report References:

<u>Title</u>	<u>Rohm and Haas TR No.</u>	<u>MRID No.</u>
Rohm and Haas Company Response to EPA CBTS May 12, 1994 Review of Residue Chemistry Data Supporting the Petition For Permanent Tolerances for Oxyfluorfen in Peanuts PP# 3F04229/FAP# 3H5674	34-95-74	Appendix 4 in 42793300
Oxyfluorfen Crop Residue Method	34-95-150	Appendix 8 in 42793300
Oxyfluorfen Residues in Peanut Nutmeat, Vine, Shell, and Hay	34-95-51	Appendix 9 in 42793300
Oxyfluorfen Residues in Apple	34-95-113	43794001
Oxyfluorfen Residues in Avocado	34-95-115	43794002
Oxyfluorfen Residues in Fig	34-95-116	43794003
Oxyfluorfen Residues in Pomegranate	34-95-117	43794004
Oxyfluorfen Residues in Kiwi	34-95-118	43794005
Oxyfluorfen Residues in Olive	34-95-119	43794006
Oxyfluorfen Residues in Artichoke	34-95-120	43794007
Oxyfluorfen Residues in Cherries	34-95-121	43794008



Oxyfluorfen Residues in Peach	34-95-114	N/A*
Oxyfluorfen Residues in Grape	34-95-104	N/A*
Oxyfluorfen Residues in Non-Dormant Grape (non-CA trials); Supplemental to TR 34-95-104	34-95-173	N/A*

\* N/A means that these reports were not submitted to the EPA as of the study completion date of this report.

Notebook References:

QZ04	Rohm and Haas No: 060662
QZ05	Rohm and Haas No: 061114
QZ06	Rohm and Haas No: 062196
QZ07	Rohm and Haas No: 062755
QZ08	Rohm and Haas No: 063646
KRK15	Rohm and Haas No: 064107

Reviewed and Edited By: Juan Chen, Rohm and Haas Company

Table 1. Apple Oxyfluorfen Fortification Summary <sup>1</sup>

RAR/ Sample No.	Fortification Level, ppm	% Recovery
94-0129-001	0.01	80.9
94-0130-001	0.01	81.1
94-0130-001	0.01	85.2
94-0130-001	0.01	88.8
95-1408-003	0.01	88.5
94-0152-001	0.01	76.5
94-0130-001	0.02	80.0
94-0130-001	0.05	90.8
94-0129-001	0.10	82.3
94-0130-001	0.10	89.7
94-0130-001	0.10	93.4
94-0152-001	0.10	40.2 <sup>2</sup>
95-1408-003	0.10	89.3
95-1408-003	0.50	91.7
94-0129-001	0.50	89.7
94-0130-001	0.50	83.5
94-0152-001	0.50	74.2
	<i>Average</i> n=16	85.4
	<i>Standard Deviation</i>	5.7

<sup>1</sup> All data reported in this table is taken from MRID 43794001.<sup>2</sup> Value outside 2 x SD, therefore not included in statistics.

Table 2. Artichoke Oxyfluorfen Fortification Summary<sup>1</sup>

RAR/ Sample No.	Fortification Level, ppm	% Recovery
94-0060-001	0.01	82.1
94-0061-001	0.01	99.0
94-0060-001	0.01	78.2
94-0060-001	0.01	68.3
94-0060-001	0.01	75.1
94-0060-001	0.02	104.0
94-0060-001	0.05	97.6
94-0060-001	0.10	104.0
94-0061-001	0.10	127.0
94-0060-001	0.50	99.4
94-0061-001	0.50	99.4
	<i>Average</i> n=11	94.0
	<i>Standard Deviation</i>	16.7

<sup>1</sup> All data reported in this table is taken from MRID 43794007.

Table 3. Avocado Oxyfluorfen Fortification Summary<sup>1</sup>

RAR/ Sample No.	Fortification Level, ppm	% Recovery
94-0141-001	0.01	51.4
94-0141-001	0.01	52.3
94-0141-001	0.01	71.7
94-0141-001	0.01	105.0
94-0141-001	0.01	108.0
94-0141-001	0.02	109.0
94-0141-001	0.02	111.0
94-0141-001	0.05	112.0
94-0141-001	0.10	100.0
94-0141-001	0.10	104.0
94-0141-001	0.50	92.8
94-0141-001	0.50	110.0
	<i>Average</i> n=12	93.9
	<i>Standard Deviation</i>	22.6

<sup>1</sup> All data reported in this table is taken from MRID 43794002.

Table 4. Cherry Oxyfluorfen Fortification Summary<sup>1</sup>

RAR/ Sample No.	Fortification Level, ppm	% Recovery
92-0041-002	0.01	73.8
92-0041-002	0.01	75.1
92-0041-002	0.01	76.1
92-0041-002	0.01	72.5
92-0042-002	0.01	75.9
92-0041-002	0.02	96.5
92-0041-002	0.05	109.0
92-0041-002	0.10	90.4
92-0042-002	0.10	101.0
92-0041-002	0.50	87.8
92-0042-002	0.50	93.6
	<i>Average</i> n=11	86.5
	<i>Standard Deviation</i>	12.6

<sup>1</sup> All data reported in this table is taken from MRID 43794008.

Table 5. Fig Oxyfluorfen Fortification Summary

RAR/ Sample No.	Fortification Level, ppm	% Recovery
94-0142-001	0.01 <sup>1</sup>	90.1
94-0142-001	0.01	88.0
94-0142-001	0.01	86.6
94-0142-001	0.01	84.1
94-0142-001	0.01	80.6
94-0142-001	0.02	103.0
94-0142-001	0.02	99.5
94-0142-001	0.05	97.8
94-0142-001	0.10	88.1
94-0142-001	0.1 <sup>1</sup>	98.5
94-0142-001	0.50	93.6
94-0142-001	0.5 <sup>1</sup>	88.8
	<i>Average</i> n=12	91.6
	<i>Standard Deviation</i>	6.9

<sup>1</sup> This data is taken from MRID 43794003.

All other data collected under 34P-95-30A.

Table 6. Grape Oxyfluorfen Fortification Summary<sup>1</sup>

Variety	RAR/ Sample No.	Fortification Level, ppm	% Recovery
Thompson	93-0012-001	0.01	70.5
	93-0012-001	0.01	71.8
	93-0012-001	0.01	74.7
	93-0012-001	0.01 <sup>1</sup>	90.7
	93-0012-001	0.01 <sup>1</sup>	91.8
	93-0012-001	0.01 <sup>1</sup>	85.6
	93-0012-001	0.01 <sup>1</sup>	84.4
	93-0012-001	0.05	76.8
	93-0012-001	0.10	81.3
	93-0012-001	0.1 <sup>1</sup>	78.1
	93-0012-001	0.1 <sup>1</sup>	84.4
	93-0012-001	0.1 <sup>1</sup>	78.1
	93-0012-001	0.1 <sup>1</sup>	85.5
	93-0012-001	0.50	83.6
	93-0012-001	0.50	88.0
	93-0012-001	0.5 <sup>1</sup>	87.8
	93-0012-001	0.5 <sup>1</sup>	85.5
	93-0012-001	0.5 <sup>1</sup>	91.1
Concord	92-0089-005	0.01 <sup>1</sup>	109.0
	92-0089-005	0.01 <sup>1</sup>	95.8
	92-0089-005	0.1 <sup>1</sup>	107.0
	92-0089-005	0.1 <sup>1</sup>	77.0
	92-0089-005	0.5 <sup>1</sup>	104.0
	92-0089-005	0.5 <sup>1</sup>	108.0

*Data continued on the next page.*

<sup>1</sup> This data is taken from Appendix 8 in MRID 42793300.

All other data on this page is taken from TR 34-95-104, under protocol 34P-95-65A, MRID not available yet.

Table 6. Grape Oxyfluorfen Fortification Summary<sup>1</sup>  
(continued)

Variety	RAR/ Sample No.	Fortification Level, ppm	% Recovery	
Concord	92-0089-005	0.01	88.8	
	92-0089-005	0.10	114.0	
Cabernet	92-0100-005	0.01	81.0	
		0.01	86.7	
		0.01	69.6	
		0.02	88.0	
		0.02	96.0	
		0.10	106.0	
Niagra	92-0072-005	0.01	86.9	
		0.01	88.0	
		0.01	93.1	
		0.02	90.5	
		0.02	108.0	
		0.05	106.0	
		<i>Average</i>	38	89.3
		<i>Standard Deviation</i>		11.6

<sup>1</sup> All other data on this page is taken from TR 34-95-173, under protocol 34P-95-65A, MRID not available yet.



Table 7. Kiwi Oxyfluorfen Fortification Summary<sup>1</sup>

RAR/ Sample No.	Fortification Level, ppm	% Recovery
94-0146-001	0.01 <sup>2</sup>	64
94-0146-001	0.01	76.4
94-0146-001	0.01	86.5
94-0146-001	0.01	87.4
94-0146-001	0.01	93.8
94-0146-001	0.02	97.5
94-0146-001	0.02	101.0
94-0146-001	0.05	103.0
94-0146-001	0.05	104.0
94-0146-001	0.10	96.9
94-0146-001	0.50	84.2
94-0146-001	0.50	91.2
	<i>Average</i> n=11	92.9
	<i>Standard Deviation</i>	8.7

<sup>1</sup> All data reported in this table is taken from MRID 43794005.

<sup>2</sup> Value outside 2 x SD, therefore not included in statistics.

Table 8. Olive Oxyfluorfen Fortification Summary

RAR/ Sample No.	Fortification Level, ppm	% Recovery
94-0172-001	0.01 <sup>1</sup>	91.5
94-0172-001	0.01	81.5
94-0172-001	0.01	75.1
94-0172-001	0.01	74.6
94-0172-001	0.01	82.4
94-0172-001	0.02	107.0
94-0172-001	0.02	106.0
94-0172-001	0.05	98.6
94-0172-001	0.10	74.5
94-0172-001	0.1 <sup>1</sup>	105.0
94-0172-001	0.50	91.6
94-0172-001	0.5 <sup>1</sup>	93.8
	<i>Average</i> n=12	90.1
	<i>Standard Deviation</i>	12.4

<sup>1</sup> This data is taken from MRID 43794006.

All other data collected under 34P-95-33A.

Table 9. Peach Oxyfluorfen Fortification Summary <sup>1</sup>

RAR/ Sample No.	Fortification Level, ppm	% Recovery
94-0117-001	0.01	70.0
94-0117-001	0.01	85.0
94-0117-001	0.01	86.1
94-0117-001	0.01	90.9
94-0117-001	0.01	92.9
94-0117-001	0.01 <sup>2</sup>	102.0
95-0196-003	0.01	88.3
94-0117-001	0.02	73.0
94-0117-001	0.02	87.5
94-0117-001	0.02	91.5
94-0117-001	0.05	78.2
94-0117-001	0.05	87.4
94-0117-001	0.05	93.4
94-0117-001	0.05 <sup>2</sup>	103.0
95-0196-003	0.05	89.8
94-0117-001	0.10	72.1
94-0117-001	0.10	86.5
94-0117-001	0.10	93.7
94-0117-001	0.10 <sup>2</sup>	88.7
94-0117-001	0.50	81.7
94-0117-001	0.50	88.5
94-0117-001	0.50	91.9
94-0117-001	0.50 <sup>2</sup>	75.2
95-0196-003	0.10	83.2
	<i>Average</i>	n=24 86.7
	<i>Standard Deviation</i>	8.5

<sup>1</sup> All data reported is taken from TR 34-95-114 (MRID not yet available), under protocols 34P-95-35A, 34P-95-51A, except those noted as below.

<sup>2</sup> This data was not collected under the protocols noted above or reported in TR 34-95-114. It was collected as an evaluation of the method.

Table 10. Peanut Nutmeat Oxyfluorfen Fortification Summary

RAR/ Sample No.	Fortification Level, ppm	% Recovery
93-0142-005	0.01 <sup>1</sup>	115.0
93-0142-005	0.01 <sup>1</sup>	131.0 <sup>2</sup>
93-0142-005	0.01 <sup>3</sup>	79.0
93-0142-005	0.01 <sup>3</sup>	86.0
93-0142-005	0.01 <sup>3</sup>	96.3
93-0142-005	0.01 <sup>3</sup>	84.8
94-0176-010	0.01 <sup>1</sup>	105.0
93-0142-005	0.02 <sup>3</sup>	102.0
93-0142-005	0.02 <sup>3</sup>	95.5
93-0142-005	0.02 <sup>3</sup>	79.5
93-0142-005	0.02 <sup>3</sup>	84.5
93-0142-005	0.05 <sup>3</sup>	113.0
93-0142-005	0.1 <sup>3</sup>	76.4
93-0157-004	0.1 <sup>1</sup>	70.1
94-0176-010	0.1 <sup>1</sup>	80.6
93-0142-005	0.5 <sup>3</sup>	90.6
93-0157-004	0.5 <sup>1</sup>	99.4
94-0176-010	0.5 <sup>1</sup>	86.6
	<i>Average</i> n=17	90.8
	<i>Standard Deviation</i>	12.9

<sup>1</sup> Data reported taken from Appendix 9 in MRID 42793300.

<sup>2</sup> Value outside 2 x SD, therefore not included in statistics.

<sup>3</sup> Data reported taken from Appendix 4 in MRID 42793300.

Table 11. Pomegranate Oxyfluorfen Fortification Summary

RAR/ Sample No.	Fortification Level, ppm	% Recovery
94-0143-001	0.01 <sup>1</sup>	76.4
94-0143-001	0.01	87.2
94-0143-001	0.01	80.1
94-0143-001	0.01	81.9
94-0143-001	0.01	70.6
94-0143-001	0.02	103.0
94-0143-001	0.02	109.0
94-0143-001	0.05	107.0
94-0143-001	0.10	92.7
94-0143-001	0.1 <sup>1</sup>	100.0
94-0143-001	0.50	104.0
94-0143-001	0.5 <sup>1</sup>	84.0
	<i>Average</i> n=12	91.3
	<i>Standard Deviation</i>	13.0

<sup>1</sup> This data is taken from MRID 43794004.

All other data collected under 34P-95-31A.

Table 12. Soybean Oxyfluorfen Fortification Summary

RAR/ Sample No.	Fortification Level, ppm	% Recovery
92-0107-014	0.01	108.0
92-0107-014	0.01	119.0
94-0136-001	0.01	119.0
94-0136-001	0.01	108.0
92-0107-014	0.10	98.5
92-0107-014	0.10	100.0
94-0136-001	0.10	86.3
94-0136-001	0.10	106.0
92-0107-014	0.50	121.0
92-0107-014	0.50	108.0
94-0136-001	0.50	91.6
94-0136-001	0.50	93.2
	<i>Average</i> n=12	104.9
	<i>Standard Deviation</i>	11.3

<sup>1</sup> This data is taken from Appendix 8 in MRID 42793300.

Table 13. Peanut Hay Oxyfluorfen Fortification Summary

RAR/ Sample No.	Fortification Level, ppm	% Recovery
93-0142-003	0.02 <sup>1</sup>	87.4
93-0142-003	0.02 <sup>2</sup>	97.5
93-0142-003	0.02 <sup>2</sup>	96.5
93-0157-002	0.02 <sup>1</sup>	111.0
93-0157-002	0.02 <sup>2</sup>	80.5
93-0157-002	0.05 <sup>2</sup>	97.0
93-0157-002	0.05 <sup>2</sup>	92.8
93-0142-003	0.05 <sup>2</sup>	99.0
93-0142-003	0.1 <sup>1</sup>	82.0
93-0157-002	0.5 <sup>1</sup>	78.0
	<i>Average</i> n=10	92.2
	<i>Standard Deviation</i>	10.2

<sup>1</sup> Data reported taken from Appendix 9 in MRID 42793300.

<sup>2</sup> Data reported taken from Appendix 4 in MRID 42793300.

Table 14. Peanut Shell Oxyfluorfen Fortification Summary

RAR/ Sample No.	Fortification Level, ppm	% Recovery
93-0142-005	0.02 <sup>1</sup>	102.0
93-0142-005	0.02 <sup>3</sup>	98.5
93-0157-004	0.02 <sup>1</sup>	65.5 <sup>2</sup>
93-0157-004	0.02 <sup>3</sup>	89.5
93-0157-004	0.05 <sup>3</sup>	95.2
93-0142-005	0.02 <sup>3</sup>	94.5
93-0142-005	0.05 <sup>3</sup>	98.0
93-0142-005	0.05 <sup>3</sup>	94.0
93-0142-005	0.1 <sup>1</sup>	89.6
93-0157-004	0.1 <sup>3</sup>	86.4
93-0157-004	0.5 <sup>1</sup>	105.0
	<i>Average</i> n=10	95.3
	<i>Standard Deviation</i>	5.8

<sup>1</sup> Data reported taken from Appendix 9 in MRID 42793300.

<sup>2</sup> Value outside 2 x SD, therefore not included in statistics.

<sup>3</sup> Data reported taken from Appendix 4 in MRID 42793300.



Table 15. Peanut Vine Oxyfluorfen Fortification Summary

RAR/ Sample No.	Fortification Level, ppm	% Recovery
93-0142-001	0.02 <sup>2</sup>	87.5
93-0142-001	0.02 <sup>2</sup>	82.5
93-0142-001	0.02 <sup>1</sup>	82.1
93-0142-001	0.02 <sup>1</sup>	87.6
93-0142-001	0.05 <sup>2</sup>	87.4
93-0142-001	0.05 <sup>2</sup>	95.2
93-0142-001	0.05 <sup>2</sup>	98.0
93-0142-001	0.1 <sup>1</sup>	73.2
93-0142-001	0.1 <sup>1</sup>	64.6
93-0142-001	0.5 <sup>1</sup>	68.2
93-0142-001	0.5 <sup>1</sup>	74.0
	<i>Average</i> n=11	81.8
	<i>Standard Deviation</i>	10.8

<sup>1</sup> Data reported taken from Appendix 9 in MRID 42793300.

<sup>2</sup> Data reported taken from Appendix 4 in MRID 42793300.

Table 16. Grape Oxyfluorfen Fortification Summary <sup>1</sup>  
Rtx-50 Confirmatory Column

Variety	RAR/ Sample No.	Fortification Level, ppm	% Recovery	
Thompson	93-0012-001	0.01	100.0	
	93-0012-001	0.01	77.9	
	93-0012-001	0.10	75.3	
	93-0012-001	0.10	97.3	
	93-0012-001	0.50	101.0	
	93-0012-001	0.50	95.2	
Concord	92-0089-005	0.01	92.5	
	92-0089-005	0.01	84.2	
	92-0089-005	0.10	97.8	
	92-0089-005	0.10	83.9	
	92-0089-005	0.50	110.0	
	92-0089-005	0.50	96.5	
		<i>Average</i>	n=12	92.6
		<i>Standard Deviation</i>		10.3

<sup>1</sup> This data is taken from Appendix 8 in MRID 42793300.

Table 17. Soybean Oxyfluorfen Fortification Summary  
Rtx-50 Confirmatory Column

RAR/ Sample No.	Fortification Level, ppm	% Recovery
92-0107-014	0.01	54 *
92-0107-014	0.01	50 *
94-0136-001	0.01	106.0
94-0136-001	0.01	113.0
92-0107-014	0.10	98.7
92-0107-014	0.10	90.8
94-0136-001	0.10	100.0
94-0136-001	0.10	124.0
92-0107-014	0.50	110.0
92-0107-014	0.50	97.7
94-0136-001	0.50	104.0
94-0136-001	0.50	95.6
	<i>Average</i> n=10	104.0
	<i>Standard Deviation</i>	9.7

<sup>1</sup> This data is taken from Appendix 8 in MRID 42793300.

\* Value outside 2 X SD, therefore not included in statistics.

Table 18. Oxyfluorfen GC/MS Confirmatory Data Summary <sup>1</sup>

Sample Type	RAR/ Sample No.	Fortification Level, ppm	% Recovery	
			GC/MS	ECD
grape	92-0100-005	0.025	93.0	99.0
grape	92-0100-005	0.05	86.0	82.9
grape	92-0100-005	0.125	108.0	91.2
		<i>Average</i>	95.7	91.0
peach	94-0117-001	0.025	85.0	91.4
peach	94-0117-001	0.05	82.0	79.4
peach	94-0117-001	0.125	102.0	89.6
		<i>Average</i>	89.7	86.8
peanut nutmeat	93-0157-004	0.025	109.0	109.0
peanut nutmeat	93-0157-004	0.05	81.0	163.0
peanut nutmeat	93-0157-004	0.125	87.0	116.0
		<i>Average</i>	92.3	129.3
peanut shell	93-0157-004	0.025	104.0	104.0
peanut shell	93-0157-004	0.05	117.0	107.0
peanut shell	93-0157-004	0.125	103.0	106.0
		<i>Average</i>	108.0	105.7

<sup>1</sup> This data was collected and reviewed under protocol 34-95-92.

LOQ=0.025 ppm


LOD=0.008 ppm

Table 19. Oxyfluorfen Data Summary by Crop

Crop	Group <sup>1</sup>	Detection	Average		
			Recovery	LOQ	LOD
apple	I	GC/ECD	85.4	0.01 ppm	0.003 ppm
artichoke	I	GC/ECD	94	0.01 ppm	0.003 ppm
avocado	I	GC/ECD	93.9	0.01 ppm	0.003 ppm
cherry	I	GC/ECD	86.5	0.01 ppm	0.003 ppm
fig	I	GC/ECD	91.6	0.01 ppm	0.003 ppm
grape	I	GC/ECD	89.3	0.01 ppm	0.003 ppm
kiwi	I	GC/ECD	92.9	0.01 ppm	0.003 ppm
olive	I	GC/ECD	90.1	0.01 ppm	0.003 ppm
peach	I	GC/ECD	86.7	0.01 ppm	0.003 ppm
peanut nutmeat	I	GC/ECD	90.9	0.01 ppm	0.003 ppm
pomegranate	I	GC/ECD	91.3	0.01 ppm	0.003 ppm
soybean	I	GC/ECD	104.9	0.01 ppm	0.003 ppm
peanut hay	II	GC/ECD	92.2	0.02 ppm	0.007 ppm
peanut shell	II	GC/ECD	95.3	0.02 ppm	0.007 ppm
peanut vine	II	GC/ECD	81.8	0.02 ppm	0.007 ppm

<sup>1</sup> Refers to section 7 crop grouping for analysis.

Figure 1. Apple Control  
95-1408-003




RUN# 1966 JUN 30, 1995 01:35:03  
**APPLE**  
 SAMPLE NAME: **95-1408** SAMPLE# 23 **Control**  
 METHOD NAME: **M:GDALEXT.MET**  
 G FINAL WT **5** VOL INJ **1** ML FINAL VOL **5** OIL FACTOR **1** MG INJ **1**

RT	HEIGHT	TYPE	WIDTH	HEIGHTX
15.348	5803	BB	.117	8.14573
18.526	42884	PB	.196	91.85424

TOTAL HEIGHT= 46687  
 MUL FACTOR=1.0000E+00

Figure 2. Apple Fortification 0.01 PPM  
95-1408-003




RUN# 1969 JUN 30, 1995 03:48:50  
**APPLE** fortification  
 SAMPLE NAME: 0.01ppm SAMPLER 26 95-1408  
 METHOD NAME: M:GOALEXT.MET  
 G FINAL WT 5 VOL INJ 1ul ML FINAL VOL 5 DIL FACTOR 1 MG INJ 1

HEIGHTX	RT	HEIGHT	TYPE	WIDTH	HEIGHTX
	15.345	1441	BB	.158	3.66331
	16.770	2657	BB	.123	6.75465
	18.520	35238	PB	.187	89.58205

TOTAL HEIGHT= 39336  
 MUL FACTOR=1.0000E+00

Figure 3. Apple Fortification 0.10 PPM  
95-1408-003



RUN# 1970 JUN 30, 1995 04:33:27  
**APPLE** fortification  
 SAMPLE NAME: 0.1ppm SAMPLE# 27 95-1408  
 METHOD NAME: M:60ALEXT.MET  
 G FINAL WT 5 VOL INJ 1ul ML FINAL VOL 5 OIL FACTOR 1 MG INJ 1

HEIGHTX	RT	HEIGHT	TYPE	WIDTH	HEIGHTX
	15.344	1374	BP	.147	2.20315
	16.767	24114	PB	.124	40.06979
	18.518	34692	PB	.105	57.64707


TOTAL HEIGHT= 60100  
 MUL FACTOR=1.0000E-00



Figure 4. Apple Fortification 0.50 PPM  
95-1408-003

RUN # 1972 JUN 30, 1995 06:02:38

START



RUN# 1972 JUN 30, 1995 06:02:38

APPLE fortification  
 SAMPLE NAME: 0.5 ppm SAMPLE# 29 95-1408  
 METHOD NAME: M:GOALEXT.MET  
 G FINAL WT. 5 UOL INJ. 10 ML FINAL UOL 5 DIL FACTOR 5 MG INJ. 0.2

RT	HEIGHT	TYPE	WIDTH	HEIGHTX
16.767	24710	PB	.124	65.70934
18.518	12895	PB	.164	34.29066

TOTAL HEIGHT= 37605  
 MUL FACTOR=1.0000E+00


### Figure 5. Artichoke Control 94-0060-001

Data file: L0666 RAR number: 94-0060  
 Method file: GOALVn Sample No: 06121 <sup>117195</sup>  
 Type: SAMPLE Component: ARTICHOKE

Sample Name: NA Cal. Curve: 07/06/95  
 Date: 6 Jul 1995 19:12 Method: GOALVn Analyst: SS <sup>52</sup>  
 Interface: 711 Cycles: 5 Channels: A

Instrument: VARIAN 3800.  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
 Prog Slope (0 or Linear): NA Inj Port Temp: 255  
 Flowrate/Gas: 3.5ml/min Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.31 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.
13.40	GOAL	0.000E+00	.000E+00	5.00	5.00


Figure 6. Artichoke Fortification 0.01 PPM  
94-0060-001

Data file: L0688 RAR number: 94-0060  
Method file: GOALVn Sample No: 001  
Type: FORTIFICATION Component: ARTICHOKE

Sample Name: NA Cal. Curve: 07/06/95  
Date: 6 Jul 1995 20:28 Method: GOALVn Analyst: AS 62  
Interface: 711 Cycles: 8 Channel: A

Instrument: VARIAN 3500,  
Column: KIX-200, 0.32mm ID, 1.00m df Column Length: 60 Meters  
Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
Prog Slope (S or Linear): NA Inj Port Temp: 255  
Flowrate/Gas: 3.5ml/min Split Ratio: NA  
Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
Plot range: 50 millivolts (-.24 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp. Wt.	ug Added
13.33	GOAL	2.770E+08	.224E+03	5.00	5.00	0.050

Figure 7. Artichoke Fortification 0.10 PPM  
94-0060-001

Data file: L0489 RAR number: 94-0060  
 Method file: GOALVA Sample No: 001  
 Type: FORTIFICATION Component: ARTICHOKE

---


Sample Name: NA Cal. Curve: 07/06/95  
 Date: 6 Jul 1995 20:53 Method: GOALVA Analyst: SS 92  
 Interface: 711 Cycles: 9 Channel: A

---

Instrument: VARIAN 1500,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
 Prog Slope (0 or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 2.5ml/minHE Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

---

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.37 uv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp wt.	ug Added
13.33	GOAL	2.960E+01	.243E+04	5.00	5.00	0.500



### Figure 9. Avocado Control 94-0141-001

```

Data file: L1985                      RAN number: 94-0141
Method file: GOALVn                   Sample No: 001
Type: SAMPLE                          Component: AVOCADO
-----
Sample Name: NA                        Cal. Curve: 07/20/98
Date: 20 Jul 1995 23:35               Method: GOALVn      Analyst: ST
Interface: 711                         Cycle#: 8           Channel#: A
-----
Instrument: VARIAN 3500,
Column: RTX-200, 0.32mm ID, 1.0um df  Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C      Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C           End Time-Temp (deg-min): 250 C-17
Prog Slope (S or Linear): NA          Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min              Split Ratio: NA
Det 1-Type & Temp: ECD/300C           Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 50 millivolts (-.24 mv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp. Wt.
5.00	GOAL	0.000E+00	.000E+00	5.00	5.00


### Figure 10. Avocado Fortification 0.01 PPM 94-0141-001

Data file: L1989  
 Method file: GOALVA  
 Type: FORTIFICATION  
 RAR number: 94-0141  
 Sample No: 001  
 Component: AVOCADO

Sample Name: NA  
 Date: 21 Jul 1995 01:16  
 Interface: 711  
 Method: GOALVA  
 Cycled: 9  
 Cal. Curve: 07/20/95  
 Analyst: SS  
 Channel: A

Instrument: VARIAN 3500,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog slope (S or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/min Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.5 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Sample Wt.	ug Added
13.31	GOAL	3.110E+00	.248E+03	5.00	5.00	0.050


Figure 11. Avocado Fortification 0.10 PPM  
94-0141-001

Data file: L118A9 RAR number: 94-0141  
 Method file: GOALVA Sample No: 001  
 Type: FORTIFICATION Component: AVOCADO

Sample Name: NA Cal. Curve: 07/13/98  
 Date: 12 Jul 1998 16:14 Method: GOALVA Analyst: SE 57  
 Interface: 711 Cycles: 9 Channel: A

Instrument: VARIAN 3500,  
 Column: RTX-200, 0.32mm ID, 1.0um d2 Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog slope (S or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/minHE Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.3 uv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Emp Wt.	ug Added
5.00	GOAL	4.240E+01	.297E+04	5.00	5.00	0.500



### Figure 12. Avocado Fortification 0.50 PPM 94-0141-001

```

Data file: 219816                      RAR number: 94-0141
Method file: GOALVM                     Sample No: 001
Type: XXXXXXXXXX                   Component: AVOCADO
-----
Sample Name: NA                          Cal. Curve: 07/20/95
Date: 21 Jul 1995 04:13 Method: GOALVM   Analyst: SS
Interface: 711                          Cycles: 16      Channel: A
-----
Instrument: VARIAN 3500.
Column: RTX-200, 0.32mm ID, 1.0um df   Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C       Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C           End Time-Temp (deg-min): 250 C-17
Prog Slope (S or Linear): NA          Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min               Split Ratio: NA
Det 1-Type & Temp: ECD/300C           Det 2-Type & Temp: NA
-----
Plot time: 0 to 30 minutes
Plot range: 50 millivolts (-.19 mv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ul)	Resp. WC.	ug Added
13.31	GOAL	1.670E+01	.137E+04	50.0	5.00	2.50






Figure 15. Cherry Fortification 0.10 PPM  
92-0041-002

Data file: L8889  
Method file: GOALVM  
Type: FORTIFICATION  
RRR number: 92-0041  
Sample No: 003  
Component: CHERRY

Sample Name: NA  
Date: 5 Jul 1995 19:50  
Interface: 711  
Method: GOALVM  
Cycles: 9  
Cal. Curve: 07/05/95  
Analyte: SE  
Channel: A

Instrument: VARIAN 3800,  
Column: MIX-200, 0.32mm ID, 1.0um dF Column Length: 60 Meters  
Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
Prog Slope (S or Linear): NA Inj Port Temp: 265  
Flowrate/Gas: 3.5ml/min Split Ratio: NA  
Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
Plot range: 50 millivolts (-.20 uv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp. Wt.	ug Added
13.33	GOAL	2.610E+01	.388E+04	5.00	5.00	0.500






### Figure 18. Fig Fortification 0.01 PPM 94-0142-001

Data file: L118A16 RAR number: 94-0142  
 Method file: GOALVA Sample No: 001  
 Type: Fortification Component: F20

Sample Name: NA Cal. Curve: 07/12/95  
 Date: 12 Jul 1995 19:10 Method: GOALVA Analyst: SG-J  
 Interface: 711 Cycles: 16 Channel: A

Instrument: VARIAN 3500.  
 Column: KIX-200, 0.22mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 280 C-17  
 Prog slope (S or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/min Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.25 uv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
13.32	GOAL	3.330E+00	.277E+03	5.00	5.00	0.050

### Figure 19. Fig Fortification 0.10 PPM 94-0142-001

```

Data file: L118A17                      RAR number: 94-0142
Method file: G0ALVn                      Sample No: 001
Type: FORTIFICATION                      Component: FIG
-----
Sample Name: NA                           Cal. Curve: 07/12/95
Date: 12 Jul 1995 19:35 Method: G0ALVn   Analyst: SS/2
Interface: 711                            Cycled: 17      Channel: A
-----
Instrument: VARIAN 3500,
Column: RTX-200, 0.32mm ID, 1.0um df   Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C       Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C           End Time-Temp (deg-min): 250 C-17
Prog Slope (S or Linear): NA          Inj Port Temp: 245
Flowrate/Gas: 3.5ml/min              Split Ratio: NA
Det 1-Type & Temp: ECD/300C           Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 50 millivolts (-.3 uv offset)
    
```




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp. Wt.	ug Added
5.00	GOAL	3.268E+01	.276E+04	5.00	5.00	0.500



### Figure 20. Fig Fortification 0.50 PPM 94-0142-001

```

Data file: L118A18                      RAR number: 94-0142
Method file: GOALVA                      Sample No: 001
Type: FORTIFICATION                   Component: FID
-----
Sample Name: NA                           Cal. Curve: 07/13/95
Date: 13 Jul 1995 20:00 Method: GOALVA    Analyst: SE-SE
Interface: 711                            Cycles: 18      Channel: A
-----
Instrument: VARIAN 3500,
Column: RTX-300, 0.32mm ID, 1.0um d2    Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C        Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C             End Time-Temp (deg-min): 250 C-17
Prog slope (0 or Linear): NA            Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min                Split Ratio: NA
Det 1-Type & Temp: ECD/100C             Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 50 millivolts (-.37 mv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp. Wt.	ug Added
13.32	GOAL	1.280E+01	.102E+04	50.0	5.00	2.50

### Figure 21. Grape Control 93-0012-001

Data file: J05J6 BAR number: 93-0012  
 Method file: GOALVn Sample No: 001  
 Type: SAMPLE Component: THOMPSON

---


Sample Name: NA Cal. Curve: 01/05/95  
 Date: 5 Jan 1995 21:29 Methods: GOALVn Analyst: SZ 52  
 Interface: 711 Cycle#: 6 Channel#: A

---

Instrument: VARIAN 3300,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
 Prog Slope (# or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/min. Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

---

Plot times: 0 to 30 minutes  
 Plot range: 360 millivolts ( 31 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.
5.26	GOAL	0.000E+00	.000E+00	5.00	5.00


Figure 22. Grape Fortification 0.01 PPM  
93-0012-001

Data file: J1337 RAR number: 93-0012  
Method file: GOALVM Sample No: 001  
Type: FORTIFICATION Component: THOMPSON

Sample Name: NA Cal. Curve: 01/14/95  
Date: 14 Jan 1995 17:49 Method: GOALVM Analyst: SZ 42  
Interface: 715 Cycle#: 7 Channel#: A

Instrument: VARIAN 3500,  
Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
Prog Slope (# or Linear): NA Inj Port Temp: 265  
Flowrate/Gas: 3.5ml/min. Split Ratio: NA  
Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
Plot range: 360 millivolts ( 31 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp. Vt.	ug Added
15.30	GOAL	3.610E+01	.231E+04	5.00	5.00	0.050

### Figure 23. Grape Fortification 0.10 PPM 93-0012-001

Data file: J05J8 RAW number: 93-0012  
 Method file: GOALVM Sample No: 001  
 Type: FORTIFICATION Component: THOMPSON

---


Sample Name: NA Cal. Curve: 01/05/95  
 Date: 5 Jan 1995 22:42 Method: GOALVM Analyst: SZ SZ  
 Interface: 711 Cycles: 8 Channel: A

---

Instrument: VARIAN 3300,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
 Prog. Slope (# or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/min. Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

---

Plot times: 0 to 30 minutes  
 Plot range: 360 millivolts ( 31 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
15.25	GOAL	2.780E+02	.160E+05	5.00	5.00	0.500

Figure 24. Grape Fortification 0.50 PPM  
93-0012-001

Data file: J05J9  
Method file: GDALVn  
Type: FORTIFICATION

RAN number: 93-0012  
Sample No: 001  
Component: THOMPSON

Sample Name: NA  
Date: 5 Jan 1995 23:19  
Interface: 711


Method: GDALVn  
Cycles: 9

Cal. Curve: 01/05/95  
Analyst: SZ JZ  
Channel: A

Instrument: VARIAN 3500,  
Column: RTX-200, 0.32mm ID, 1.0um df  
Start Temp-Time (deg-min): 215 C  
Program Rate (deg/min): 5C  
Flowrate/Gas: 3.5ml/min.Det 1-Type & Temp: ECD/300C

Column Length: 60 Meters  
Ramp Hold (deg-min): 1  
End Time-Temp (deg-min): 250 C-25  
Inj Port Temp: 265  
Split Ratio: NA  
Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
Plot range: 360 millivolts ( 31 mv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume * Samp (ml)	Samp Wt.	ug Added
15.25	GDAL	1.590E+02	.902E+04	50.0	5.00	2.50

### Figure 25. Kiwi Control 94-0146-001

```

Data file: L26JAS                      RAR number: 94-0146
Method file: GOALVn                     Sample No: 001
Type: SAMPLE                          Component: KIWI
-----
Sample Name: NA                          Cal. Curve: 07/22/95
Date: 22 Jul 1995 10:28 Method: GOALVn  Analyst: JB
Interface: 711                          Cycle#: 5      Channel#: A
-----
Instrument: VARIAN 3800,
Column: RTX-200, 0.32mm ID, 1.0um dF   Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C       Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C           End Time-Temp (deg-min): 250 C-17
Prog Slope (S or Linear): NA          Inj Port Temp: 255
Flowrate/Gas: 3.5ml/min               Split Ratio: NA
Det 1-Type & Temp: ECD/300C           Det 2-Type & Temp: NA
-----
Plot time: 0 to 30 minutes
Plot range: 50 millivolts (-.25 uv offset)
    
```




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.
13.33	GOAL	0.000E+00	.000E+00	5.00	5.00

### Figure 26. Kiwi Fortification 0.01 PPM 94-0146-001

```

Data file: 130726                      RAN number: 94-0146
Method file: GOALVn                     Sample No: 001
Type: FORTIFICATION                   Component: KIWI
-----
Sample Name: NA                          Cal. Curve: 07/22/95
Date: 22 Jul 1995 18:54                 Method: GOALVn      Analyst: JB
Interface: 711                           Cycle#: 6           Channel#: A
-----
Instrument: VARIAN 3500,
Column: RTX-200, 0.32mm ID, 1.0um dE   Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C       Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C            End Time-Temp (deg-min): 250 C-17
Prog slope (0 or Linear): NA           Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min               Split Ratio: NA
Det 1-Type & Temp: ECD/300C           Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 50 millivolts (-.25 uv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Resp Wt.	ug Added
15.32	GOAL	1.390E+00	.281E+03	5.00	5.00	0.050








### Figure 29. Olive Control 94-0172-001

```

Data file:  L07813                      RAR number: 94-0172
Method file: GOALVn                      Sample No: 0013
Type:       SAMPLE                       Component: OLIVE
-----
Sample Name: NA                           Cal. Curve: 07/07/95
Date: 7 Jul 1995 20:52   Method: GOALVn   Analyst: SE SZ
Interface: 711           Cycles: 13       Channels: A
-----
Instrument: VARIAN 1500.
Column: RTX-200, 0.32mm ID, 1.0um d2   Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C       Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C            End Time-Temp (deg-min): 250 C-17
Prog Slope (# or Linear): NA           Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min               Split Ratio: NA
Det 1-Type & Temp: ECD/300C            Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 50 millivolts (-.25 mv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.
5.00	GOAL	0.000E+00	.000E+00	5.00	5.00






Figure 32. Olive Fortification 0.50 PPM  
94-0172-001

Data file: 107218 RAR number: 94-0172  
Method file: GOALVn Sample No: 001  
Type: FORTIFICATION Component: OLIVE

Sample Name: NA Cal. Curve: 07/07/98  
Date: 7 Jul 1995 12:58 Method: GOALVn Analyst: ASZ  
Interface: 711 Cycle#: 18 Channel#: A


Instrument: VARIAN 3500.  
Column: RTX-300, 0.32mm ID, 1.0um df Column Length: 60 Meters  
Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
Prog Slope (0 or Linear): NA Inj Port Temp: 250  
Flowrate/Gas: 3.5ml/min Split Ratio: NA  
Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
Plot range: 50 millivolts (-.25 uv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp Wt.	ug Added
13.33	GOAL	1.320E+01	.106E+04	50.0	5.00	2.50

Figure 33. Peach Control  
94-0117-001




RUN# 1989 JUL 7, 1995 13:47:03  
FRUIT  
SAMPLE NAME: 94-0117-001 SAMPLE# 6 CONTROL  
METHOD NAME: M:GOALEXT.MET  
G FINAL WT\_\_\_ UOL INJ\_\_\_ ML FINAL UOL\_\_\_ DIL FACTOR\_\_\_ MG INJ\_\_\_

RT	HEIGHT	TYPE	WIDTH	HEIGHTX
15.334	4164	BB	.118	7.55470
17.250	812	PU	.143	1.47320
18.500	50142	PB	.199	90.97210

TOTAL HEIGHT= 55118  
MUL FACTOR=1.0000E+00

Figure 34. Peach Fortification 0.01 PPM  
94-0117-001




RUN# 1992 **Fortification** JUL 7, 1995 16:01:09  
**FRUIT** SAMPLE NAME: **0.01ppm** SAMPLE# 9 **94-0117-001**  
 METHOD NAME: M:GOALEXT.ME  
 G FINAL WT **5** UOL INJ **1** ML FINAL UOL **5** DIL FACTOR **1** MG INJ **1**

RT	HEIGHT	TYPE	WIDTH	HEIGHTX
15.337	4040	PB	.115	7.08685
16.759	3154	PB	.122	5.53265
17.249	765	PB	.129	1.34194
18.510	49048	PB	.199	86.03853

TOTAL HEIGHT= 57007  
 MUL FACTOR=1.0000E+00

Figure 35. Peach Fortification 0.10 PPM  
94-0117-001

RUN # 1993 JUL 7, 1995 16:45:46  
START



STOP

RUN# 1993 JUL 7, 1995 16:45:46  
**FRUIT** **Fortification**  
 SAMPLE NAME: 0.10ppm SAMPLE# 10 **94-0117-001**  
 METHOD NAME: M:GOALEXT.MET  
 G FINAL WT 5 UOL INJ 1ul ML FINAL VOL 5 DIL FACTOR 1 MG INJ 1

HEIGHTX	RT	HEIGHT	TYPE	WIDTH	HEIGHTX
	15.336	3362	BB	.133	4.00271
	16.759	28139	PB	.125	33.50160
	17.250	751	PB	.122	.89412
	18.509	51741	PB	.203	61.60157


TOTAL HEIGHT= 83993  
MUL FACTOR=1.0000E+00



Figure 36. Peach Fortification 0.50 PPM  
94-0117-001

RUN # 1996 JUL 7, 1995 18:59:35

START



RUN# 1996 JUL 7, 1995 18:59:35

FRUIT

Fortification

SAMPLE NAME: 0.5 ppm SAMPLE# 13 94-0117-001

METHOD NAME: M:GOREXT.MET

G FINAL WT 5 VOL INJ 1ml ML FINAL VOL 5 DIL FACTOR 5 MG INJ 0.2

HEIGHTX	RT	HEIGHT	TYPE	WIDTH	HEIGHTX
	16.757	27553	PB	.124	59.21939
	18.505	18974	PB	.173	40.78062


TOTAL HEIGHT= 46527

MUL FACTOR=1.0000E+00

### Figure 37. Peanut Nutmeat Control 94-0176-010

```

Data file: Y1086                      RAR number: 94-0176
Method file: GOALVM                   Sample No: 010
Type: SAMPLE                          Component: NUT
-----
Sample Name: NA                        Cal. Curve: 05/16/95
Date: 16 May 1995 16:23 Method: GOALVM Analyst: SS 43
Interface: 707 Cycles: 6 Channel: A
-----
Instrument: VARIAN 3500.
Column: RTX-200, 0.23mm ID, 1.0um df Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1
Program Rate (deg/min): 5C End Time-Temp (deg-min): 350 C-25
Prog slope (0 or Linear): NA Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min. Split Ratio: NA
Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 750 millivolts ( 70 mv offset)
    
```




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.
15.63	GOAL	0.000E+00	.000E+00	5.00	5.00

### Figure 38. Peanut Nutmeat Fortification 0.01 PPM 93-0142-005

```

Data file: Y042A9                      RAR number: 93-0142
Method file: GOALVn                     Sample No: 005
Type: FORTIFICATION                  Component: NUT
-----
Sample Name: NA                          Cal. Curve: 05/12/95
Date: 12 May 1995 19:29 Method: GOALVn   Analyst: SE LS
Interface: 707                          Cycle#: 9      Channel#: A
-----
Instrument: VARIAN 3500.
Column: RTX-200, 0.32mm ID, 1.0um df   Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C       Ramp Hold (deg-min): 1
Program Rate (deg/min): 5C             End Time-Temp (deg-min): 250 C-25
Prog slope (# or Linear): NA           Inj Port Temp: 265
Flowrate/Gas: 1.5ml/min.              Split Ratio: NA
Det 1-Type & Temp: ECD/300C           Det 2-Type & Temp: NA
-----
Plot time: 0 to 30 minutes
Plot range: 750 millivolts ( 70 mv offset)
    
```




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
5.62	GOAL	8.520E+01	.680E+04	5.00	5.00	0.050

### Figure 39. Peanut Nutmeat Fortification 0.10 PPM 94-0176-010

```

Data file: Y1028                      RAR number: 94-0176
Method file: GOALVA                   Sample No: 010
Type: FORTIFICATION                 Component: NUT
-----
Sample Name: NA                        Cal. Curve: 05/16/95
Date: 14 May 1995 17:37 Method: GOALVA Analyst: SE (J)
Interface: 707                         Cycles: 8      Channel: A
-----
Instrument: VARIAN 3500.
Column: RTX-200, 0.32mm ID, 1.0um dE Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C      Ramp Hold (deg-min): 1
Program Rate (deg/min): 5C            End Time-Temp (deg-min): 250 C-25
Prog Slope (S or Linear): NA          Inj Port Temp: 255
Flowrate/Gas: 3.5ml/min.              Split Ratio: NA
Det 1-Type & Temp: ECD/300C           Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 750 millivolts ( 70 mv offset)
    
```




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
15.63	GOAL	2.790E+02	.224E+05	10.0	5.00	0.500

### Figure 40. Peanut Nutmeat Fortification 0.50 PPM 93-0157-004

```

Data file: Y048A11                      RAR number: 93-0157
Method file: GOALVn                      Sample No: 084
Type: FORTIFICATION                   Component: NUT
-----
Sample Name: NA                           Cal. Curve: 05/12/98
Date: 13 May 1998 20:43 Method: GOALVn   Analyst: SZ LJ
Interface: 707                           Cycles: 11      Channel: A
-----
Instrument: VARIAN 3500.
Column: RTX-200, 0.22mm ID, 1.0um df   Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C       Ramp Hold (deg-min): 1
Program Rate (deg/min): 5C             End Time-Temp (deg-min): 250 C-25
Prog Slope (S or Linear): NA           Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min.              Split Ratio: NA
Det 1-Type & Temp: ECD/300C           Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 750 millivolts ( 70 uv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
15.61	GOAL	1.570E+02	.287E+05	50.0	5.00	2.50


### Figure 41. Pomegranate Control 94-0143-001

Data file: L0788 RAR number: 94-0143  
 Method file: GOALVn Sample No: 0013  
 Type: SAMPLE Component: POMEGRANATE

Sample Name: NA Cal. Curve: 07/07/98  
 Date: 7 Jul 1998 17:31 Method: GOALVn Analyst: SS  $\zeta$  E  
 Interface: 711 Cycles: 3 Channels: A

Instrument: VARIAN 3800,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog slope (0 or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 1.5ml/min Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 10 minutes  
 Plot range: 50 millivolts (-.28 uv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp Wt.
5.00	GOAL	0.000E+00	.000E+00	5.00	5.00

Figure 42. Pomegranate Fortification 0.01 PPM  
94-0143-001

```

Data file: L0788                      RAR number: 94-0143
Method file: GOALVM                   Sample No: 001
Type: FORTIFICATION                 Component: POMEGRANATE
-----
Sample Name: MA                        Cal. Curve: 07/07/98
Date: 7 Jul 1995 18:46                Method: GOALVM      Analyst: BS
Interface: 711                         Cycles: 8 62       Channel: A
-----
Instrument: VARIAN 3500,
Column: RTX-200, 0.32mm ID, 1.0um df  Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C      Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C           End Time-Temp (deg-min): 280 C-17
Prog Slope (S or Linear): MA          Inj Port Temp: 265
Flowrate/Gas: 3.5ml/minR             Split Ratio: MA
Det 1-Type & Temp: ECD/300C           Det 2-Type & Temp: MA
-----
Plot times: 0 to 30 minutes
Plot ramps: 50 millivolts (-.26 mv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ul)	Comp Wt.	ug Added
13.23	GOAL	2.660E+00	.220E+03	5.00	5.00	.0050






### Figure 44. Pomegranate Fortification 0.50 PPM 94-0143-001

```

Data file: 107810                      RAW number: 94-0143
Method file: GOALPa                     Sample No: 001
Type: FORTIFICATION                     Component: POMEGRANATE
-----
Sample Name: NA                          Cal. Curve: 07/07/98
Date: 7 Jul 1995 19:37                   Method: GOALPa      Analyst: SS / J
Interface: 711                            Cycle#: 10          Channel#: A
-----
Instrument: VARIAN 3500.
Column: RTX-200, 0.32mm ID, 1.0um df    Column Length: 50 Meters
Start Temp-Time (deg-min): 115 C        Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C             End Time-Temp (deg-min): 250 C-17
Prog Slope (S or Linear): NA            Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min                 Split Ratio: NA
Det 1-Type & Temp: ECD/300C             Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 50 millivolts (-.27 mv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
13.33	GOAL	1.180E+01	.945E+03	50.0	5.00	2.50


### Figure 45. Soybean Control 92-0107-014

Data files: X209949                      RAR numbers: 92-0107  
 Method files: NQDALP1m                Sample No: 014  
 Type: SAMPLE                        Component: SOYBEAN

Sample Name: SOYBEAN                      Cal. Curves: 11-9-94  
 Date: 9 Nov 1994 12:27                  Method: NQDALP1m                  Operator: DAN  
 Interfaces: 711                              Cycles: 9                              Channel: A

Instrument: VARIAN 3300  
 Column: Rtx-200                              Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215              Ramp Hold (deg-min): 17  
 Program Rate (deg/min): 10                  End Time-Temp (deg-min): 250  
 Prog Slope (# or Linear):                      Inj Port Temp: 265  
 Flowrate/Gas: 3.5                              Split Ratio:  
 Det 1-Type & Temp:                              Det 2-Type & Temp: ECD/270

Plot times: 0 to 22 minutes  
 Plot ranges: 50 millivolts (-.038 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Sample No.
16.10	GOAL	0.000E+00	.000E+00	5.00	5.00


Figure 46. Soybean Fortification 0.01 PPM  
92-0107-014

Data files: X0099410                      RAR number: 92-0107  
 Method files: HQCALP1m                  Sample No: 014  
 Type: FORTIFICATION                      Component: SOYBEAN

Sample Name: SOYBEAN                      Cal. Curves: 11-9-94  
 Date: 9 Nov 1994 12:52                  Methods: HQCALP1m                  Operator: DAN  
 Interfaces: 711                      Cycles: 10                      Channel: A

Instrument: VARIAN 3300                      Column Length: 60 Meters  
 Column: RTX-200                      Ramp Hold (deg-min): 17  
 Start Temp-Time (deg-min): 215                  Program Rate (deg/min): 10                  End Time-Temp (deg-min): 250  
 Prog Slope (S or Linear):                  Inj Port Temp: 265  
 Flowrate/Gas: 3.5                      Split Ratio:  
 Det 1-Type & Temp:                      Det 2-Type & Temp: ECD/Z70

Plot times: 0 to 22 minutes  
 Plot range: 50 millivolts (-.092 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Sample Vol.	ug Added
16.09	GDAL	5.260E+00	.164E+03	5.00	5.00	0.050


Figure 47. Soybean Fortification 0.10 PPM  
92-0107-014

Data files: XH099412                      BAR number: 92-0107  
 Method files: HQDALP1a                  Sample No: 014  
 Types: FORTIFICATION                    Component: SOYBEAN

Sample Name: SOYBEAN                      Cal. Curves: 11-9-94  
 Date: 9 Nov 1994 13:42                  Method: HQDALP1a                  Operator: DAN  
 Interfaces: 711                              Cycles: 12                              Channel: A

Instrument: VARIAN 3300  
 Column: RTX-200                              Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215                  Ramp Hold (deg-min): 17  
 Program Rate (deg/min): 10                  End Time-Temp (deg-min): 250  
 Prog Slope (# or Linear):                      Inj Port Temp: 265  
 Flowrate/Gas: 3.5                              Split Ratio:  
 Det 1-Type & Temp:                              Det 2-Type & Temp: ECD/270

Plot times: 0 to 22 minutes  
 Plot ranges: 50 millivolts ( .001 mv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
16.08	GOAL	5.410E+01	.160E+04	5.00	5.00	0.500

### Figure 48. Soybean Fortification 0.50 PPM 92-0107-014

```

Data file: X2099414          RAR number: 92-0107
Method file: HQCALP1m       Sample No: 014
Type: FORTIFICATION         Component: SOYBEAN
-----
Sample Name: SOYBEAN        Cal. Curve: 11-9-94
Date: 9 Nov 1996 16:32     Method: HQCALP1m      Operator: DAN
Interface: 711              Cycle#: 16              Channel#: A
-----
Instrument: VARIAN 3300
Column: RTX-200            Column Length: 60 Meters
Start Temp-Time (deg-min): 215   Ramp Hold (deg-min): 17
Program Rate (deg/min): 10       End Time-Temp (deg-min): 250
Prog Slope (S or Linear):          Inj Port Temp: 265
Flowrate/Gas: 3.5           Split Ratio:
Det 1-Type & Temp:          Det 2-Type & Temp: ECD/Z70
    
```

Plot times 0 to 22 minutes  
Plot range: 50 millivolts (-.038 mv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Sample Wt.	ug Added
16.09	GOAL	2.890E+01	.897E+03	50.0	5.00	2.50

### Figure 49. Peanut Hay Control 93-0142-003

```

Data file: Y1257                      RAW number: 93-0142
Method file: GOALYN                    Sample No: 003
Type: SAMPLE                           Component: HAY
-----
Sample Name: NA                          Cal. Curve: 05/15/95
Date: 16 May 1995 11:37                 Method: GOALYN      Analyte: GS (3)
Interface: 707                           Cycles: 7          Channel: A
-----
Instrument: VARIAN 3500.
Column: RTX-200, 0.32mm ID, 1.0um dF   Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C        Ramp Hold (deg-min): 1
Program Rate (deg/min): 5C              End Time-Temp (deg-min): 250 C-15
Prog Slope (S or Linear): NA            Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min.               Split Ratio: NA
Det 1-Type & Temp: ECD/300C            Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 750 millivolts ( 70 uv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ul)	Comp. Wt.
15.63	GOAL	1.870E+06	.140E+03	10.0	5.00


**Figure 50. Peanut Hay Fortification 0.02 PPM  
93-0157-002**

Data file: Y1512 RAR number: 93-0157  
 Method file: GOALPa Sample No: 003  
 Type: FORTIFICATION Component: HAY

Sample Name: HA Cal. Curve: 05/16/93  
 Date: 16 May 1995 14:41 Method: GOALPa Analyst: SS  
 Interface: 787 Cycles: 12 Channel: A

Instrument: VARIAN 2800,  
 Column: RTX-200, 0.32mm ID, 1.0um dF Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
 Prog Slope (S or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/min. Split Ratio: NA  
 Det 1-Type & Temp: ECD/200C Det 2-Type & Temp: NA

Plot time: 0 to 30 minutes  
 Plot range: 780 millivolts ( 70 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
18.61	GOAL	7.760E+01	.896E+04	10.0	5.00	0.100







### Figure 53. Peanut Shell Control 93-0157-004


```

Data file: Y1554                      RRR number: 93-0157
Method file: GOALVA                   Sample No: 004
Type: SAMPLE                           Component: SHELL

-----
Sample Name: NA                        Cal. Curve: 05/17/95
Date: 17 May 1995 13:12               Method: GOALVA      Analyst: SE 52
Interface: 707                        Cycles: 0           Channel: A

-----
Instrument: VARIAN 3500.
Column: RTX-100, 0.32mm ID, 1.0um df  Column Length: 60  Retore
Start Temp-Time (deg-min): 215 C      Ramp Hold (deg-min): 1
Program Rate (deg/min): 5C            End Time-Temp (deg-min): 250 C-25
Prog Slope (0 or Linear): NA          Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min.             Split Ratio: NA
Det 1-Type & Temp: SCD/300C           Det 2-Type & Temp: NA

-----
Plot times: 0 to 30 minutes
Plot range: 750 millivolts ( 70 uv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Resp. (uv)
15.64	GOAL	0.000E+00	.000E+00	10.0	1.00








### Figure 57. Peanut Vine Control 93-0142-001

```

Data file: Y11E7                      RAR number: 93-0142
Method file: GOALVn                   Sample No: 001
Type: SAMPLE                           Component: VINE
-----
Sample Name: NA                        Cal. Curve: 05/15/95
Date: 15 May 1995 11:17              Method: GOALVn      Analyst: JS
Interface: 787                        Cycles: 7           Channel: A
-----
Instrument: VARIAN 3500,
Column: RTX-200, 0.32mm ID, 1.0um dF  Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C      Ramp Hold (deg-min): 1
Program Rate (deg/min): 5C            End Time-Temp (deg-min): 250 C-25
Prog Slope (S or Linear): NA          Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min.              Split Ratio: NA
Det 1-Type & Temp: ECD/300C           Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 750 millivolts ( 70 uv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp Wt.
15.62	GOAL	0.000E+00	.000E+00	10.0	1.00

### Figure 58. Peanut Vine Fortification 0.02 PPM 93-0142-001

Data file: Y11214 RAR number: 93-0142  
 Method file: GOALVA Sample No: 001  
 Type: ~~XXXXXXXXXXXX~~ Component: VINE

---


Sample Name: NA Cal. Curve: 05/15/95  
 Date: 15 May 1995 15:33 Method: GOALVA Analyst: JH  
 Interface: 707 Cycled: 14 Channels: A

---

Instrument: VARIAN 3500,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
 Prog Slope (S or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.6ml/min. Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

---

Plot times: 0 to 30 minutes  
 Plot range: 750 millivolts ( 70 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp. Wt.	ug Added
10.03	GOAL	6.550E+01	.507E+04	10.0	1.00	0.100





### Figure 60. Peanut Vine Fortification 0.50 PPM 93-0142-001

Data file: Y11218 RAR number: 93-0142  
 Method file: GOALVn Sample No: 001  
 Type: FORTIFICATION Component: VINE

---


Sample Name: NA Cal. Curve: 05/15/95  
 Date: 15 May 1995 18:00 Method: GOALVn Analyst: JS  
 Interface: 707 Cycle: 18 Channel: A

---

Instrument: VARIAN 3500.  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
 Prog slope (S or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/min. Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

---

Plot times: 0 to 30 minutes  
 Plot range: 750 millivolts ( 70 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp. Wt.	ug Added
18.62	GOAL	2.678E+02	.212E+05	50.0	5.00	2.50


### Figure 61. Grape Control (Rtx-50 confirmatory)

Data file: N153A8 RAR number: 93-0012  
 Method file: GOALS0Vm Sample No: 001 N/A  
 Type: SAMPLE Component: FRUIT

Sample Name: NA Cal. Curve: 11/18/94  
 Date: 19 Nov 1994 03:23 Method: GOALS0Vm Analyst: EZ 4  
 Interface: 707 Cycles: 8 Channel: A

Instrument: VARIAN 3500, GC 14  
 Column: RTX-50, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 5C End Time-Temp (deg-min): 250 C-25  
 Prog Slope (# or Linear): NA Inj Port Temp: 265C  
 Flowrate/Gas: 3.5ML/MIN Split Ratio: NA  
 Det 1-Type & Temp: ECD/270C Det 2-Type & Temp: NA

Plot times: 0 to 40 minutes  
 Plot range: 40 millivolts (-1.2 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Comp Wt.
22.50	GOAL	0.000E+00	.000E+00	5.00	5.00

Figure 62. Grape Fortification 0.01 PPM  
(Rtx-50 confirmatory)

Data file: N155A11 RAR number: 93-0012  
 Method file: GOALS01a Sample No: 001 57K 15  
 Type: FORTIFICATION Component: FRUIT

---


Sample Name: NA Cal. Curves: 11/18/94  
 Date: 19 Nov 1994 05:13 Method: GOALS01a Analyst: SZ42  
 Interface: 707 Cycles: 11 Channel: A

---

Instrument: VARIAN 3500, GC 14  
 Column: RTX-50, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 5C End Time-Temp (deg-min): 230 C-25  
 Prog Slope (S or Linear): NA Inj Port Temp: 265C  
 Flowrate/Gas: 3.5ML/MIN Split Ratio: NA  
 Det 1-Type & Temp: ECD/270C Det 2-Type & Temp: NA

---

Plot time: 0 to 40 minutes  
 Plot range: 40 millivolts (-1.1 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	Wt Added
22.56	GOAL	2.760E+00	.139E+03	5.00	5.00	0.050





### Figure 65. Soybean Control (Rtx-50 confirmatory)

Data files: N0859	RAR number: 92-0107
Method files: GOALS0Vm	Sample No: 014
Type: <u>SAMPLE</u>	Component: SOYSEAN

---


Sample Name: NA	Cal. Curve: 11/21/94	
Date: 21 Nov 1994 16:54	Method: GOALS0Vm	Analyst: SZ SZ
Interface: 707	Cycle#: 9	Channel#: A

---

Instrument: VARIAN 3300, GC 14

Column: RTX-50, 0.32mm ID, 1.0um df.	Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C	Ramp Hold (deg-min): 1
Program Rate (deg/min): 5C	End Time-Temp (deg-min): 250 C-25
Prog Slope (S or Linear): NA	Inj Port Temp: 265C
Flowrate/Gas: 3.5ML/MIN	Split Ratio: NA
Det 1-Type & Temp: ECD/270C	Det 2-Type & Temp: NA

Plot times: 0 to 40 minutes  
 Plot range: 40 millivolts (-1.2 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.
22.78	GOAL	1.750E+00	.955E+02	5.00	5.00



Figure 67. Soybean Fortification 0.10 PPM  
(Rtx-50 confirmatory)

Data file: N08213 RAR number: 92-0107  
 Method file: GOAL50Vn Sample No: 016  
 Type: FORTIFICATION Component: SOYBEAN

---


Sample Name: NA Cal. Curves: 11/21/94  
 Date: 21 Nov 1994 19:21 Method: GOAL50Vn Analyst: SZ 42  
 Interfaces: 707 Cycles: 13 Channel: A

---

Instrument: VARIAN 3300, GC 14  
 Column: RTX-50, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 5C End Time-Temp (deg-min): 230 C-25  
 Prog Slope (\$ or Linear): NA Inj Port Temp: 265 C  
 Flowrate/Gas: 3.5ML/NIN Split Ratio: NA  
 Det 1-Type & Temp: ECD/270C Det 2-Type & Temp: NA

---

Plot times: 0 to 40 minutes  
 Plot range: 40 millivolts (-1.1 mv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Sample Wt.	ug Added
22.60	GOAL	2.580E+01	.121E+06	5.00	5.00	0.500



### Figure 68. Soybean Fortification 0.50 PPM (Rtx-50 confirmatory)

Data files: N08S14		RAR number: 92-0107
Method file: GOALS0VM		Sample No: 014
Type: <u>FORTIFICATION</u>		Component: SOYBEAN
<hr/>		
Sample Name: NA	Cal. Curves: 11/21/94	
Date: 21 Nov 1994 19:58	Method: GOALS0VM	Analyst: SFC
Interface: 707	Cycle#: 14	Channel#: A
<hr/>		
Instrument: VARIAN 3500, GC 14		
Column: RTX-50, 0.32mm ID, 1.0um df	Column Length: 60 Meters	
Start Temp-Time (deg-min): 215 C	Ramp Hold (deg-min): 1	
Program Rate (deg/min): 5C	End Time-Temp (deg-min): 250 C-25	
Prog Slope (S or Linear): NA	Inj Port Temp: 265C	
Flowrate/Gas: 3.5NL/MIN	Split Ratio: NA	
Det 1-Type & Temp: ECD/270C	Det 2-Type & Temp: NA	
<hr/>		
Plot times: 0 to 40 minutes		
Plot range: 40 millivolts (-1.1 mv offset)		



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
22.59	GOAL	1.620E+01	.785E+03	50.0	5.00	2.50







Figure 72. Standard 0.10 µg/ml (GC/ECD)

Data file: L8782 Type: STANDARD

---


Sample Name: NA Cal. Curve: 07/07/98  
 Date: 7 Jul 1998 16:15 Method: GOALV6 Analyze: 02 1-2  
 Interface: 711 Cycles: 1 Channel: A

---

Instrument: VARIAN 3800.  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog slope (S or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 1.5ml/minE Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

---


Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.25 mv offset)



Retention Time	Compound Name	FW Injected	Area	Height
13.34	GOAL	0.100	2.810E+01	2.250E+02



Figure 74. Oxyfluorfen Calibration Curve (GC/ECD)




Concentrations in report are calculated from equation:  
 $HEIGHT = B0 + B1(UG/ML) + B2(UG/ML)^2$   
obtained by least-squares fit of standard injection data.

Figure 75. Grape Control (GC/MS)

File : C:\HPCHEM\1\DATA\N08S9.D  
 Operator : sz  
 Acquired : 9 Nov 95 10:41 pm using AcqMethod GOAL1  
 Sample Name: 92-0100-005  
 Misc Info : blk SW=5g FV=2.5ml  
 Vial Number: 9  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.37  
 Concentration: 12.33  
 Pk # and Type: 1 RTEINT used  
 # Qualifiers Not Satisfied



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.37	19.16	114	sys def
Q1	361.00	0.0	0.0- 22.0	0.00	to	0	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def







Figure 76. Grape Fortification 0.025 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N08S10.D  
 Operator : SZ  
 Acquired : 9 Nov 95 11:19 pm using AcqMethod GOAL1  
 Sample Name: 92-0100-005  
 Misc Info : spka SW=5g FV=2.5ml 0.025ppm  
 Vial Number: 10  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 46.28  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	4397	sys def
Q1	361.00	13.5	0.0- 22.0	19.64	to	592	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 77. Grape Fortification 0.05 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N08S11.D  
 Operator : sz  
 Acquired : 9 Nov 95 11:57 pm using AcqMethod GOAL1  
 Sample Name: 92-0100-005  
 Misc Info : spkb SW=5g FV=2.5ml 0.05ppm  
 Vial Number: 11  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 85.98  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	9405	sys def
Q1	361.00	11.5	0.0- 22.0	19.65	to	1084	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 78. Grape Fortification 0.125PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N08S12.D  
 Operator : sz  
 Acquired : 10 Nov 95 12:35 am using AcqMethod GOAL1  
 Sample Name: 92-0100-005  
 Misc Info : spkc SW=5g FV=2.5ml 0.125ppm  
 Vial Number: 12  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 269.02  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.16	32497	sys def
Q1	361.00	11.4	0.0- 22.0	19.65	to	3689	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 79. Peach Control (GC/MS)

File : C:\HPCHEM\1\DATA\N08S5.D  
 Operator : sz  
 Acquired : 9 Nov 95 8:08 pm using AcqMethod GOAL1  
 Sample Name: 94-0117-001  
 Misc Info : blk SW=5g FV=2.5ml  
 Vial Number: 5  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 20.06  
 Concentration: 12.30  
 Pk # and Type: 1 RTEINT used  
 # Qualifiers Not Satisfied



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		20.06	19.16	110	sys def
Q1	361.00	0.0	0.0- 22.0	0.00	to	0	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 80. Peach Fortification 0.025 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N08S6.D  
 Operator : sz  
 Acquired : 9 Nov 95 8:46 pm using AcqMethod GOAL1  
 Sample Name: 94-0117-001  
 Misc Info : spka SW=5g FV=2.5ml 0.025ppm  
 Vial Number: 6  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 42.26  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	3890	sys def
Q1	361.00	10.5	0.0- 22.0	19.65	to	410	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 81. Peach Fortification 0.05 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N08S7.D  
 Operator : SZ  
 Acquired : 9 Nov 95 9:24 pm using AcqMethod GOAL1  
 Sample Name: 94-0117-001  
 Misc Info : spkb SW=5g FV=2.5ml 0.05ppm  
 Vial Number: 7  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 82.39  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	8952	sys def
Q1	361.00	12.4	0.0- 22.0	19.65	to	1110	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 82. Peach Fortification 0.125 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N08S8.D  
 Operator : sz  
 Acquired : 9 Nov 95 10:03 pm using AcqMethod GOAL1  
 Sample Name: 94-0117-001  
 Misc Info : spkc SW=5g FV=2.5ml 0.125ppm  
 Vial Number: 8  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 255.18  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	30752	sys def
Q1	361.00	11.3	0.0- 22.0	19.65	to	3485	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 83. Peanut Nutmeat Control (GC/MS)

File : C:\HPCHEM\1\DATA\N09S3.D  
 Operator : SZ  
 Acquired : 10 Nov 95 4:39 pm using AcqMethod GOAL1  
 Sample Name: 93-0157-004 nutmeat  
 Misc Info : blk SW=5g, FV=2.5ml  
 Vial Number: 3  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.60  
 Concentration: 27.22  
 Pk # and Type: 1 RTEINT used  
 # Qualifiers Not Satisfied



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.60	19.16	197	sys def
Q1	361.00	0.0	0.0- 22.0	0.00	to	0	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def







Figure 84. Peanut Nutmeat Fortification 0.025 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N09S7.D  
 Operator : sz  
 Acquired : 10 Nov 95 7:12 pm using AcqMethod GOAL1  
 Sample Name: 93-0157-004 nutmeat  
 Misc Info : spka SW=5g FV=2.5ml 0.025ppm  
 Vial Number: 7  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 54.38  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	7784	sys def
Q1	361.00	11.5	0.0- 22.0	19.64	to	899	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 85. Peanut Nutmeat Fortification 0.05 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N09S11.D  
 Operator : SZ  
 Acquired : 10 Nov 95 9:44 pm using AcqMethod GOAL1  
 Sample Name: 93-0157-004 nutmeat  
 Misc Info : spkb SW=5 FV=2.5 0.05ppm  
 Vial Number: 11  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 80.70  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	15137	sys def
Q1	361.00	11.2	0.0- 22.0	19.65	to	1700	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 86. Peanut Nutmeat Fortification 0.125 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N09S13.D  
 Operator : sz  
 Acquired : 10 Nov 95 11:01 pm using AcqMethod GOAL1  
 Sample Name: 93-0157-004 nutmeat  
 Misc Info : spk SW=5 FV=2.5 0.125ppm  
 Vial Number: 13  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 217.28  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.16	53290	sys def
Q1	361.00	11.5	0.0- 22.0	19.65	to	6113	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 87. Peanut Shell Control (GC/MS)

File : C:\HPCHEM\1\DATA\N09S4.D  
 Operator : sz  
 Acquired : 10 Nov 95 5:17 pm using AcqMethod GOAL1  
 Sample Name: 93-0157-004 shell  
 Misc Info : blk SW=5g, FV=2.5ml  
 Vial Number: 4  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.79  
 Concentration: 28.01  
 Pk # and Type: 1 RTEINT used  
 # Qualifiers Not Satisfied



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.79	19.16	418	sys def
Q1	361.00	0.0	0.0- 22.0	0.00	to	0	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 88. Peanut Shell Fortification 0.025 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N09S8.D  
 Operator : SZ  
 Acquired : 10 Nov 95 7:50 pm using AcqMethod GOAL1  
 Sample Name: 93-0157-004 shell  
 Misc Info : spka SW=5g FV=2.5ml 0.025ppm  
 Vial Number: 8  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 52.23  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	7183	sys def
Q1	361.00	11.3	0.0- 22.0	19.65	to	814	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 89. Peanut Shell Fortification 0.05 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N09S12.D  
 Operator : sz  
 Acquired : 10 Nov 95 10:23 pm using AcqMethod GOAL1  
 Sample Name: 93-0157-004 shell  
 Misc Info : spkb SW=5 FV=2.5 0.05ppm  
 Vial Number: 12  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 117.28  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	25354	sys def
Q1	361.00	11.3	0.0- 22.0	19.65	to	2876	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 90. Peanut Shell Fortification 0.125 PPM (GC/MS)

File : C:\HPCHEM\1\DATA\N09S14.D  
 Operator : SZ  
 Acquired : 10 Nov 95 11:39 pm using AcqMethod GOAL1  
 Sample Name: 93-0157-004 she#11 *RE RB 11/15/95*  
 Misc Info : spk# SW=5 FV=2.5 0.125ppm  
 Vial Number: 14 *RE RB 11/15/95*  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 256.92  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.16	64363	sys def
Q1	361.00	11.5	0.0- 22.0	19.66	to	7428	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 91. Standard 0.025 µg/ml (GC/MS)

File : C:\HPCHEM\1\DATA\N08S4.D  
 Operator : sz  
 Acquired : 9 Nov 95 7:30 pm using AcqMethod GOAL1  
 Sample Name: std 0.025ug/ml  
 Misc Info : C103195-5  
 Vial Number: 4  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 30.45  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.16	2400	sys def
Q1	361.00	12.6	0.0- 22.0	19.65	to	303	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def







Figure 92. Standard 0.05 µg/ml (GC/MS)

File : C:\HPCHEM\1\DATA\N08S3.D  
 Operator : sz  
 Acquired : 9 Nov 95 6:51 pm using AcqMethod GOAL1  
 Sample Name: std 0.05ug/ml  
 Misc Info : C103195-4  
 Vial Number: 3  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 79.04  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.22	5254	sys def
Q1	361.00	12.6	0.0- 22.0	19.65	to	662	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.21	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 93. Standard 0.10 µg/ml (GC/MS)

File : C:\HPCHEM\1\DATA\N08S2.D  
 Operator : sz  
 Acquired : 9 Nov 95 6:13 pm using AcqMethod GOAL1  
 Sample Name: std 0.1ug/ml  
 Misc Info : C103195-3  
 Vial Number: 2  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 119.45  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.22	11274	sys def
Q1	361.00	11.5	0.0- 22.0	19.65	to	1294	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.21	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 94. Standard 0.50 µg/ml (GC/MS)

File : C:\HPCHEM\1\DATA\N08S1.D  
 Operator : sz  
 Acquired : 9 Nov 95 5:35 pm using AcqMethod GOAL1  
 Sample Name: std 0.5ug/ml  
 Misc Info : C103195-2  
 Vial Number: 1  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 682.50  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.22	95162	sys def
Q1	361.00	11.1	0.0- 22.0	19.65	to	10593	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.21	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def





Figure 95. Oxyfluorfen Calibration Curve (GC/MS)



Method Name: C:\HPCHEM\1\METHODS\GOAL1.M  
Calibration Table Last Updated: Tue Nov 14 09:59:48 1995

## **Appendix 1**

**GC/MS Method Validation Protocol  
Certificate of Analysis (RH-32915)  
GLP Statements  
Quality Assurance Statements**

I. Protocol Number: 34P-95-92

Project Number: 3107.14

II. Study Title

Method Validation- GC/MS Confirmatory Detection in Oxyfluorfen (Goal®) Crop Residue Analytical Method, TR 34-94-150

III. Study Objective

The objective is to validate the Oxyfluorfen (GOAL®) Crop Residue Analytical Method, TR 34-94-150, using GC/MS as the confirmatory detection method. The crops that will be used are peanut nutmeat, peanut shell, peach, and grape\*. The data will be used in the reissue of this residue method as the Tolerance Enforcement Method, TR 34-95-111.

\*These samples will be taken from the residue analysis covered in subprotocol 34P-94-65A.

IV. Sponsor

Rohm and Haas Company  
727 Norristown Road  
Spring House, PA 19477-0904

V. Testing Laboratory

Centre Analytical Laboratories  
3048 Research Drive  
State College, PA 16801

VI. Personnel

5.1 Study Director:\* Dennis A. Martin  
Rohm and Haas Company  
215-619-5576

5.3 Principal Investigator: Shaozhi Zheng  
Centre Analytical Laboratories

- 5.2 Contract Laboratory Coordinator:\*\* Roberta C. Regetta  
Rohm and Haas Company  
215-619-5588
- 5.4 RMEF Manager: Dr. Stanley S. Stavinski  
Rohm and Haas Company

\*Primary contact for technical issues

\*\*Primary contact for non-scientific issues

#### VII. Dates

- Proposed Analytical Start Date: October 1995
- Proposed Experimental Completion: November 1995
- Proposed Final Report Date: November 1995

#### VIII. Identification and Justification For Selection of Test System

The test system for this study is defined as sets of untreated crop samples collected from field trials. Samples will be identified by the Rohm and Haas sample numbers defined later in this protocol. The samples are selected to support magnitude of residue studies. The crops specifically to be used in this study are peanut nutmeat, grape, peanut shell, and peach.

#### IX. Experimental Design and Methods to Control Bias

Untreated crop samples will be fortified with oxyfluorfen and run through the method to validate the GC/MS method as the confirmatory detection technique. Sample analysis will follow documented standard operating procedures. Methods to control bias include analysis of untreated control samples, fortification recovery data, and generation of standard curves.

#### X. Method

All analyses and preparation of stock fortification solutions will be performed using the residue analytical method TR 34-94-150, "Oxyfluorfen (GOAL®) Crop Residue Analytical Method." However, there is the addition of GC/MS detection, developed by Centre

Analytical Laboratories, as the mode of quantitation in place of GC/ECD. (Data acceptability - SOP RMEF.RES.05.02)

**XI. Test and Reference Substances**

1. Test and Reference Material

<u>Standard</u>	<u>Lot</u>	<u>Purity/%</u>	<u>Appearance</u>	<u>Expiration Date</u>
RH-32915 * oxyfluorfen	RPO8674FP	99.4	Pale yellow powder	August 2, 2000

Analytical grade compounds will be used as both a test substance, for fortifications, and a reference substance, to generate calibration curves. The reference material has been GLP analyzed and possesses a certificate of analysis. The characterization of the reference substance is detailed in Rohm and Haas report, APR-95-334.

\* RH-32915 and RH-2915 identify the same standard. The former appears on the certificate of analysis, the latter is used in the residue database.

2. Chemicals- The grade and supplier of chemicals are defined in the method.
3. Hazard Information- Current MSD sheets for the standards will be available by the testing laboratory.

**XII. Treatment of Test System**

The samples used are untreated control samples from field studies. The RMEF RAR code system will be used to identify all samples. The samples have been processed according to SOP SAMP.02.06; and, they are maintained in frozen storage until analysis.

The following samples will be analyzed. Each set will consist of one control and three fortifications. The fortifications will be 1 at the LOQ, 1 at 2 X LOQ, and 1 at 5 X LOQ. The grape samples are taken from subprotocol 34P-94-65A.

<u>RAR Code</u>	<u>Sample No.</u>	<u>Variety</u>	<u>Component</u>
93-0157	004	peanut	nutmeat
93-0157	004	peanut	shell
94-0117	001	peach	fruit
92-0100	005	Cabernet grape	fruit



**XIII. Records**

Originals of the items listed below, and any other raw data generated will be provided to the Study Director.

**1. Chromatograms**

All chromatograms will contain the following:

- A. Sample identification, date, final sample volume, sample weight, arrows indicating the area of interest, the analyst's initials, and injection number corresponding to the run sheet.
- B. Fortifications will additionally include the  $\mu\text{g}$  added and the sample number of the sample that was fortified.
- C. Standard chromatograms will additionally include the  $\mu\text{g}/\text{ml}$  concentration.
- D. The first standard chromatogram in a daily run should also include temperatures, flow rates, column parameters, gases, instrument parameters, and instrument type if any of these differ from the method.

- 2. Storage History and Chain of Custody Records
- 3. Method Run Sheets
- 4. Diskettes of Results
- 5. Analyst's Bench Sheet

**XIV. Reporting****1. Data Reporting**

The testing laboratory will provide the original version of run sheets, standard curves, and any other raw data generated. Rohm and Haas will provide the original processing records. Work performed will be certified by the signature of all personnel involved in the analysis. Signed GLP Compliance and Quality Assurance Statements should also be provided.

## 2. Final Report

The analytical data generated under this protocol will be used in the reissue of the Oxyfluorfen Crop Residue Analytical Method, TR 34-94-150, as the Oxyfluorfen Crop Tolerance Enforcement Method, TR 34-95-111. The final method report will be written by the Study Director upon completion of this analysis. The report will include data generated under other protocols covering residue analysis of samples used in the reregistration of oxyfluorfen. The report will include a description of the materials and methods used, the residue analytical data generated, along with signed and dated Study Compliance and Quality Assurance Statements.

## XV. GLP Statement

This analysis will be conducted in compliance with FIFRA's Good Laboratory Practice Standard (40 CFR 160). All raw data generated in this study will be archived by the Rohm and Haas Company. The Study Director will be responsible for the overall conduct of the study and the final report. The Principal Investigator is the Study Director's delegate responsible for this protocol, oversight of the analysis, the analytical experimentation, and will review the protocol, procedures, Good Laboratory Practice and safety guidelines with all support personnel prior to their participation in this study. Any deviations will be recorded. Changes in plans will be approved by the Study Director and attached to this protocol.

## XVI. Quality Assurance

The Centre Analytical Laboratories, Inc's QA Unit will conduct in-progress inspections, audit the raw data to ensure the integrity of the study, and will report findings, noting any problem found and corrective actions taken, to site management, the study director, Rohm and Haas QAU, and Rohm and Haas Management.

**XVII. Signatures**

Study Director:

Dennis A. Martin 10/19/95  
Dennis A. Martin Date

RMEF Manager:

Stan Stavinski Oct 19/95  
Dr. Stanley S. Stavinski Date

Principal Investigator:

Shaozhi Zheng 10/24/95  
Shaozhi Zheng Date

Reviewed by QAU

Rhonda S. Krick 10/19/95  
Rhonda S. Krick Date

**PROTOCOL MODIFICATION FORM**

**SUBPROTOCOL/**  
**PROTOCOL No.** 34P-95-92      **STUDY DIRECTOR** Dennis A. Martin

**STUDY TITLE** Method Validation - GC/MS Confirmatory Detection

**MODIFICATION:** in Oxyfluorfen (GowE) Crop Residue Analytical Method (TR 34-94-150)

① Analyze crop samples generated by GC/ECD as well as GC/MS, using method found in TR 34-94-150.

**JUSTIFICATION:**

② Results are needed on both detection methods to use for comparison.

**EFFECT ON STUDY:**

None.

**Originator:**

Dennis A. Martin      11-14-95  
 signature      date

**Study Director Approval:**

Dennis A. Martin      11-14-95  
 signature      date

## PROTOCOL MODIFICATION FORM

SUBPROTOCOL/

PROTOCOL No. 34-95-92STUDY DIRECTOR D.A. MartinSTUDY TITLE Method Validation- GC/MS Confirmatory Detection inoxyfluorfen (60a1E) crop Residue Analytical  
Method, TR 34-94-150**MODIFICATION:**

The LOQ for the confirmatory method using GC/MS is 0.025 ppm, instead of 0.01 ppm as the primary method dictates.

**JUSTIFICATION:**

The sensitivity of the confirmatory detection method does not need to be the same as the primary method. The LOQ of 0.025 ppm results in the best chromatography for this method.

**EFFECT ON STUDY:**

There will be no effect on the study.

Originator:

Jenniss A. Martin 12/18/95  
signature date

Study Director Approval:

Jenniss A. Martin 12/18/95  
signature date



**PROTOCOL/METHOD MODIFICATION FORM**

Subprotocol/Protocol No: 34P-95-92Study Director: Dennis MartinCAL Project Number: 002-165 Sponsor Study Number: 3107.14\*  
34P-95-92Study Title: method validation - GC/MS confirmatory detection in  
oxyfluorfen (Goal) Crop Residue Analytical method, TR34-94-150.MODIFICATION:

1. In step 6.6, 40 ml of 50:50 (V/V) ethyl acetate/hexane was used instead of 25 ml of 5:95 (V/V) ethyl acetate/hexane.
2. The final volume was changed from 5 ml to 2.5 ml.

JUSTIFICATION:

1. ~~2 modifications were issued.~~ \* 11/14/95
1. ensure the good recovery
2. ensure the fortification \* ~~low~~ sample within the standard range.

EFFECT ON STUDY:

1. none
2. none.

Originator: Shaoshi Zheng11/10/95  
DateStudy Director: Dennis A. Martin1/10/96  
Date



3048 Research Drive  
Phone: (814) 231-8032

State College, PA 16801  
Facsimile: (814) 231-1253

Page 1 of 1

PROTOCOL/METHOD MODIFICATION FORM

Subprotocol/Protocol No: 34P-95-92  
Study Director: Dennis A. Martin

CAL Project Number: 002-165 Sponsor Study Number: 3107.14

Study Title: method validation - GC/MS confirmatory detection in  
Oxyfluorfen (Goal) Crop Residue Analytical method, TR 34-94-150

MODIFICATION:

1. 1000ug/ml stock solution was made instead of 200ug/ml.

JUSTIFICATION:

1. It is laboratory practice to prepare stock solution concentration of 1000ug/ml in order to use the solution for different projects.

EFFECT ON STUDY:

1. none.

Originator: Shaozhi Zheng

1/2/96  
Date

Study Director: Dennis A. Martin

1/10/96  
Date



**PROTOCOL/SOP/METHOD DEVIATION**

Deviation Number: 1  
Date of Occurrence: 10/31/95

CAL Project Number: 002-165 Sponsor Study Number: 3107.14

**DESCRIPTION OF DEVIATION**

Is this deviation a  PROTOCOL, or  SOP, or  METHOD deviation?

standard (RH-32915) purity is 99.7% and expiration date is 7/12/2000.

**ACTIONS TAKEN**

i.e., amendment issued, SOP revision, etc...

protocol deviation was issued.

Recorded By/Date: Shawzhi Zheng 1/3/96

**IMPACT ON THE STUDY**

None.

Dennis A. Martens 1/10/96  
Study Director Signature Date

CAL QAU Review GW 1/5/96



TR 34 95 111

INDEPENDENCE MALL WEST PHILADELPHIA, PA. 19105, U.S.A. TELEPHONE (215) 592-3000  
CABLE ADDRESS: ROHMHAAS TELEX 845-247 TWX 710-670-5335 TELECOPIER (215) 592-3377

REPLY TO:  
RESEARCH LABORATORIES  
727 NORRISTOWN ROAD  
P.O. BOX 904  
TRING HOUSE, PA 19477-0904  
(215) 641-7000  
(215) CH 2-0400

000152



## CERTIFICATE OF ANALYSIS ANALYTICAL STANDARD

**RH NUMBER:** 32915  
**LOT NUMBER:** RPO8674FP  
**MATERIAL NAME:** GOAL  
**PRODUCTLINE:** GOAL  
**PURITY:** 99.4 %  
**GLP ANALYSIS:** YES  
**STORAGE CONDITIONS:** ROOM TEMP. OR BELOW  
**APPEARANCE:** PALE YELLOW POWDER  
**EXPIRATION DATE:** 02-aug-2000  
**DATE OF LAST ANALYSIS:** 02-aug-1995

QUERIED BY Denise Milligan DATE 28-sep-1995  
(Signature)

Supporting documentation for this information is retained at the  
Rohm and Haas Company in the AgChem Analytical Standard Archive.

## GLP COMPLIANCE STATEMENT

CAL Study Number 002-165 "Method Validation - GC/MS Confirmatory Detection in Oxyfluorfen (Goal®) Crop Residue Analytical Method, TR 34-94-150" was performed in compliance with EPA Good Laboratory Practice Standards (40 CFR, Part 160; FR. 8/17/89) by Centre Analytical Laboratories, Inc.

Shaozhi Zheng  
Shaozhi Zheng  
Principal Investigator

4/4/96  
Date

**QUALITY ASSURANCE STATEMENT**

CAL Study Number 002-165 was reviewed by Centre Analytical Laboratories' Quality Assurance Unit. All reviewed phases were reviewed for conduct according to Centre Analytical Laboratories' Standard Operating Procedures, the Study Protocol and all applicable Good Laboratory Practice Standards. All findings were reported to the study director and management.

<u>Phase</u>	<u>Date Inspected</u>	<u>Date Reported to Principal Investigator and CAL management</u>	<u>Date Reported to Study Director and Sponsor management</u>
1. Protocol Review	11/01/95	11/30/95	12/18/95
2. Extraction	11/02/95	11/30/95	12/18/95
3. Extraction	11/09/95	11/30/95	12/18/95
4. Raw Data Review	12/06/95	01/10/96	02/02/96

Gaya Wickremesinhe  
Gaya Wickremesinhe  
Quality Assurance Officer

4/4/96  
Date

The GLP and QA statements that follow are for additional fortifications submitted after completion of the reports, as listed below:

<u>TR #</u>	<u>MRID #</u>	<u>Protocol #</u>	<u>Crop</u>
34-95-116	43794003	34P-95-30A	fig
34-95-119	43794006	34P-95-33A	olive
34-95-117	43794004	34P-95-31A	pomegranate

**GLP COMPLIANCE STATEMENT**

CAL Study Number 002-146 "GOAL® 1.6E Field Residue Trials on Figs" was performed in compliance with EPA Good Laboratory Practice Standards (40 CFR, Part 160; FR. 8/17/89) by Centre Analytical Laboratories, Inc.

Shaozhi Zheng

Shaozhi Zheng  
Principal Investigator

3/19/96

Date

**QUALITY ASSURANCE STATEMENT**

CAL Study Number 002-146 was reviewed by Centre Analytical Laboratories' Quality Assurance Unit. All reviewed phases were reviewed for conduct according to Centre Analytical Laboratories' Standard Operating Procedures, Study Protocol and all applicable Good Laboratory Practice Standards. All findings were reported to the study director and management.

<u>Phase</u>	<u>Date Inspected</u>	<u>Date Reported to Principal Investigator and CAL management</u>	<u>Date Reported to Study Director and Sponsor management</u>
1. Protocol Review	07/10/95	07/20/95	08/18/95
2. Extraction	07/11/95	08/16/95	08/18/95
3. Raw Data Review	08/07-08/95	08/16/95	08/18/95
4. Raw Data Review	08/18/95	09/25/95	12/18/95

Gail L. Keller

Gail L. Keller  
Quality Assurance Officer

3/19/96

Date

**GLP COMPLIANCE STATEMENT**

CAL Study Number 002-151 "GOAL® 1.6E Field Residue Trials on Olives" was performed in compliance with EPA Good Laboratory Practice Standards (40 CFR, Part 160; FR. 8/17/89) by Centre Analytical Laboratories, Inc.

Shaozhi Zheng

Shaozhi Zheng  
Principal Investigator

3/19/96

Date

### QUALITY ASSURANCE STATEMENT

CAL Study Number 002-151 was reviewed by Centre Analytical Laboratories' Quality Assurance Unit. All reviewed phases were reviewed for conduct according to Centre Analytical Laboratories' Standard Operating Procedures, Study Protocol and all applicable Good Laboratory Practice Standards. All findings were reported to the study director and management.

<u>Phase</u>	<u>Date Inspected</u>	<u>Date Reported to Principal Investigator and CAL management</u>	<u>Date Reported to Study Director and Sponsor management</u>
1. Protocol Review	07/07/95	07/20/95	08/18/95
2. Extraction	07/07/95	08/16/95	08/18/95
3. Raw Data Review	08/08/95	08/16/95	08/18/95
3. Raw Data Review	08/18/95	09/25/95	12/18/95

Gail L. Keller

Gail L. Keller  
Quality Assurance Officer

3/19/96

Date



**GLP COMPLIANCE STATEMENT**

CAL Study Number 002-147 "GOAL® 1.6E Field Residue Trials on Pomegranates" was performed in compliance with EPA Good Laboratory Practice Standards (40 CFR, Part 160; FR. 8/17/89) by Centre Analytical Laboratories, Inc.

shaozhi zheng

Shaozhi Zheng  
Principal Investigator

3/19/96

Date

## QUALITY ASSURANCE STATEMENT

CAL Study Number 002-147 was reviewed by Centre Analytical Laboratories' Quality Assurance Unit. All reviewed phases were reviewed for conduct according to Centre Analytical Laboratories' Standard Operating Procedures, Study Protocol and all applicable Good Laboratory Practice Standards. All findings were reported to the study director and management.

<u>Phase</u>	<u>Date Inspected</u>	<u>Date Reported to Principal Investigator and CAL management</u>	<u>Date Reported to Study Director and Sponsor management</u>
1. Protocol Review	07/10/95	07/20/95	08/18/95
2. Extraction	07/07/95	08/16/95	08/18/95
3. Raw Data Review	08/08/95	08/16/95	08/18/95
3. Raw Data Review	08/18/95	09/25/95	12/18/95

Gail L. Keller

Gail L. Keller  
Quality Assurance Officer

3/19/96

Date

The GLP and QA statements that follow are for peach fortifications which will be reported, as listed below:

<u>TR #</u>	<u>MRID #</u>	<u>Protocol #</u>	<u>Crop</u>
34-95-114	not yet assigned as of study completion	34P-95-35A 34P-95-51A	peach

GLP Compliance Statement

We, the undersigned, hereby certify that the data contained in this raw data package for Rohm & Haas Sub Protocol Number 34P-95-35A, "GOAL® 1.6E Field Residue Trials on Peaches," was generated in compliance with Good Laboratory Practice Standards (40 CFR Part 160) applicable to analytical testing facilities.

Kathleen Koltay<sup>2</sup> for Melinda Lalko MLC  
Melinda Lalko, Residue Lab Supervisor  
Principal Investigator

04AUG 95  
[Date]

Diana Price  
Diana Price, Quality Assurance Supervisor

04 Aug 95  
[Date]

Kathleen Lacey  
Kathleen Lacey, Laboratory Director

04AUG 95  
[Date]

Quality Assurance Statement

This data package for Rohm & Haas Sub Protocol Number, 34P-95-35A, "GOAL@ 1.6E Field Residue Trials on Peaches," has been reviewed by the Quality Assurance Unit of McKenzie Laboratories, Inc. It has been found to accurately describe and/or identify the methods and protocol followed in the conduct of the study. In addition, the Quality Assurance Unit conducted the following inspections of this study and has submitted written reports of these inspections to the Study Director, their management, and to the management of McKenzie Laboratories, Inc.

## Date reported to:

<u>Date Inspected</u>	<u>Phase</u>	<u>Study Director and Management</u>	<u>McKenzie Management</u>
06 Jul 95	<ul style="list-style-type: none"> <li>• Fortification</li> <li>• Extraction</li> <li>• Partitioning</li> <li>• Concentration</li> <li>• Protocol/Amendments</li> <li>• Test System Receipt Documentation</li> <li>• Analytical Reference Substance   Receipt Documentation</li> <li>• Chain of Custody</li> <li>• Study Notebook</li> <li>• Training Files</li> <li>• Standard Preparation Documentation</li> </ul>	17 Jul 95*	12 Jul 95
07 Jul 95	<ul style="list-style-type: none"> <li>• Column Clean-up</li> <li>• Centrifuge</li> <li>• Equipment/Instrument Use Logs</li> <li>• Data Packages</li> </ul>	17 Jul 95*	12 Jul 95
04 Aug 95	<ul style="list-style-type: none"> <li>• Data Package Review</li> <li>• Review of Study Package</li> </ul>	04 Aug 95*	04 Aug 95

\* Date inspection report was sent to the Study Director.

  
 Diana L Price, Quality Assurance Supervisor

04 Aug 95  
 [Date]

TR 34 95 111

000165

**GOAL® 1.6E Field Residue Trials on Peaches  
Protocol Number: 34P-95-51A****GLP Compliance Statement**

We, the undersigned, hereby certify that the data contained in this report was generated in compliance with Good Laboratory Practice Standards (40 CFR Part 160) applicable to analytical testing facilities. All original raw data (or certified copies thereof) have been provided to Rohm and Haas Company and copies (or facility-related originals) have been retained at McKenzie Laboratories, Inc.

   
\_\_\_\_\_  
Diana Price, Quality Assurance Supervisor

15 Dec 95  
[Date]

   
\_\_\_\_\_  
Melinda Lalko, Principal Investigator  
Residue Laboratory Supervisor

15 Dec 95  
[Date]

   
\_\_\_\_\_  
Kathryn Koltavy  
Residue Laboratory Manager

15 Dec 95  
[Date]

   
\_\_\_\_\_  
Kathleen Lacey, Laboratory Director

15 DEC 95  
[Date]

Quality Assurance Statement

This data package for Sub Protocol Number: 34P-95-51A (Goal @ 1.6E Field Residue Trials on Peaches), has been reviewed by the Quality Assurance Unit of McKenzie Laboratories, Inc. It has been found to accurately describe and/or identify the methods and protocol followed in the conduct of the study. In addition, the Quality Assurance Unit conducted the following inspections of this study and has submitted written reports of these inspections to the Study Director, their management, and to the management of McKenzie Laboratories, Inc.

<u>Date Inspected</u>	<u>Phase</u>	<u>Date reported to:</u>	
		<u>Study Director and Management</u>	<u>McKenzie Management</u>
21 Nov 95	<ul style="list-style-type: none"> <li>• Fortification</li> <li>• Extraction</li> <li>• Partitioning</li> <li>• Protocol/Amendments</li> <li>• Study Notebook</li> </ul>	*11 Dec 95	07 Dec 95
22 Nov 95	<ul style="list-style-type: none"> <li>• Training Files</li> </ul>	*11 Dec 95	07 Dec 95
28 Nov 95	<ul style="list-style-type: none"> <li>• Column Clean-up</li> <li>• Chain-of-Custody</li> <li>• Gas Chromatography</li> <li>• Test System Receipt Documentation</li> </ul>	*11 Dec 95	07 Dec 95
29 Nov 95	<ul style="list-style-type: none"> <li>• Maintenance Logs</li> <li>• Equipment/Instrument Use Logs</li> </ul>	*11 Dec 95	07 Dec 95

\* Date inspection report was sent to the Study Director/ Study Director Management.

  
 Diana L Price, Quality Assurance Supervisor

14 Dec 95  
 [Date]

## **Appendix II**

**Representative Analytical Run (GC/ECD)**  
**Representative Analytical Run (GC/MS)**



GC RUN SHEET

Date 7/28/95  
 Analyst Julie Burton  
 Instrument Varian 3500 G C14  
 Analytical Method TR34-94-150  
 Notebook Ref. 002-146  
 Crop Fla

Compound 70AL  
 Data Reduction Method 70ALV  
 Raw Data File L263A  
 Date Extracted 7/25-26/95  
 Date Injected 7/28/95, 7/29/95

Inj. No	STD (µg/mL)	RAR No.	Spl. No	Sample Component	Sample Weight (g)	Final Volume (mL)	Fortification Level (µg)	Notes
1	0.005							
2	0.01							
3	0.1							
4	0.1							
5		94-0142	001	Fla	5.0	5.0		
6		↓	↓	↓	↓	↓	0.05	Spk A.
7		↓	↓	↓	↓	↓	↓	Spk A.
8		↓	↓	↓	↓	↓	↓	Spk A.
9		↓	↓	↓	↓	↓	↓	Spk A.
10	0.01							
11	0.05							
12		94-0142	001	Fla	5.0	5.0	0.1	Spk B.
13		↓	↓	↓	↓	↓	↓	Spk B.
14		↓	↓	↓	↓	↓	0.25	
15		↓	↓	↓	↓	↓	0.5	
16		↓	↓	↓	↓	100.0	2.5	
17	0.005							
18	0.05							
19	0.15							
20	0.15							
21								
22								
23								
24								
25								
26								

NOTES: Samples stored in freezer 10A prior to analysis. 7/28/95

Signature Julie Burton

Data file: L26JA1 Type: STANDARD

---

Sample Name: NA Cal. Curve: 07/28/95  
 Date: 28 Jul 1995 16:22 Method: GOALVA Analyst: JS  
 Interface: 707 Cycles: 1 Channel: A


---

Instrument: VARIAN 3500.  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog Slope (0 or Linear): NA Inj Port Temp: 255°C  
 Flowrate/Gas: 3.5ml/minE Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

---

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.13 uv offset)

Intake JS  
 Date 8/7/95  
 Run# 1 to 20



Retention Time	Compound Name	PPM Injected	Area	Height
13.36	GOAL	0.005	1.460E+00	1.270E+02

Data file: L367A1      Type: STANDARD

---

Sample Name: MA      Cal. Curve: 07/28/95

Date: 28 Jul 1995 16:47    Method: GOALVA    Analyst: JB

Interface: 707      Cycles: 2      Channel: A

---

Instrument: VARIAN 3500,

Column: RTX-200, 0.32mm ID, 1.0um df    Column Length: 60 Meters

Start Temp-Time (deg-min): 215 C    Ramp Hold (deg-min): 1

Program Rate (deg/min): 10C    End Time-Temp (deg-min): 250 C-17

Prog Slope (S or Linear): NA    Inj Port Temp: 265


Flowrate/Gas: 3.5ml/minute    Split Ratio: NA

Det 1-Type & Temp: ECD/300C      Det 2-Type & Temp: NA

---

Plot times: 0 to 30 minutes

Plot range: 50 millivolts (-.12 mv offset)



Retention Time	Compound Name	FWW Injected	Area	Height
12.37	GOAL	0.010	1.220E+00	2.660E+02

TR 34 95 111

Data file: L267A3      Type: STANDARD

-----

Sample Name: NA      Cal. Curve: 07/28/95

Date: 28 Jul 1995 17:12    Method: GOALVA    Analyst: JS

Interface: 707      Cycle#: 3      Channel#: A

-----

Instrument: VARIAN 3500,

Column: RTX-200, 0.32mm ID, 1.0um df    Column Length: 60 Meters

Start Temp-Time (deg-min): 215 C      Ramp Hold (deg-min): 1

Program Rate (deg/min): 10C    End Time-Temp (deg-min): 250 C-17

Prog Slope (0 or Linear): NA      Inj Port Temp: 265


Flowrate/Gas: 3.5ml/minE      Split Ratio: NA

Det 1-Type & Temp: ECD/300C      Det 2-Type & Temp: NA

-----

Plot times: 0 to 30 minutes

Plot range: 50 millivolts (-.13 mv offset)




Retention Time	Compound Name	PPM Injected	Area	Height
13.38	GOAL	0.100	3.090E+01	2.590E+03

Data file: L267M Type: STANDARD

Sample Name: NA Cal. Curve: 07/28/95  
 Date: 28 Jul 1995 17:37 Method: GOALVa Analyst: JS  
 Interface: 707 Cycles: 4 Channel: A


Instrument: VARIAN 3500,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog Slope (\$ or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/min Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.15 uv offset)



Retention Time	Compound Name	PW Injected	Area	Height
13.37	GOAL	0.100	3.060E+01	2.590E+03

Data file: L26JA10 Type: STANDARD  
 -----  
 Sample Name: MA Cal. Curve: 07/28/95  
 Date: 28 Jul 1995 20:08 Method: GOALVM Analyst: JB  
 Interface: 707 Cycles: 10 Channel: A  
 -----  
 Instrument: VARIAN 3500,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog Slope (# or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/min Split Ratio: NA  
 Det 1-Type & Temp: BCD/300C Det 2-Type & Temp: NA  
 -----  
 Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.094 mv offset)




Retention Time	Compound Name	FW Injected	Area	Height
13.37	GOAL	0.010	3.370E+00	2.760E+02

Data file: L26JAL1 Type: STANDARD  
 -----  
 Sample Name: NA Cal. Curve: 07/26/95  
 Date: 28 Jul 1995 20:33 Method: GOALVs Analyst: JB  
 Interface: 707 Cycles: 11 Channel: A  
 -----


Instrument: VARIAN 3500,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog slope (\$ or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/minEE Split Ratio: NA  
 Det 1-Type & Temp: BCD/300C Det 2-Type & Temp: NA  
 -----

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.14 uv offset)



Retention Time	Compound Name	FWW Injected	Area	Height
13.36	GOAL	0.050	1.680E+01	1.380E+03

Data file: L26JA17 Type: STANDARD  
 -----  
 Sample Name: NA Cal. Curve: 07/28/95  
 Date: 28 Jul 1995 23:05 Method: GOALWA Analyst: JB  
 Interface: 707 Cycle#: 17 Channel#: A  
 -----  
 Instrument: VARIAN 1500,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog Slope (S or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 1.5ml/minHE Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA  
 -----  
 Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.15 mv offset)



Retention Time	Compound Name	PPM Injected	Area	Height
13.36	GOAL	0.005	1.540E+00	1.320E+02






Data file: L267A19 Type: STANDARD  
 -----  
 Sample Name: NA Cal. Curve: 07/28/95  
 Date: 28 Jul 1995 23:55 Method: GOALVA Analyst: JB  
 Interface: 707 Cycled: 19 Channel: A  
 -----

Instrument: VARIAN 3500,  
 Column: RTX-300, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog Slope (S or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/minHE Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA  
 -----

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.15 mv offset)




Retention Time	Compound Name	PFW Injected	Area	Height
13.36	GOAL	0.150	5.270E+01	4.880E+03

Data file: L26JA30 Type: STANDARD


Sample Name: MA Cal. Curve: 07/26/95  
Date: 29 Jul 1995 00:21 Method: GOALVa Analyst: JB  
Interface: 707 Cycle#: 20 Channel#: A

Instrument: VARIAN 3500,  
Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
Prog Slope (0 or Linear): MA Inj Port Temp: 265  
Flowrate/Gas: 1.5ml/minHE Split Ratio: MA  
Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: MA

Plot times: 0 to 30 minutes  
Plot range: 50 millivolts (-.16 mv offset)



Retention Time	Compound Name	PPM Injected	Area	Height
13.36	GOAL	0.150	5.250E+01	4.640E+03



Concentrations in report are calculated from equation:  
 $HEIGHT = B0 + B1(UG/ML) + B2(UG/ML)^2$   
 obtained by least-squares fit of standard injection data.




Data file: L267A6 RAR number: 94-0142  
 Method file: GOALVn Sample No: 001  
 Type: FORTIFICATION Component: FIG

Sample Name: MA Cal. Curve: 07/20/95  
 Date: 20 Jul 1995 16:27 Method: GOALVn Analyst: JB  
 Interface: 707 Cycles: 6 Channel: A

Instrument: VARIAN 1500,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog Slope (0 or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/min Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.17 mv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Temp. WT.	ug Added
13.37	GOAL	2.970E+00	.249E+03	5.00	5.00	0.050

Data file: L267A7 RAR number: 94-0142  
 Method file: GOALVa Sample No: 001  
 Type: FORTIFICATION Component: FIG

Sample Name: MA Cal. Curve: 07/28/95  
 Date: 28 Jul 1995 18:52 Method: GOALVa Analyst: JB  
 Interface: 707 Cycle#: 7 Channel#: A

Instrument: VARIAN 3500,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog Slope (S or Linear): NA Inj Port Temp: 255  
 Flowrate/Gas: 3.5ml/minHE Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot time: 0 to 30 minutes  
 Plot range: 50 millivolts (-.12 mv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Temp Wt.	ug Added
5.00	GOAL	2.990E+00	.346E+03	5.00	5.00	0.050

Data file: L267A8 RAR number: 94-0142  
 Method file: GOALVa Sample No: 001  
 Type: FORTIFICATION Component: FIG

Sample Name: MA Cal. Curve: 07/28/95  
 Date: 28 Jul 1995 19:18 Method: GOALVa Analyst: JB  
 Interface: 707 Cycle#: 8 Channel#: A

Instrument: VARIAN 3500,  
 Column: RTX-300, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog Slope (S or Linear): MA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/min Split Ratio: MA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: MA

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.14 mv offset)




Ret. Time	Compound Name	Peak Area	Peak Weight	Volume (ml)	Samp Wt.	ug Added
5.00	GOAL	2.920E+00	.341E+03	5.00	5.00	0.050



```

Data file: L367AS                      RAW number: 94-0142
Method file: GOALVs                     Sample No: 001
Type: PORTIFICATION                  Component: FIG
-----
Sample Name: MA                          Cal. Curve: 07/28/95
Date: 28 Jul 1995 19:43                 Method: GOALVs      Analyst: JS
Interface: 707                          Cycle#: 9           Channel#: A
-----
Instrument: VARIAN 3500,
Column: RTX-200, 0.32mm ID, 1.0um df   Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C       Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C            End Time-Temp (deg-min): 250 C-17
Prog Slope ($ or Linear): NA           Inj Port Temp: 265
Flowrate/Gas: 3.5ml/min                Split Ratio: NA
Det 1-Type & Temp: ECD/300C            Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 80 millivolts (-.15 mv offset)
    
```



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Resp Wt.	ug Added
13.37	GOAL	2.770E+00	.234E+03	5.00	5.00	0.050

Data file: L26JAL2  
 Method file: GOALVA  
 Type: FORTIFICATION

RAR number: 94-0142  
 Sample No: 001  
 Component: FIG

Sample Name: MA  
 Date: 28 Jul 1995 20:59  
 Interface: 707

Cal. Curve: 07/28/95

Method: GOALVA


Analyst: JB *JB*

Cycles: 12

Channel: A

Instrument: VARIAN 3500.  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog slope (# or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/minHE Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.095 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Sampl. Wt.	ug Added
13.37	GOAL	6.140E+00	.499E+03	5.00	5.00	0.100

Data file: L26JA13 RAR number: 94-0142  
 Method file: GOALVA Sample No: 001  
 Type: FORTIFICATION Component: FIG

---


Sample Name: MA Cal. Curve: 07/20/95  
 Date: 28 Jul 1995 21:24 Method: GOALVA Analyst: JB  
 Interface: 707 Cycle#: 13 Channel#: A

---

Instrument: VARIAN 3500,  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog Slope (S or Linear): NA Inj Port Temp: 365  
 Flowrate/Gas: 3.5ml/minHE Split Ratio: NA  
 Det 1-Type & Temp: BCD/300C Det 2-Type & Temp: NA

---

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.061 mv offset)




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
5.00	GOAL	5.780E+00	.482E+03	5.00	5.00	0.100

Data file: L26J114 RAR number: 94-0142  
Method file: GOALVA Sample No: 001  
Type: PORTIFICATION Component: FIG

-----  
Sample Name: MA Cal. Curve: 07/28/95  
Date: 28 Jul 1995 21:49 Method: GOALVA Analyst: JB  
Interface: 707 Cycle#: 14 Channel#: A

-----  
Instrument: VARIAN 3500,  
Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
Prog Slope (0 or Linear): NA Inj Port Temp: 265  
Flowrate/Gas: 3.5ml/min Split Ratio: NA  
Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA


-----  
Plot times: 0 to 30 minutes  
Plot range: 50 millivolts (-.11 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
5.00	GOAL	1.440E+01	.118E+04	5.00	5.00	0.250

```

Data file: L26JA15                      RAR number: 94-0142
Method file: GOALVa                      Sample No: 001
Type: FORTIFICATION                      Component: FIG
-----
Sample Name: MA                          Cal. Curve: 07/28/95
Date: 28 Jul 1995 22:15                 Method: GOALVa      Analyst: JB
Interface: 707                           Cycle#: 15          Channel#: A
-----
Instrument: VARIAN 3500,
Column: RTX-300, 0.32mm ID, 1.0um df   Column Length: 60 Meters
Start Temp-Time (deg-min): 215 C       Ramp Hold (deg-min): 1
Program Rate (deg/min): 10C            End Time-Temp (deg-min): 250 C-17
Prog Slope (S or Linear): MA           Inj Port Temp: 265
Flowrate/Gas: 1.5ml/minHE              Split Ratio: MA
Det 1-Type & Temp: ECD/300C            Det 2-Type & Temp: NA
-----
Plot times: 0 to 30 minutes
Plot range: 50 millivolts (-.093 uv offset)
    
```




Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Temp. Wt.	ug Added
5.00	GOAL	2.760E+01	.232E+04	5.00	5.00	0.500

Data file: L36JA15 RAR number: 94-0142  
 Method file: GOALVa Sample No: 001  
 Type: FORTIFICATION Component: FIG

Sample Name: MA Cal. Curve: 07/28/95  
 Date: 28 Jul 1995 22:40 Method: GOALVa Analyst: JB  
 Interface: 707 Cycle#: 16 Channel#: A

Instrument: VARIAN 3500.  
 Column: RTX-200, 0.32mm ID, 1.0um df Column Length: 60 Meters  
 Start Temp-Time (deg-min): 215 C Ramp Hold (deg-min): 1  
 Program Rate (deg/min): 10C End Time-Temp (deg-min): 250 C-17  
 Prog Slope (0 or Linear): NA Inj Port Temp: 265  
 Flowrate/Gas: 3.5ml/minHE Split Ratio: NA  
 Det 1-Type & Temp: ECD/300C Det 2-Type & Temp: NA

Plot times: 0 to 30 minutes  
 Plot range: 50 millivolts (-.15 mv offset)



Ret. Time	Compound Name	Peak Area	Peak Height	Volume (ml)	Samp Wt.	ug Added
13.36	GOAL	6.850E+00	.561E+03	100	5.00	2.50

Calibration Standards  
07/28/95

Concentration of GOAL	Peak Height
0.150	4580.0
0.150	4640.0
0.100	2590.0
0.100	2590.0
0.050	1380.0
0.050	1290.0
0.010	266.0
0.010	276.0
0.005	127.0
0.005	132.0

Least Squares Calibration Equation  
Standard Curve Number 07/28/95

$$\text{Height} = 75.758 + 19035 \cdot (\text{Ug/ML}) + 72899 \cdot (\text{Ug/ML})^2$$

94-0141 Fortifications

Sample Num.	Comp.	Peak Height	Ug/ML GOAL Found	Final Volume (ML)	Control Corr. (Ug)	PPM GOAL Found	PPM GOAL Added	Pct. Recov.
001	FIG	249.0	0.0088	5.00	0.00	0.0088	0.010	88.
001	FIG	246.0	0.00866	5.00	0.00	0.00866	0.010	86.6
001	FIG	261.0	0.00841	5.00	0.00	0.00841	0.010	84.1
001	FIG	224.0	0.00806	5.00	0.00	0.00806	0.010	80.6
001	FIG	499.0	0.0206	5.00	0.00	0.0206	0.020	103.
001	FIG	482.0	0.0199	5.00	0.00	0.0199	0.020	99.5
001	FIG	1180.0	0.0489	5.00	0.00	0.0489	0.050	97.8
001	FIG	2320.0	0.0881	5.00	0.00	0.0881	0.100	88.1
001	FIG	861.0	0.0234	100	0.00	0.468	0.500	93.6

94-0141 Samples

Sample Num.	Comp.	Peak Height	Ug/ML GOAL Found	Final Volume (ML)	Ug GOAL Found	Sample Weight (G)	Recov Fact.	PPM GOAL Found
001	FIG	0.0	0.00	5.00	0.00	5.00	1.00	0.00

GC/MS RUN SHEET

Date 11/9/95  
 Analyst Shaozhi zheng  
 Instrument GC/MS #4  
 Analytical Method TR 34-94-150  
 Notebook Ref. 002-165  
 Crop fruit (Grape, peach)

Compound Goal  
 Data Reduction Method NA  
 Raw Data File 11085  
 Date Extracted 11-9-95  
 Date Injected 11/9-10/95  
SDS 11/18/95

The preceding run(s) apply to the following analysis

Inj. No	STD (ug/mL)	RAR No.	Spl. No	Sample Component	Sample Weight (g)	Final Volume (mL)	Fortification Level (ug)	Notes
1	0.5							
2	0.1							
3	0.05							
4	0.025							
5		94-0117	001	fruit	5	2.5		BIK
6							0.125	SPKA
7							0.25	SPKB
8							0.625	SPKC
9		92-0100	005					BIK
10							0.125	SPKA
11							0.25	SPKB
12							0.625	SPKC
13	WASH							
14	0.025							
15	0.05							
16	0.1							
17	0.5							
18								
19								
20								
21								
22								
23								
24								
25								
26								


NOTES: Samples were stored in freezer #10A before injected.  
 \*\* 21 is not included in the standard curve. SD 12/18/95  
 94-0117 (peach) Amis Martin  
 92-0100 (grape) 3/19/96  
 Signature Shaozhi zheng




Thu Nov 09 10:41:22 1995

C:\HPCHEM\1\5972\TARGET.U

TR 34 95 111



Scan: 10.00 - 700.00 Samples: 8 Thresh: 100 Step: 0.10  
 300 peaks Base: 218.95 Abundance: 840448



Mass	Abund	Rel Abund	Iso Mass	Iso Abund	Iso Ratio
69.00	773760	100.00	70.00	8631	1.12
218.95	840448	108.62	219.95	37408	4.45
501.95	6749	0.87	502.95	737	10.92

TARGET MASS:	69	131	219	502
DYNAMIC ENT OFFSET:	7.0	10.3	14.6	10.3
TARGET ABUND(%) :	100.0	74.0	115.0	0.9
ACTUAL TUNE ABUND(%) :	100.0	77.6	108.6	0.9

Analytical Data Report for GC-MS Confirmation of Goal in Grape and Peach

Protocol No.:  
Matrix:  
Extraction Date:

34P-95-92  
Grape and Peach  
11/8/95

CAL. Study No.:  
Analyte:  
Analysis Date:

002-165  
Goal  
11/9-10/95


52 12/18/95

Peak	Height	Area	Height	Area	Height	Area
2400	2465	4496	5254	9344	11274	61774

F(x) = 126 \* Amt - 1440  
Corr Coef = 0.999

R & H	Samples ID	CAL No.	Final Vol. ml	Dilution Factor	Sample Vol. (ul)	Peak Area	ng/ml Found	ug Found	ppm Found	ppm Added	% Recovery
94-0117-001*	952309 Btk		2.5	1	5.0	0	0.00	0.000	0.000	NA	NA
94-0117-001	952309 Spk A		2.5	1	5.0	3890	42.26	0.106	0.021	0.025	85
94-0117-001	952309 Spk B		2.5	1	5.0	8952	82.39	0.206	0.041	0.050	82
94-0117-001	952309 Spk C		2.5	1	5.0	30752	255.18	0.638	0.128	0.125	102
92-0100-005*	954668 Btk		2.5	1	5.0	0	0.00	0.000	0.000	NA	NA
92-0100-005	954668 Spk A		2.5	1	5.0	4397	46.28	0.116	0.023	0.025	93
92-0100-005	954668 Spk B		2.5	1	5.0	9405	85.98	0.215	0.043	0.050	86
92-0100-005	954668 Spk C		2.5	1	5.0	32497	269.02	0.673	0.135	0.125	108

\* 94-0117 (peach)  
92-0100 (grape) Deminor-Master  
3/19/96




Method Name: C:\HPCHEM\1\METHODS\GOAL1.M  
Calibration Table Last Updated: Fri Nov 10 08:35:09 1995

TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S1.D  
Operator : sz  
Acquired : 9 Nov 95 5:35 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: std 0.5ug/ml  
Misc Info : C103195-2  
Vial Number: 1

000195



GC-MSD INSTRUMENTATION

- a. Instrumentation: Hewlett-Packard model 5890 Series II Gas Chromatograph/model 5972 mass selective detector
- b. Column: Rtx-50, (0.25 mm i.d. x 30 m, 0.5 µm film thickness)
- c. Oven Temperature: Hold at 100°C for 2 min., then 100 to 270°C at 10°C/min., then hold 15 min.
- d. Injector Temperature: 280°C
- e. Detector Temperature: 300°C
- f. Carrier Gas: Helium
- g. Carrier Gas Flow Rate: 2.18 ml/min.
- h. Head Pressure: 20 Psi
- i. Injection Mode: Splitless
- j. Injection Liner: Silanized Single Taper
- k. Injector Purge Delay: 1.0 min.
- l. Septum Purge: 2 ml/min.
- m. Injection Volume: 2 µl
- n. Ionization Potential: 70 eV
- o. Electron Multiplier Voltage: 1900 to 2318 V
- p. Dwell Time: 100 msec

Initials sz  
Date 11/17/95  
Run # N08S1 to N08S17

---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S1.D  
Operator : SZ  
Acquired : 9 Nov 95 5:35 pm using AcqMethod GOAL1  
Sample Name: std 0.5ug/ml  
Misc Info : C103195-2  
Vial Number: 1  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding:  
\* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---

Number of Compounds in Database: 1

## Time Reference Peaks:

Compound	Expected RT	Actual RT
-----	-----	-----


TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S1.D  
 Operator : sz  
 Acquired : 9 Nov 95 5:35 pm using AcqMethod GOAL1  
 Sample Name: std 0.5ug/ml  
 Misc Info : C103195-2  
 Vial Number: 1  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 682.50  
 Pk # and Type: 1 RTEINT used

19.47	19.96	19.47	19.96	19.47	19.96
Target		Q1		Merged	


	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.22	95162	sys def
Q1	361.00	11.1	0.0- 22.0	19.65	to	10593	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.21	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S2.D  
Operator : sz  
Acquired : 9 Nov 95 6:13 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: std 0.1ug/ml  
Misc Info : C103195-3  
Vial Number: 2

000198



---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S2.D  
Operator : sz  
Acquired : 9 Nov 95 6:13 pm using AcqMethod GOAL1  
Sample Name: std 0.1ug/ml  
Misc Info : C103195-3  
Vial Number: 2  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding:  
\* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---

Number of Compounds in Database: 1


## Time Reference Peaks:

Compound	Expected RT	Actual RT
-----	-----	-----




File : C:\HPCHEM\1\DATA\N08S2.D  
 Operator : sz  
 Acquired : 9 Nov 95 6:13 pm using AcqMethod GOAL1  
 Sample Name: std 0.1ug/ml  
 Misc Info : C103195-3  
 Vial Number: 2  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 119.45  
 Pk # and Type: 1 RTEINT used




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.22	11274	sys def
Q1	361.00	11.5	0.0- 22.0	19.65	to	1294	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.21	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S3.D  
Operator : sz  
Acquired : 9 Nov 95 6:51 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: std 0.05ug/ml  
Misc Info : C103195-4  
Vial Number: 3

000201



---

Information from Data File:  
File : C:\HPCHEM\1\DATA\N08S3.D  
Operator : SZ  
Acquired : 9 Nov 95 6:51 pm using AcqMethod GOAL1  
Sample Name: std 0.05ug/ml  
Misc Info : C103195-4  
Vial Number: 3  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

Quantitation Settings:

Reference Peak Window:	10.00 Percent
Non-Reference Peak Window:	5.00 Percent
Correlation Window:	0.02 Minutes
Default Multiplier:	1.00
Default Sample Concentration:	0.00
Peak Type Decoding:	* -> Time Reference Peak m -> Manually re-integrated 00X -> Compound references ISTD X

---

Number of Compounds in Database: 1

---

Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

TR 34 95 111


000203

File : C:\HPCHEM\1\DATA\N08S3.D  
 Operator : sz  
 Acquired : 9 Nov 95 6:51 pm using AcqMethod GOAL1  
 Sample Name: std 0.05ug/ml  
 Misc Info : C103195-4  
 Vial Number: 3  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 79.04  
 Pk # and Type: 1 RTEINT used




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.22	5254	sys def
Q1	361.00	12.6	0.0- 22.0	19.65	to	662	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.21	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S4.D  
Operator : SZ  
Acquired : 9 Nov 95 7:30 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: std 0.025ug/ml  
Misc Info : C103195-5  
Vial Number: 4

000204



---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S14.D  
Operator : sz  
Acquired : 10 Nov 95 1:52 am using AcqMethod GOAL1  
Sample Name: std 0.025ug/ml  
Misc Info : C103195-5  
Vial Number: 14  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding:  
\* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---

Number of Compounds in Database: 1


## Time Reference Peaks:

Compound	Expected RT	Actual RT
-----	-----	-----


TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S4.D  
 Operator : sz  
 Acquired : 9 Nov 95 7:30 pm using AcqMethod GOAL1  
 Sample Name: std 0.025ug/ml  
 Misc Info : C103195-5  
 Vial Number: 4  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 30.45  
 Pk # and Type: 1 RTEINT used




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.16	2400	sys def
Q1	361.00	12.6	0.0- 22.0	19.65	to	303	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S5.D  
Operator : sz  
Acquired : 9 Nov 95 8:08 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: 94-0117-001  
Misc Info : blk SW=5g FV=2.5ml  
Vial Number: 5

000207





---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S5.D  
Operator : SZ  
Acquired : 9 Nov 95 8:08 pm using AcqMethod GOAL1  
Sample Name: 94-0117-001  
Misc Info : blk SW=5g FV=2.5ml  
Vial Number: 5  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding: \* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---

Number of Compounds in Database: 1


## Time Reference Peaks:

Compound	Expected RT	Actual RT
-----	-----	-----


TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S5.D  
 Operator : sz  
 Acquired : 9 Nov 95 8:08 pm using AcqMethod GOAL1  
 Sample Name: 94-0117-001  
 Misc Info : blk SW=5g FV=2.5ml  
 Vial Number: 5  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 20.06  
 Concentration: 12.30  
 Pk # and Type: 1 RTEINT used  
 # Qualifiers Not Satisfied




	Signal	Ratios	Limits	RT	Limits	Resp	Integ	Type
Tgt	252.00	100.0%		20.06	19.16	110		sys def
Q1	361.00	0.0	0.0- 22.0	0.00	to	0		sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0		sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0		sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S6.D  
Operator : SZ  
Acquired : 9 Nov 95 8:46 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: 94-0117-001  
Misc Info : spka SW=5g FV=2.5ml 0.025ppm  
Vial Number: 6

000210



---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S6.D  
Operator : SZ  
Acquired : 9 Nov 95 8:46 pm using AcqMethod GOAL1  
Sample Name: 94-0117-001  
Misc Info : spka SW=5g FV=2.5ml 0.025ppm  
Vial Number: 6  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding: \* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---


Number of Compounds in Database: 1

## Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

File : C:\HPCHEM\1\DATA\N08S6.D  
 Operator : sz  
 Acquired : 9 Nov 95 8:46 pm using AcqMethod GOAL1  
 Sample Name: 94-0117-001  
 Misc Info : spka SW=5g FV=2.5ml 0.025ppm  
 Vial Number: 6  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 42.26  
 Pk # and Type: 1 RTEINT used




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	3890	sys def
Q1	361.00	10.5	0.0- 22.0	19.65	to	410	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S7.D  
Operator : sz  
Acquired : 9 Nov 95 9:24 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: 94-0117-001  
Misc Info : spkb SW=5g FV=2.5ml 0.05ppm  
Vial Number: 7

000213



## Information from Data File:

File : C:\HPCHEM\1\DATA\N08S7.D  
Operator : SZ  
Acquired : 9 Nov 95 9:24 pm using AcqMethod GOAL1  
Sample Name: 94-0117-001  
Misc Info : spkb SW=5g FV=2.5ml 0.05ppm  
Vial Number: 7  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding:  
\* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

Number of Compounds in Database: 1


## Time Reference Peaks:

Compound	Expected RT	Actual RT
-----	-----	-----


TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S7.D  
 Operator : sz  
 Acquired : 9 Nov 95 9:24 pm using AcqMethod GOAL1  
 Sample Name: 94-0117-001  
 Misc Info : spkb SW=5g FV=2.5ml 0.05ppm  
 Vial Number: 7  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 82.39  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	8952	sys def
Q1	361.00	12.4	0.0- 22.0	19.65	to	1110	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def






TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S8.D  
Operator : sz  
Acquired : 9 Nov 95 10:03 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: 94-0117-001  
Misc Info : spkc SW=5g FV=2.5ml 0.125ppm  
Vial Number: 8

000216



---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S8.D  
Operator : sz  
Acquired : 9 Nov 95 10:03 pm using AcqMethod GOAL1  
Sample Name: 94-0117-001  
Misc Info : spkc SW=5g FV=2.5ml 0.125ppm  
Vial Number: 8  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding:  
\* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---


Number of Compounds in Database: 1

## Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

File : C:\HPCHEM\1\DATA\N08S8.D  
 Operator : SZ  
 Acquired : 9 Nov 95 10:03 pm using AcqMethod GOAL1  
 Sample Name: 94-0117-001  
 Misc Info : spkc SW=5g FV=2.5ml 0.125ppm  
 Vial Number: 8  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 255.18  
 Pk # and Type: 1 RTEINT used




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	30752	sys def
Q1	361.00	11.3	0.0- 22.0	19.65	to	3485	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S9.D  
Operator : sz  
Acquired : 9 Nov 95 10:41 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: 92-0100-005  
Misc Info : blk SW=5g FV=2.5ml  
Vial Number: 9

000219



---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S9.D  
Operator : SZ  
Acquired : 9 Nov 95 10:41 pm using AcqMethod GOAL1  
Sample Name: 92-0100-005  
Misc Info : blk SW=5g FV=2.5ml  
Vial Number: 9  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding: \* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---


Number of Compounds in Database: 1

## Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

File : C:\HPCHEM\1\DATA\N08S9.D  
 Operator : sz  
 Acquired : 9 Nov 95 10:41 pm using AcqMethod GOAL1  
 Sample Name: 92-0100-005  
 Misc Info : blk SW=5g FV=2.5ml  
 Vial Number: 9  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.37  
 Concentration: 12.33  
 Pk # and Type: 1 RTEINT used  
 # Qualifiers Not Satisfied




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.37	19.16	114	sys def
Q1	361.00	0.0	0.0- 22.0	0.00	to	0	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S10.D  
Operator : sz  
Acquired : 9 Nov 95 11:19 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: 92-0100-005  
Misc Info : spka SW=5g FV=2.5ml 0.025ppm  
Vial Number: 10

000222



---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S10.D  
Operator : SZ  
Acquired : 9 Nov 95 11:19 pm using AcqMethod GOAL1  
Sample Name: 92-0100-005  
Misc Info : spka SW=5g FV=2.5ml 0.025ppm  
Vial Number: 10  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding: \* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---

Number of Compounds in Database: 1


## Time Reference Peaks:

Compound	Expected RT	Actual RT
-----	-----	-----




File : C:\HPCHEM\1\DATA\N08S10.D  
 Operator : sz  
 Acquired : 9 Nov 95 11:19 pm using AcqMethod GOAL1  
 Sample Name: 92-0100-005  
 Misc Info : spka SW=5g FV=2.5ml 0.025ppm  
 Vial Number: 10  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 46.28  
 Pk # and Type: 1 RTEINT used




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	4397	sys def
Q1	361.00	13.5	0.0- 22.0	19.64	to	592	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S11.D  
Operator : sz  
Acquired : 9 Nov 95 11:57 pm using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: 92-0100-005  
Misc Info : spkb SW=5g FV=2.5ml 0.05ppm  
Vial Number: 11

000225



Information from Data File:

File : C:\HPCHEM\1\DATA\N08S11.D  
 Operator : sz  
 Acquired : 9 Nov 95 11:57 pm using AcqMethod GOAL1  
 Sample Name: 92-0100-005  
 Misc Info : spkb SW=5g FV=2.5ml 0.05ppm  
 Vial Number: 11  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Quantitation Settings:

Reference Peak Window: 10.00 Percent  
 Non-Reference Peak Window: 5.00 Percent  
 Correlation Window: 0.02 Minutes  
 Default Multiplier: 1.00  
 Default Sample Concentration: 0.00  
 Peak Type Decoding:  
 \* -> Time Reference Peak  
 m -> Manually re-integrated  
 00X -> Compound references ISTD X


Number of Compounds in Database: 1

Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

File : C:\HPCHEM\1\DATA\N08S11.D  
 Operator : sz  
 Acquired : 9 Nov 95 11:57 pm using AcqMethod GOAL1  
 Sample Name: 92-0100-005  
 Misc Info : spkb SW=5g FV=2.5ml 0.05ppm  
 Vial Number: 11  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 85.98  
 Pk # and Type: 1 RTEINT used




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	9405	sys def
Q1	361.00	11.5	0.0- 22.0	19.65	to	1084	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S12.D  
Operator : sz  
Acquired : 10 Nov 95 12:35 am using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: 92-0100-005  
Misc Info : spkc SW=5g FV=2.5ml 0.125ppm  
Vial Number: 12

000228



---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S12.D  
Operator : SZ  
Acquired : 10 Nov 95 12:35 am using AcqMethod GOAL1  
Sample Name: 92-0100-005  
Misc Info : spkc SW=5g FV=2.5ml 0.125ppm  
Vial Number: 12  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding: \* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---


Number of Compounds in Database: 1

## Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

File : C:\HPCHEM\1\DATA\N08S12.D  
 Operator : sz  
 Acquired : 10 Nov 95 12:35 am using AcqMethod GOAL1  
 Sample Name: 92-0100-005  
 Misc Info : spkc SW=5g FV=2.5ml 0.125ppm  
 Vial Number: 12  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 269.02  
 Pk # and Type: 1 RTEINT used




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.16	32497	sys def
Q1	361.00	11.4	0.0- 22.0	19.65	to	3689	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S13.D  
Operator : sz  
Acquired : 10 Nov 95 1:14 am using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: wash  
Misc Info : na  
Vial Number: 13

000331





---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S13.D  
Operator : sz  
Acquired : 10 Nov 95 1:14 am using AcqMethod GOAL1  
Sample Name: wash  
Misc Info : na  
Vial Number: 13  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding:  
\* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---

Number of Compounds in Database: 1

## Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

TR 34 95 111


000233

File : C:\HPCHEM\1\DATA\N08S13.D  
 Operator : SZ  
 Acquired : 10 Nov 95 1:14 am using AcqMethod GOAL1  
 Sample Name: wash  
 Misc Info : na  
 Vial Number: 13  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.85  
 Concentration: 12.42  
 Pk # and Type: 1 RTEINT used  
 # Qualifiers Not Satisfied




	Signal	Ratios	Limits	RT	Limits	Resp	Integ	Type
Tgt	252.00	100.0%		19.85	19.16	125		sys def
Q1	361.00	0.0	0.0- 22.0	0.00	to	0		sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0		sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0		sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S14.D  
Operator : SZ  
Acquired : 10 Nov 95 1:52 am using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: std 0.025ug/ml  
Misc Info : C103195-5  
Vial Number: 14

000234



---

Information from Data File:  
File : C:\HPCHEM\1\DATA\N08S4.D  
Operator : sz  
Acquired : 9 Nov 95 7:30 pm using AcqMethod GOAL1  
Sample Name: std 0.025ug/ml  
Misc Info : C103195-5  
Vial Number: 4  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

Quantitation Settings:  
Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding:  
\* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---

Number of Compounds in Database: 1


---

Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

File : C:\HPCHEM\1\DATA\N08S14.D  
 Operator : SZ  
 Acquired : 10 Nov 95 1:52 am using AcqMethod GOAL1  
 Sample Name: std 0.025ug/ml  
 Misc Info : C103195-5  
 Vial Number: 14  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 30.97  
 Pk # and Type: 1 RTEINT used




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.16	2465	sys def
Q1	361.00	13.3	0.0- 22.0	19.65	to	328	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.14	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S15.D  
Operator : sz  
Acquired : 10 Nov 95 2:30 am using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: std 0.05ug/ml  
Misc Info : C103195-4  
Vial Number: 15

000237



---

Information from Data File:  
File : C:\HPCHEM\1\DATA\N08S15.D  
Operator : sz  
Acquired : 10 Nov 95 2:30 am using AcqMethod GOAL1  
Sample Name: std 0.05ug/ml  
Misc Info : C103195-4  
Vial Number: 15  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

Quantitation Settings:

Reference Peak Window:	10.00 Percent
Non-Reference Peak Window:	5.00 Percent
Correlation Window:	0.02 Minutes
Default Multiplier:	1.00
Default Sample Concentration:	0.00
Peak Type Decoding:	* -> Time Reference Peak m -> Manually re-integrated 00X -> Compound references ISTD X

---

Number of Compounds in Database: 1


---

Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

File : C:\HPCHEM\1\DATA\N08S15.D  
 Operator : SZ  
 Acquired : 10 Nov 95 2:30 am using AcqMethod GOAL1  
 Sample Name: std 0.05ug/ml  
 Misc Info : C103195-4  
 Vial Number: 15  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 73.95  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.22	4496	sys def
Q1	361.00	13.3	0.0- 22.0	19.65	to	596	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.21	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def






TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S16.D  
Operator : sz  
Acquired : 10 Nov 95 3:08 am using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: std 0.1ug/ml  
Misc Info : C103195-3  
Vial Number: 16

60000



---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S16.D  
Operator : SZ  
Acquired : 10 Nov 95 3:08 am using AcqMethod GOAL1  
Sample Name: std 0.1ug/ml  
Misc Info : C103195-3  
Vial Number: 16  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding: \* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---


Number of Compounds in Database: 1

## Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

File : C:\HPCHEM\1\DATA\N08S16.D  
 Operator : SZ  
 Acquired : 10 Nov 95 3:08 am using AcqMethod GOAL1  
 Sample Name: std 0.1ug/ml  
 Misc Info : C103195-3  
 Vial Number: 16  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.64  
 Concentration: 106.49  
 Pk # and Type: 1 RTEINT used




	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.64	19.22	9344	sys def
Q1	361.00	12.2	0.0- 22.0	19.64	to	1138	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.21	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



TR 34 95 111

File : C:\HPCHEM\1\DATA\N08S17.D  
Operator : sz  
Acquired : 10 Nov 95 3:46 am using AcqMethod GOAL1  
Instrument : GC/MS # 4  
Sample Name: std 0.5ug/ml  
Misc Info : C103195-2  
Vial Number: 17

000243



---

Information from Data File:

File : C:\HPCHEM\1\DATA\N08S17.D  
Operator : sz  
Acquired : 10 Nov 95 3:46 am using AcqMethod GOAL1  
Sample Name: std 0.5ug/ml  
Misc Info : C103195-2  
Vial Number: 17  
CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

---

## Quantitation Settings:

Reference Peak Window: 10.00 Percent  
Non-Reference Peak Window: 5.00 Percent  
Correlation Window: 0.02 Minutes  
Default Multiplier: 1.00  
Default Sample Concentration: 0.00  
Peak Type Decoding: \* -> Time Reference Peak  
m -> Manually re-integrated  
00X -> Compound references ISTD X

---


Number of Compounds in Database: 1

## Time Reference Peaks:


Compound	Expected RT	Actual RT
-----	-----	-----

File : C:\HPCHEM\1\DATA\N08S17.D  
 Operator : SZ  
 Acquired : 10 Nov 95 3:46 am using AcqMethod GOAL1  
 Sample Name: std 0.5ug/ml  
 Misc Info : C103195-2  
 Vial Number: 17  
 CurrentMeth: C:\HPCHEM\1\METHODS\GOAL1.M

Compound: Goal  
 Ret Time: 19.65  
 Concentration: 458.40  
 Pk # and Type: 1 RTEINT used



	Signal	Ratios	Limits	RT	Limits	Resp	Integ Type
Tgt	252.00	100.0%		19.65	19.22	61774	sys def
Q1	361.00	11.2	0.0- 22.0	19.65	to	6923	sys def
Q2	0.00	0.0	0.0- 0.0	0.00	20.21	0	sys def
Q3	0.00	0.0	0.0- 0.0	0.00		0	sys def



CAL STUDY NO. <u>002-165</u>
SPONSOR PROTOCOL NO. <u>34P-95-92</u>

DAILY STANDARD INJECTIONS ARE COMPARED AGAINST THE CURRENT CALIBRATION FILE. THIS FILE IS GENERALLY BASED ON THE PREVIOUS RUN'S STANDARDS. THE CALIBRATION FILE IS UPDATED BASED ON THE RUN'S STANDARD, AND SAMPLE ANALYSES ARE QUANTITATED USING THE UPDATED CALIBRATION FILE.

*shanzhi zheng 1/4/96*