



Q&A #4

Date: October 2, 2014

Project: Lighting Retrofit / T12 to T8 Fluorescent Fixtures Conversion

Bidders must make sure that their bids are based on the latest version of the tender documents published and take into consideration the following amendments and information, including any information provided on amendments or Q&As previously published for this project.

Bidders that do not comply with this requirement will be discarded.

Question 10: Concerning Addendum #2 (2014/09/30)

At the beginning of the addendum you modify the quantity of 2 tubes fixtures from 603 to 528, therefore 75 less.

Later in the addendum, at page 4, you indicate 70 fixtures must be replaced instead rather than be retrofitted.

My question is, should we consider there are 145 fixtures less? Or the 70 fixtures to replaced are included in the 75 less? Please clarify.

Answer 10: There are 528 2-tubes fixtures to convert to T8.
There are 5 2-tubes fixtures to convert to LED (refrigerated rooms).
There are 70 2-tubes fixtures to replace.

ALL OTHER TERMS AND CONDITIONS REMAIN UNCHANGED