



Annex A

Statement of Work (SOW)

Service Agreement for Medical Device Inspection Program

1. OBJECTIVE

Correctional Service of Canada (CSC) intends to contract a qualified firm to provide expert inspection services for the medical equipment located at Correctional Service of Canada (CSC) federal correctional facilities in the Atlantic Region (Nova Scotia and New Brunswick).

This requirement involves the provision of technical support to Correctional Service of Canada (CSC) at the following sites:

- Dorchester Complex Minimum Sector, Dorchester (NB)
- Dorchester Complex Medium Sector, Dorchester (NB)
- Shepody Healing Center, Dorchester (NB)
- Atlantic Institution, Renous (NB)
- Springhill Institution, Springhill (NS)
- Nova Institution, Truro (NS)
- Regional Pharmacy, Moncton (NB)
- Regional Headquarters, Moncton (NB)

2. BACKGROUND

2.1 Project Background

The Health Services Sector of CSC has developed a Medical Device Inspection Program (MDIP) which outlines the processes and responsibilities for the effective inspection, upgrading, and replacement of medical devices and equipment used in CSC Health Services clinics, Hospitals, Pharmacies and Regional Head Quarters (RHQ).

The MDIP ensures that medical devices and equipment are safe and functional. It also provides guidance to identify and address potential problems with medical devices and equipment that may result in injury to staff or patients and requires a follow up on medical equipment or devices involved in an adverse event.



To support the MDIP, CSC requires the services of expert resources to conduct an annual inspection of all relevant medical equipment that falls within the inspection profile of the MDIP (hereinafter referred to as “Equipment”).

2.2 Objectives

As part of its mandate, Correctional Service Canada (CSC) has built and maintains a large number of secure sites with a diverse security profile. In order to provide the requisite health services to the inmate population within these sites, each Institution has a wide variety of medical systems and components that, in all aspects, constitutes a mission critical infrastructure that must be maintained.

To that end, CSC has a requirement to contract technical support for the inspection of the medical systems and associated equipment described in Annexes C, D, E, F, G, H, I, J.

3. SCOPE

3.1 Nature of Work

The work involves the provision of qualified resources to conduct the inspection of the Equipment identified by CSC in its federal correctional facilities (hereinafter referred to as “Institutions”). The Contractor must perform the Work, as described in this Statement of Work (SOW).

3.2 Nature of Services

The Contractor will provide technical support services designed and delivered to ensure the quality and effectiveness of CSC Equipment and associated components.

3.3 Institutional Medical Equipment

The medical devices (Equipment) within an Institution’s Health Services clinic can be broken down into two distinct categories:

1. Equipment that should be repaired when it ceases to function; and
2. Equipment that should be replaced when it ceases to function.



3.4 Equipment at End of Life

In all of the above cases, the useful service life of Equipment is established by the recommendations of the manufacturer and captured by Health Services in its Lifecycle Management Program. The Technician must assess the lifecycle status of each piece of Equipment based on the manufacturer's recommended service life as well as through the performance of a physical examination. If the Equipment is at or near the end of its useful service life, the Technician must note that in his or her report to the chief of health or project authority at RHQ.

4. INSPECTION PERSONNEL

4.1 General

The Contractor must perform this work using specialized resources with expertise in the inspection of medical equipment (hereinafter referred to as "Technicians").

The Contractor's Technicians will be responsible for:

1. Travelling to the required Institutions
2. Providing an annual report to the chief of health care, Regional Pharmacist or Manager of Clinical Services and the project authority at RHQ.

4.2 Technician Qualifications

Each Technician performing the Work must be an employee of the Contractor or an employee of a subcontractor.

Each Technician must have graduated from a recognized Biomedical Electronics Technician program.

Each Technician must at least five (5) years experience performing inspection on medical equipment.

4.3 Technician Certifications

Many manufacturers of medical equipment provide certification programs for the inspection and service of its product lines (including the inspection thereof). The Contractor must provide Technician(s) to complete the Work that have the

necessary certifications and training to inspect all of the Equipment appearing in Annexes C, D, E, F, G, H, I, J of this document.



It is reasonable that the Contractor may provide more than one Technician to an Institution at disparate times to ensure that all of the Equipment is maintained by someone with the proper certifications and training to perform that part of the Work.

5. TASKS AND ACTIVITIES

5.1 Test Equipment and Tools

The provision, inspection, management and control of all required test equipment will be the responsibility of the Contractor. A complete list of required testing equipment shall be presented to the chief of health prior to each site visit.

It should be noted that CSC will have final approval of all proposed test equipment before it is deployed on-site. This is due to a possible situation where some forms of equipment qualify as contraband within the precincts of an Institution and are therefore not permitted on-site.

5.2 Travel to CSC Institutions

The Contractor's Technician(s) must travel to Institutions identified by CSC as a facility that has Equipment that requires Inspection. The Contractor must arrange for the travel and accommodations for the Technician(s) and will remain entirely responsible for the travel requirements of its Technicians throughout the project.

The Contractor must ensure that the target Institution receives prior notification of the planned Inspection activities at least ten (10) working days prior to the site visit and the Technician has undergone a CPIC check in advance of that date.

The Technician will be accompanied by an individual assigned by the Institution to escort the Technician throughout the Institution. The Contractor must ensure that the Technician will abide by the instructions regarding the escort given by the Institutional Head or its delegate to the Technician.

6. CONSTRAINTS

The following constraints must be taken into account in the performance of the Work as described in this document:



6.1 Security Requirements

Any person entering a correctional Institution must undergo a clearance performed by the Canadian Police Information Centre (CPIC). This mandatory "CPIC Check" must be initiated at least two (2) weeks in advance of the scheduled date of entry into the Institution and completed without issue before the person is permitted to enter.

Delays caused by the CSC security clearance requirement will not incur any additional costs to CSC.

7. DELIVERABLES

The contractor shall provide the following services with regards to medical devices listed in Annexes C, D, E, F, G, H, I, J.

- a. Incoming on-site inspections to be conducted on an annual basis. Time of inspection will be determined after contract award between the Chief Health Services, Regional pharmacist or Manager of Clinical Services and the contractor. Each inspection shall include industry standard verifications for each piece of equipment listed in Annexes C, D, E, F, G, H, I, J.

Equipment that undergoes inspection must be tested to ensure it meets quality control standards as per the manufacturer's guidelines.

- b. Annual performance testing, inspection and calibration of electro medical devices is required in order to comply with standards set forth by the respective provinces.
- c. During inspection, assist in identifying requirements. This will be in regards to replacing existing equipment or purchasing new equipment.
- d. During inspection, repair all issues or problems with the equipment which do not require any parts.
- e. The report must be provided in writing **within 30 days of the site inspection**, and must detail all pieces of equipment inspected, any concerns identified and/or repairs recommended and or completed.

List of equipment

From time to time, equipment may be added or removed from the lists of equipment when required during the period of the contract. These additions or deletions may come from the respective site authority (Chief of Health Services) or from the Project Authority at the Regional office.



Any changes greater than 10% in the total number of items to be inspected as part of the contract or any changes greater than 30% in the number of items at anyone specific site will require a review of the inspection fees to take into consideration the additional or reduced workload. Should this situation occur, the Project Authority will request an updated quote from the contractor for work to be done.

NOTE: This clause excludes the addition or removal of diagnostic imaging systems and equipment. Any addition or removal of such equipment would require a new quote.

Response Time:

The Contractor shall provide these services Monday to Friday from 8:00 a.m. to 4:30 p.m.

Due to the nature of the environment in which the programs will be delivered, the scheduled dates of inspections/repairs will be subject to change due to urgent institutional circumstances such as institutional lockdowns, closures due to inclement weather, security requirements, or other emergency situations at no cost to the CSC. In these cases, the Project Authority or his/her representative will advise the Contractor and work will be scheduled.

Travel and Living Expenses

Travel and Living Expenses will be the responsibility of the contractor to and from CSC work sites. Vehicle transportation, driver's permits, insurances are also the responsibility of the contractor.

Project Authority

All work will be performed subject to the inspection and acceptance by the Project Authority or the on-site Chief of Health Services.

Project Authority name and phone number will be provided in the contract document.

The Project Authority or his designated representative will be responsible for monitoring the progress of the work and will be responsible for the technical requirements; the acceptance and approval of the deliverables. Any proposed changes to the scope of work may be discussed with the Project Authority, but any resulting changes can only be authorized by an amendment issued by the Contracting Authority.



APPENDIX 'A' – Glossary of Terms and Acronyms

Term	Definition
Contractor	A private business entity that is awarded the contract to perform the Work for Canada.
Equipment	Medical equipment at CSC Institutions that falls within the inspection profile of the Medical Device Maintenance Program.
Institution	Correctional facilities under the care of Correctional Service Canada where incarcerated offenders serve criminal sentences of two years or more as inmates. Institutions are secure facilities not broadly open to the public.
Technician	A specialized resource provided by the Contractor that is a trained Biomedical Technician with the requisite expertise and experience in the inspection of medical equipment.



APPENDIX 'B' – List of CSC Institutions and Addresses

ATLANTIC REGION

Atlantic Institution: Attention – Chief of Health Care

13175 Route 8, P.O. Box 102
Renous, New Brunswick
E9E 2E1

Dorchester Complex Minimum Sector: Attention – Chief of Health Care

4902 Main Street
Dorchester, New Brunswick
E4K 2Y9

Dorchester Complex Medium Sector: Attention – Chief of Health Care

4902 Main Street
Dorchester, New Brunswick
E4K 2Y9

Nova Institution: Attention – Chief of Health Care

180 James Street
Truro, Nova Scotia
B2N 6R8

Springhill Institution: Attention – Chief of Health Care

330 McGee Street, P.O. Box 2140
Springhill, Nova Scotia
B0M 1X0

Shepody Healing Centre: Attention – Nursing Supervisor

4902 Main Street
Dorchester, New Brunswick
E4K 2Y9

Correctional Service Canada: Attention – Regional Pharmacist or Manager

Clinical Services Health Services (RHQ – Atlantic)

1045 Main Street
Moncton, New Brunswick
E1C 1H1

Appendix C - Equipment List Dorchester Minimum

Item	~~Equipment Code~~Code d'équipement~~	~~Institution~~Établissement~~	~~Location~~Localisation~~	~~Manufacturer~~Nom du fabricant~~	~~Model~~Modèle~~	~~Serial Number~~Numéro de série~~
Aspirator	221-006303	Westmorland Institution	Dr room	Gomco	300	L10395
Blanket warmer	221-002028	Westmorland Institution	Medication Room	Pedigo	P-2010S	545084-000
Blood drawing chair	221-002015	Westmorland Institution	Nurses Office	Winco	574	574100103
Camscope angled wand	221-002023	Westmorland Institution	Dr room	Sometech	DCS-104T (NTSC)	DSCN-11B-061
Camscope receiver	221-002023	Westmorland Institution	Dr room	Sometech	DCSM (XGA)	DCSM-11C-029
Camscope straight wand	221-002023	Westmorland Institution	Dr room	Sometech	P101	DCSM-11C-115
Caulk Dentsply	221-0005000	Westmorland Institution	Dentists rm	Dentsply International	VM-D	71970
Centrifuge	221-002014	Westmorland Institution	Nurses Office	Thermo Scientific	CL-2	42619248
Dental Chair	221-009078	Westmorland Institution	Dentists rm	Belmont	X-Caliber	AN09A0073
Compressor	221-009082	Westmorland Institution	storage rm	Midmark	277OCE50-660A	
Curing Light		Westmorland Institution	Dentists rm	Coltene/Whaledent INC.	C7970100115	110507017
Dental cart equip	221-002027	Westmorland Institution	Dentists rm	Aidec	PACI	C912329
Dental chair equip		Westmorland Institution	Dentists rm	Belmont	BDS0050	13344
Dental X-Ray Controler	221-009079	Westmorland Institution	Dentists rm	Belmont	96	EX08K0779
Dental X-Ray Head	221-009079	Westmorland Institution	Dentists rm	Belmont	096-H	EH08K0779
Electrocardiograph	221-002006	Westmorland Institution	Dr room	Welch Allyn Canada Ltd.	CP10	105068614163
Grinder / polisher	221-000535	Westmorland Institution	Dentists rm	Red-wing	26A	816-227 84
Interferential Ultrasound	221-002007	Westmorland Institution	Dr room	Dynatronics Research Corp.	D850PLUS	860726
Examining Light	221-002009	Westmorland Institution	Dr room	Welch Allyn Canada Ltd.	LS200	ref # 74227

Examining Light	221-002008	Westmorland Institution	Dr room	Welch Allyn Canada Ltd.	48100	
Monitor - NIBP	221-002013	Westmorland Institution	Dr room	Critikon / Johnson & Johnson	100V2	043M0299009
O2 sat monitor	221-002017	Westmorland Institution	Nurses Office	Masimo Corporation	RAD5	500255
Polycom Practitioner Cart	221-001623	Westmorland Institution	Dr room	Polycom	PPCE-120BP	L100614525
Prophylaxis unit	221-000499	Westmorland Institution	Dentists rm	Dental Equipment	2002	103-05362
O2 sat monitor		Westmorland Institution	Dr room	Pulse Oximeer	M50E09CS0	9060429
Cast Saw	221-002021	Westmorland Institution	Dr room	Stryker Canada	840	92111273
Sterilizer	221-002024	Westmorland Institution	Dentists rm	L&R	T-14B	Z85 83120
Stretcher	221-002020	Westmorland Institution	Nurses Office	Midmark	K520	TAV003996
Physio Treatment Table	221-002005	Westmorland Institution	Dr room	Cardon Rehabilitation & medical equipment	R28535	039-08011
Vibration machine		Westmorland Institution	Dentists rm			G-10-1-4-V
Dental X-ray Arm	221-009079	Westmorland Institution	Dentist rm	Belmont	096-A	EA08K0779
IV Pump	221-009028	Westmorland Institution	Nurses office	Plum A+	84095270002	13898345
Mister	221-006303	Westmorland Institution	Dr. room	Schuco	5711102	188532
Aerosol Compressor		Westmorland Institution	Storage rm	Airlife	RT500	169547
Monitor - NIBP	221-002026	Westmorland Institution	Nurses office	Critikon/Johnson & Johnson	DCP400N-CE	2019205-001
Ophthalmoscope	221-002018	Westmorland Institution	Nurses office	Welch Allyn Canada Ltd.	74710	92192
Ophthalmoscope		Westmorland Institution	Dr room	Welch Allyn Canada Ltd.	767	5290
Temperature Monitor	221-002018	Westmorland Institution	Nurses office	Welch Allyn Canada Ltd.	690	4106441
Temperature Monitor	221-002010	Westmorland Institution	Dr room	Welch Allyn Canada Ltd.	690	09361171
Curing Light		Westmorland Institution	Dentist rm	3M ESPE	Elipar S10	939112017014
Aerosol Compressor		Westmorland Institution	Storage rm	Invacare	IRC1196	99H328019
Geriatric chair		Westmorland Institution	Storage Room	AMG Medical Inc.	770-050	10860-0609-0003
Vaccine Refrigerator		Westmorland Institution	Inf Dis Office	Panasonic	SR-L6111W-PA	121101390
Dental Assistant Chair		Westmorland Institution	Dentist rm	Belmont	093G	36911
Dental Assistant Chair		Westmorland Institution	Dentist rm	Belmont	092Q	36910
Dental Light		Westmorland Institution	Dentist rm	Belmont	AL501T5CR3	AV09A0040
Aerosol Compressor		Westmorland Institution	Storage rm	Airlife	RT500	169547

Appendix D – Equipment List Dorchester Medium

Item	~~Equipment Code~~Code d'équipement~~	~~Institution~~Établissement~~	~~Location~~Localisation~~	~~Manufacturer~~Nom du fabricant~~	~~Model~~Modèle~~	~~Serial Number~~Numéro de série~~
AED	220-009019	Dorchester Institution	Treatment room 1	Zoll	AED+	X08D152026
Aspirator 690	220-009018	Dorchester Institution	Treatment room 1	Contemporary Products	6260	AS21745
Assembly cam-14		Dorchester Institution	Treatment room 1	GE Medical Systems	CAM-14	RNC11015386PA
Battery pack for head lamp	220-009033	Dorchester Institution	X-ray room	Heine	M Pack	2157354
Bladder scan	220-009012	Dorchester Institution	Treatment room 2	Bladder Scan	0570-0188	PA011179
Blanket Warmer		Dorchester Institution	Medical Supply Room, 1st floor	Pedigo	P-20105	545088-000
Blood pressure/Pulse/O2 sat monitor	220-009181	Dorchester Institution	Storage room 3rd floor	Critikon	Dynamap Compact BP monitor	ref. 117207
Blood Pressure/Temp/O2 Sat Monitor	220-0090004	Dorchester Institution	Treatment room 1	GE Medical Systems	Dynamap DPC400N-CE	2019205-001
Camscope box	220-009022	Dorchester Institution	Nursing supervisor's office	Sometec	DCSM(XGA)	DCSM-10G-167
Camscope wand (angled)	220-009022	Dorchester Institution	Nursing supervisor's office	Sometec	DCS-104T (NTSC)	DCSN-10E-217
Camscope wand (straight)	220-009022	Dorchester Institution	Nursing supervisor's office	Sometec	P101	DCSN-11C-116
Cassette autoclave	220-006702	Dorchester Institution	Dentist's lab	SciCan	1102	228756A
Centrifuge	220-007436	Dorchester Institution	Treatment room 1	Thermo IEC	CENTRA CL-2	42613405
Chart Projector	220-009030	Dorchester Institution	X-ray room	Marco	--	22143
Compressed Air		Dorchester Institution	Dentist's office	Belmont	4000	20213
Compressor		Dorchester Institution	Dentist's office	Dominic Hunter	DHDRD10	3701300004
Controls for X-ray machine	220-009036	Dorchester Institution	X-ray room	GE Medical System	2320993	07115877
Curing Light		Dorchester Institution	Dentist's office	Coltène/Whaledent	Coltolux 75	060128666
Dental Cart	220-009024	Dorchester Institution	Treatment room 1	--	447192	C912329
Dental Chair	220-009023	Dorchester Institution	Dentist's office	Takara Belmont	447193	012057301

Dental Chair Instruments		Dorchester Institution	Dentist's office	Belmont		12057301
Dental light	220-009025	Dorchester Institution	Dentist's office	Belmont	AL501T	T040461
Dental X-ray processor	220-009026	Dorchester Institution	Dentist's office	Dell	Studio Hybrid 140g	18HVBMI
Diagnostic X-ray wallstand	220-009035	Dorchester Institution	X-ray room	GE Hualun Medical System	600-0301	50997HL6
Digitizer	220-009034	Dorchester Institution	X-ray room	AGFA	CR30-X	
Drill - Console - Surgical	220-009024	Dorchester Institution	Dentist's office	--	--	--
Eclipse Proteus Collimator	220-009038	Dorchester Institution	X-ray room	GE Hualun Medical System	2379827	53339HL8
Electro therapy unit	220-009015	Dorchester Institution	Treatment room 2	Chattanooga Canada	Intellect Adv	4978
Electrocardiograph	220-009002	Dorchester Institution	Treatment room 1	GE Medical Systems	MAC5500	11-41583
Emergency Cart	220-009042	Dorchester Institution	Storage room 3rd floor	Harloff Co.		46269
Head lamp	220-009033	Dorchester Institution	X-ray room	Heine	--	1130545
Ice maker		Dorchester Institution	B-6 End of hallway	Whirlpool	GI15NDXTB1	EW1501128
Laser module	220-009037	Dorchester Institution	Treatment room 2	Chattanooga Canada	2762CC	2033
Lensmeter		Dorchester Institution	X-ray room	Marco	201	M2474
Lift - Patient - Floor		Dorchester Institution	Range	BHM Medical	Ergolift	ERLI2746
Lift - Patient - Floor		Dorchester Institution	Range		0102-8851	96511054100
Light - Examining	220-009014	Dorchester Institution	Treatment room 2	Welch Allyn Canada Ltd	44200	--
Monitor		Dorchester Institution	Nursing supervisor's office	Samsung	LN26B460B2D	Z1HC3CXZ105915A
Monitor - O2 sat	220-009003	Dorchester Institution	Treatment room 2	Masimo Corporation	RAD5	532646
Monitor - O2 sat	220-009008	Dorchester Institution	Treatment room 1	Masimo Corporation	RAD5	530779
Ophthaomalogy Chair	220-009031	Dorchester Institution	X-ray room	Power Tilt	FX920L	93010707007
Practitioner cart		Dorchester Institution	Nursing supervisor's office	Polycom	PPCE-120BP	N091215032
Projector - Ophthalmic	220-009030	Dorchester Institution	X-ray room	Marco	50	M22143
Proteus elevating table	220-009038	Dorchester Institution	X-ray room	GE Hualun Medical System	2259988	50508HL9
Proteus System Console	220-009038	Dorchester Institution	X-ray room	GE Hualun Medical System	2259976	50529HL7
Proteus XR-A System	220-009038	Dorchester Institution	X-ray room	GE Hualun Medical System	2259973-2	50567HL7
Pump - Infusion	220-006718	Dorchester Institution	Medical supply room 1st floor	Abbott Laboratories	PLUMXL 11555-04-01	11204399

IV Pump - Infusion	220-009001	Dorchester Institution	Treatment room 1	Hospira	1239-104-03	13898441
IV Pump - Infusion	220-009028	Dorchester Institution	Medical supply room 1st floor	Hospira	1239-104-03	13898345
RAD-14X-ray-tube	220-009038	Dorchester Institution	X-ray room	GE Hualun Medical System	2259981	95296-7U
Refractor	220-0004032	Dorchester Institution	X-ray room	Marco	RT-I	10975
Refrigerator	220-007096	Dorchester Institution	Pharmacy	Maytag	MTB1956BEW	11101861WD
Ritter automatic sterilizer		Dorchester Institution	Dentist's lab	Midmark	M9 UltraClave M9-022	V546354
Saw - Cast	220-009016	Dorchester Institution	Treatment room 2	Stryker	8127	E34984
Solid state ultrasonic sterilizer	220-000201	Dorchester Institution	Dentist's lab	L & R Manufacturing	T-14B	150
Stand - optometrist equipment	220-009032	Dorchester Institution	X-ray room	Reliance Medical Products	7900	79010702014
Stretcher	220-006439	Dorchester Institution	Treatment room 1	Midmark	520	
Table - Physio-treatment	220-009013	Dorchester Institution	Treatment room 2	Chattanooga Canada	Triton DTS TRT-600	2214
Tonometer		Dorchester Institution	X-ray room	Clement Clarke	Perkins tonometer MK2	1999T2073
Traction unit		Dorchester Institution	Treatment room 2	Chattanooga Canada	2841/3336	
Trolley		Dorchester Institution	Treatment room 1	GE Medical Systems	MAC 5000	SCD11051932PA
Ultra sound applicator		Dorchester Institution	Treatment room 2	Chattanooga Canada	27335	37018
UV inspection lamp		Dorchester Institution	Treatment room 1	Luxo	0316040	89705
Vaccine Refridgerator	220-007643	Dorchester Institution	ID Nurse's office	Symphony	128934	B89237-35605
Vital signs monitor	220-009011	Dorchester Institution	Treatment room 2	Welch Allyn	Spot Vital Signs Lxi	20110201548
X-ray Controller	220-009027	Dorchester Institution	Dentist's office	Belmont	097-CM	EX10J0141
X-ray Head	220-009027	Dorchester Institution	Dentist's office	Belmont	097-H	EH10J0510
X-ray printer	220-009039	Dorchester Institution	X-ray room	Sharp AGFA	Drystar 5300	
X-ray Regular Cone	220-009027	Dorchester Institution	Dentist's office	Belmont	097-R	ER10J0510
Stryker Stretcher	220-009087	Dorchester Institution		Stryker	0747-000	1105034705
Stryker Stretcher	220-009088	Dorchester Institution		Stryker	0747-000	1105034706
Slit Lamp	221-006503	Dorchester Institution	Dr room	Top Con PS-20A	SL-2E	6211049
Refrigerator - Pharmacy 17.2 cu ft	220-07992	Dorchester Institution	Secondary Pharmacy		SAN-MPR514	12010029
Refrigerator - Pharmacy	220-07993	Dorchester Institution	IDN Office		SAN-SRL6111W	120200257
Aspirator / Suction Machine		Dorchester Institution	Trauma Room	Gomco by Allied	405	20080331005

IV Pump		Dorchester Institution	Trauma Room	Hospira	1239-104-07	17446358
Integrated Diagnostic System		Dorchester Institution	Trauma Room	Welch Allyn	GS 777 WA 7791-MX	

Appendix E – Equipment List Shepody Healing Center

Item	~~Equipment Code~~Code d'équipement~~	~~Institution~~Établissement~~	~~Location~~Localisation~~	~~Manufacturer~~Nom du fabricant~~	~~Model~~Modèle~~	~~Serial Number~~Numéro de série~~
MONITOR SPO2	BMD-7967	Shepody Healing Centre	Acute Unit	O-Two Systems International Inc	01MN1000	
MONITOR SPO2	BMD-7968	Shepody Healing Centre	Acute Unit	O-Two Systems International Inc	01MN1000	
MONITOR NIBP SPO2 TEMP	BMD-7963	Shepody Healing Centre	Acute Unit	Welch Allyn Canada Ltd	42NTB	
MONITOR TEMP	BMD-7964	Shepody Healing Centre	Acute Unit	Welch Allyn Canada Ltd	690	
MONITOR TEMP	BMD-7966	Shepody Healing Centre	Acute Unit	Welch Allyn Canada Ltd	690	
MONITOR TEMP TYMPANIC	BMD-7969	Shepody Healing Centre	Acute Unit	Welch Allyn Canada Ltd	PRO4000	
SPHYGMOMANOMETER ANALOG	BMD-7970	Shepody Healing Centre	Acute Unit	Heine	GAMMA G5	
Vaccine fridge Acute	BMR-11881	Shepody Healing Centre	Acute Unit	Haier (LG Electronics)		NSA36LACG
MONITOR NIBP SPO2 TEMP	BMD-7973	Shepody Healing Centre	PSR Unit	Welch Allyn Canada Ltd	42NTB	
SCALE PATIENT x2	BMD-7976	Shepody Healing Centre	PSR Unit	Health-O-Meter	402KL	4020079413
MONITOR TEMP	BMD-5216	Shepody Healing Centre	PSR Unit	Welch Allyn Canada Ltd	690	
MONITOR SPO2	BMD-7975	Shepody Healing Centre	PSR Unit	O-Two Systems International Inc	01MN1000	
JAMAR DYNAMETER		Shepody Healing Centre	PSR Unit			
Stretcher	BMD-11882	Shepody Healing Centre	PSR Unit	Techlem	500	TMRR1303029
Vaccine fridge PSR	BMD-11883	Shepody Healing Centre	PSR Unit	LG	GC151RCW	

Appendix F – Equipment List Atlantic Institution

Item	Equipment Code Code d'équipement	Institution Établissement	Location Localisation	Manufacturer Nom du fabricant	Model Modèle	Serial Number Numéro de série
FRIDGE	200-004118	Atlantic Institution	Pharmacy		428730	EW3317486
Fridge- Vaccine	231-004141	Atlantic Institution	Pharmacy	Panasonic	SR-L6111W-PA	120700476
Auto Reader	231-3335	Atlantic Institution	D306	3M Attest	305807	208585
Defibrillator, Automated, E	231-004502	Atlantic Institution	Hall		403018	X08B147213
Electrocardiograph	231-003476	Atlantic Institution	D303	Welch Allyn Canada Ltd.	CD200	20004343
TONOMETER	231-000285	Atlantic Institution	D104	Topcon (Koenigkramer)	SL4E	643280
Light-examining	231-000273	Atlantic Institution	D303	Heraeus	55002286	100005
Monitor - SP02		Atlantic Institution		Nonin	8500	500257365
Monitor - SP02	231-004527	Atlantic Institution	Unit 3	Nonin	8500	500714515
Monitor NIBP-SP02-Temp- Heart Rate	231-003761	Atlantic Institution		General Electric Mesystems Limited	2009284001	AAW080502335
Pump - Infusion	231-001972	Atlantic Institution	D303	Abbot Laboratories	XL	11203144
Pump - Infusion	231-004501	Atlantic Institution	D303	Hospira Worldwide Inc.	123910403	13898447
Strecher Ferno	231-3270	Atlantic Institution	Hall	FERNO	287857	785076
Fridge- Vaccine	231-003427	Atlantic Institution	Pharmacy	Follett	128934	B89235-35605
Warmer, Blanket	231-3582	Atlantic Institution	D306	Pedigo	398442	545085-00
MONITOR-BP-ECG-NIBP-REC-SPO2	231-004512	Atlantic Institution	Unit 3	Dinamap	V100-PBXA-AXAA-DX	SH6111400455A
OPHTHALMIC INSTRUMENT STAND	231-004504	Atlantic Institution	D104	Topcon (Koenigkramer)	IS2000	351541
CLEANER-ULTRASONIC	231-004519	Atlantic Institution	D308	L & R	T14B	Q87122194T
STERILIZER-STEAM-TABLE-TOP	231-001988	Atlantic Institution	D308	STATim	1102	228426A
CHAIR-DENTAL	231-003832	Atlantic Institution	D308	DENTECH	ULTIMA2	CM2315-3-K1

XRAY-UNIT	231-003766	Atlantic Institution	D308	Sinclair Dental	BELRAY 096	EX08K0769
TABLE-PHYSIO-TREATMENT	231-004508	Atlantic Institution	D304	Centry Medical Dist	TRT600	2245
TRACTION-UNIT	231-004509	Atlantic Institution	D304	Triton-Chattanooga	MP1	3358
TABLE-EXAMINING	231-002886	Atlantic Institution	NOT IN USE - In Basement(Tunnel)		056-605	2702002001
BED-ELECTRIC	231-003076	Atlantic Institution	D2017 - Dry Cell	Hill Room	P8400D3867	E086AL6139
LENSOMETER	231-000281	Atlantic Institution	D104		322056	322056
ASPIRATOR-PORTABLE	231-004510	Atlantic Institution	Unit 3		L190-GR	NONE
THERAPY SYSTEM - Physio	231-004506	Atlantic Institution	D304	Intelect Advanced	2762CC	4879
PUMP-COMPRESSOR		Atlantic Institution	D303	Nedquip	RT500	170682
ASPIRATOR	231-004505	Atlantic Institution	D303	Advantage	L190-GR	NONE
DEFIB-AED	231-004503	Atlantic Institution	Hall	AEDPLUS	X08D152305	7478539
MONITOR-SPO2	231-004513	Atlantic Institution	Unit 3		8500	500327436
Optometrist Chair	231-004576	Atlantic Institution	D104	TOPCON		
Autoscaler		Atlantic Institution	Room D308	South East- Henry Schein		
Fridge	CSC-004229	Atlantic Institution		Danby	2306	
Dr. Camscope - Oscope		Atlantic Institution	D112	Sometech Inc	DCS-104 (NTSC)	DCSM-118-059
Dr. Camscope - Intra Oral Camera		Atlantic Institution	D112	Sometech Inc.	DCS-P101	DCSN-11C-118
Dr. Camscope - C.V.P. System		Atlantic Institution	D112	Sometech Inc.	DCSM (XGA)	DCSM-11C-030
AMD-2500		Atlantic Institution	D112	AMD Global Telemedicine	AMD-2500	
AMD-3700		Atlantic Institution	D112	AMD Global Medecine		21151
AMD-400		Atlantic Institution	D112	AMD Global Telemedicine		
Hot Water Bath	231-308	Atlantic Institution		Eagle		
Lathe	231-000312	Atlantic Institution	D308	BDMI		816-22-1287
Bed, Examination	23100-004036	Atlantic Institution	D112	Global Medco Ltd.	9105-30	
Bed, Examination	23100-004035	Atlantic Institution	Unit 3	Global Medco Ltd.	9105-30	
Stretcher	231- 003579	Atlantic Institution	D303	Stryker	1007	
Amalgamator	231-24502	Atlantic Institution	D308	Caulk		

Light - Examining (Optical)	231-0001697	Atlantic Institution	D104	Hein		
Phoropter (Optical)	231-001698	Atlantic Institution	D104	Reichert	11625 AR	10806-2
KEROMETER (Optical)	231-284	Atlantic Institution	D104	Bausch & Lomb	B & L 71-21-35	18256
Microscope Double Zoom (Optical)	231-000284	Atlantic Institution	D104	Amtek	CP 40	105072
Dry Imager	23100-004521	Atlantic Institution	D304	AGFA	Drystar 5300	
Proteus XR/A	23100-004515	Atlantic Institution	D304	Ge	AL01F	02246
X-Ray Table	23100-04516	Atlantic Institution	D304	GE		
Digitizer	23100-004522	Atlantic Institution	D304	Agfa	CR 30-X	
Chest X-Ray Machine	231-289	Atlantic Institution	D304		260-076-R	1462
Sterilizer	23100-004140	Atlantic Institution	D308	Statim	Statim 2000 G4	101212j00026
Refrigerator, Laboratory	23100-004141	Atlantic Institution	Pharmacy	Panasonic	SR-L6111W-PA	120700476
Console Xray Unit		Atlantic Institution	D304	GE	2259976	
Light Ultraviolet		Atlantic Institution		DSP	RADII PLUS	
Monitor Temp Tympanic		Atlantic Institution		Welch Allyn	Braun 6021	
Tube - XRay Unit		Atlantic Institution	D304	GE	2269646	
Monitor SPO2	23100 004169	Atlantic Institution	Unit 5	NONIN	SPECTRO2	4047289
MONITOR-BP-ECG-NIBP-REC-SPO2	23100-004170	Atlantic Institution	Unit 5	Welch Allyn Canada Ltd.		20140201184
MONITOR-BP-ECG-NIBP-REC-SPO2	23100-004172	Atlantic Institution	Unit 5	Welch Allyn Canada Ltd.		20140201183
Otto/Ophthalmoscope-Handheld	23100-004171	Atlantic Institution	Unit 5	Welch Allyn Canada Ltd.		05259-M
STRETCHER-FERNO	23100-004142	Atlantic Institution	Hallway	Ferno	35X proflex	12N232511
Monitor SPO2		Atlantic Institution	Unit 5	Nonin	8000AA	500716022
Stretcher	23001-004142	Atlantic Institution	Hallway	FERNO		238-1493-01B
Centrifuge		Atlantic Institution	Pharmacy	Unico	Power Spin FX	250-520-02
Stretcher		Atlantic Institution	Basement-For use in Unit 5	FERNO	287784	BBB113616

Appendix G – Equipment List Springhill Institution

Item	~~Equipment Code~~Code d'équipement~~	~~Institution~~Établissement~~	~~Location~~Localisation~~	~~Manufacturer~~Nom du fabricant~~	~~Model~~Modèle~~	~~Serial Number~~Numéro de série~~
AMTER KERATOMETER	210-003190	Springhill Institution	Optometrist	Topcon	238714	01074G
Light Head - HEINE OMEGA 180 INDIRECT OP	210-003429	Springhill Institution	Optometrist	Heine	268299	0/11848
TOPCON COUNTERBALANCED MIDD	210-003194	Springhill Institution	Optometrist	Topcon	238725	203100017B
TOPCON HALOGEN CHART PROJEC	210-003193	Springhill Institution	Optometrist	Topcon	238723	262933G
TOPCON MOTORIZED OPLTHALMIC Chair	210-003191	Springhill Institution	Optometrist	Topcon	238717	222466
TOPCON POHTHALMIC STAND	210-003192	Springhill Institution	Optometrist	Topcon	238721	350236
TOPCON VISION TESTER	210-003189	Springhill Institution	Optometrist	Topcon	238712	3617864
TONOMETER	210-001322	Springhill Institution	Optometrist	TONO-PENXL		T5932
Dina care Procure (BP temp, O2 sats)	210-005006	Springhill Institution	Treatment	GE	2009284001	AAW080502158A
ECG Machine	210-004364	Springhill Institution	Treatment	GE	MAC5500	SCD10087696PA
Gomco Suction Aspirator	210-005009	Springhill Institution	Treatment	Gomco	4005	20071218026
Portable Suction-Aspirator	210-005004	Springhill Institution	Treatment	Devil Bless	7305	PD235889
HILL ROM STRETCHER	210-003427	Springhill Institution	Treatment	Hill Rom	268261	F093AW7612
IV Pump	210-005008	Springhill Institution	Treatment	Hospira PlumA+	12910405	13898186
Blanket Warmer	210-005008	Springhill Institution	Treatment	Pedigo	P-2010S	545087
Pulse Oximeter Nellcor NPB4	210-003983	Springhill Institution			365367	G00842987
Belmont Light		Springhill Institution	dentist	Belmont 501T	105000 Dental Light	H71903
Curing Light-Densply		Springhill Institution	dentist	Bobcat	506 QHL75	506-30012
Speed Clave Sterilizer		Springhill Institution	dentist	Midmark	M7011 Speed Clave	MH003986
Ultra Sonic Unit		Springhill Institution	dentist	DENSPLY	Cavitron	11546484
Vacuum, water spurt		Springhill Institution	dentist	Henry Schien		07344-563850
Dental Chair		Springhill Institution	dentist	Belmont	X-Calibur	73203603

X-Ray Head	210-004318	Springhill Institution	dentist	Belmont	428384	EH07I0201
Amalgamator		Springhill Institution	dentist	Ultra Mat 2		71903
Ultrasound (physio)	210-005015	Springhill Institution		Chattanooga Canada	2762CC Intellect Advanced	4295
Hot Wax machine Inva care		Springhill Institution		Chattanooga Canada	Para Care	24753
Steam pack heating unit Hydro collator	210-005014	Springhill Institution		Chattanooga Canada	Hydrocollator -2102C	2328C
X-Ray System	210-000434	Springhill Institution	X-ray	GE	5892919	4034
X-Ray Unit	210-000435	Springhill Institution	X-ray	GE	600-030	46335HL6
X-Ray Table	210-00436	Springhill Institution	X-ray	GE	X171137 Picker	B25863
X-Ray Digital Panoramic System	210-004213	Springhill Institution	X-ray	GENDEX	Orthoralix 8500	35-1858143DP
Medication Refridgerator		Springhill Institution	Medication Room	Sanyo	MPR721	11110815
Physio Treatment Table		Springhill Institution	Doctor's Office	Chattanooga Canada	TRT400	11717
Monitor SPO2		Springhill Institution		Nellcor Puritan Bennett	NPB40	G00842991
Dr. Camscope Intra Oral Camera Telemedicine		Springhill Institution	Back office	Sometech Inc	DCSN	DCSN-11C-117
Dr. Camscope Otoscope - Telemedicine		Springhill Institution	Back Office	SomeTech	DCS-104(NTSC)	DCSN-11B-058
Dr. Camscope C.V.P. System - Telemedicine Control Box		Springhill Institution	Back Office		DCSM(XGA)	DCSN-11B-002
Ophthalmoscope		Springhill Institution	Dr Office	Welch Allyn Canada	767	
Monitor-Temp		Springhill Institution		Welch Allyn	690 Suretemp	743382
Monitor-Temp		Springhill Institution		Welch Allyn	679 Suretemp	3137764
Ophthalmoscope		Springhill Institution		Welch Allyn	767	
Dental Ultra Sonic		Springhill Institution	Dentist	L&R Medical	Quantrex140	196231
Dental Cart		Springhill Institution	Dentist	ADEC		
PUMP Compressor		Springhill Institution	Treatment Room	Devilbiss	5650C Pulmoaide	C170892
Examining Light		Springhill Institution	Treatment Room	Welch Allyn	44300/LS135	

Pump Compressor		Springhill Institution	Treatment	Invacare Corp	IRC1196 - Passport	59400000890
Pump Compressor		Springhill Institution	Treatment	Invacare Corp	IRC1196 - PASSPORT	59400000908
Pump Compressor		Springhill Institution	Treatment	Invacare Corp	IRC1196 - PASSPORT	59400000897
Monitor-NIBP-SP02-Temp		Springhill Institution		Welch Allyn	42MTB Spot Vital	201103450
Ophthalmoscope		Springhill Institution		Welch Allyn	767	
Lensometer		Springhill Institution		Bausch & Lomb	712655	512045
Vaccine Refrigerator - RRC Panasonic Laboratory Refrigerator		Springhill Institution	Nursing Office - RRC	Panasonic	SR-L6111W	120700510
Monitor-NIBP-SP02-TEMP		Springhill Institution	RRC - Nursing Office	Welch Allyn	420 Series	201216120
Monitor-NIBP-SP02-TEMP - Blood Pressure Machine		Springhill Institution	Dr. Office	Welch Allyn	420 Series	201216780
Pulse Oximeter Nellcor N-65		Springhill Institution	Dr. Office	Nellcor	N-65P	P12827521
Vaccine Fridge - Panasonic Vaccine Fridge - Treatment Room		Springhill Institution	Treatment Room	Panasonic		121001043
Cast Saw		Springhill Institution	Treatment Room	BSN Medical Inc.	0295-200	28330
PortaCount	004382	Springhill Institution	Fit Testing Office	Portacount Pro+	8038	80389433
Particle Generator (part of fit tester)	004382	Springhill Institution	Fit Testing Office	Particle Generator	8026	80263697
AGFA		Springhill Institution	X-Ray room	AGFA	CF30X	1380
AGFA		Springhill Institution	X-ray room	AGFA	Drystar 5300	
Pedigo Stretcher		Springhill Institution	Treatment Room	Pedigo	PP7500-W-Spec	0199025
Blood pressure/O2Stat and Temp		Springhill Institution	Treatment Room/Portable	Welch Allyn	42NTB	201403441
Monitor-NIBP-SP02-Temp		Springhill Institution		Welch Allan	42NTB	

MONITOR-TEMP-TYMPANIC		Springhill Institution		Welch Allan	6021	
MONITOR-TEMP-TYMPANIC		Springhill Institution		Welch Allan	6022	
Pump Compressor		Springhill Institution		Innovocorp	IRC1196	
Scale - Patient		Springhill Institution		DET	6437	
Tonometer		Springhill Institution		CLC	MK2	
Wheelchair		Springhill Institution		AMG	700-631	
Wheelchair		Springhill Institution		AMG	700-625	
Wheelchair		Springhill Institution		AMG	700-621	
Wheelchair		Springhill Institution		AMG	700-627	
Wheelchair		Springhill Institution		AMG	700-632	

Appendix H – Equipment List Nova Institution for Women

Item	~~Equipment Code~~Code d'équipement~~	~~Institution~~Établissement~~	~~Location~~Localisation~~	~~Manufacturer~~Nom du fabricant~~	~~Model~~Modèle~~	~~Serial Number~~Numéro de série~~
Amalgam separator	250-001441	Nova Institution	Dental suite (pump room)	Henry Schein	Solmetex/Vacstar/Internat	25669/89824/5858009
Amalgamator		Nova Institution	Dental suite	Silamat	S5	
Aspirator	250-001576	Nova Institution	Tx room	Gomco	4005	20071218030
Canal measurement tool	250-001442	Nova Institution	Dental suite	Sinclair Dental	Root2XII	
Curing light	250-001445	Nova Institution	Dental suite	Demetron	A.2	
Dental chair	250-001449	Nova Institution	Dental suite	Xcalibur		94492407
Dental examination lamp	250-001446	Nova Institution	Dental suite	ALT	508T12N	CN07L0139
Digital scale	250-001572	Nova Institution	Tx room	Health-o-Meter	2842KL	7E620450000069
Doppler	250-001560	Nova Institution	Ax room	Medisonics		FB06949
Drill station	250-001452	Nova Institution	Dental suite	Adec Pac	3420	C888511

ECG	250-001458	Nova Institution	Tx room	Welch Allyn	CP2AS	20004345
Exam table	250-001568	Nova Institution	Secure Unit	Midmark	Ritter 104	
Examination lamp	250-001578	Nova Institution	Tx room	Luxo		
Examination lamp	250-001566	Nova Institution	Secure Unit		LS-135	
Heat source	250-001443	Nova Institution	Dental suite	System B		
Hot pack warmer	250-001118	Nova Institution	Tx room		E-1	137706
ID Fridge	250-001467	Nova Institution	Med Room	Follett	Symphony series	
ID Fridge	250-001557	Nova Institution	Med Room	Follett	REF series	
Infusion pump	250-001577	Nova Institution	Tx room	Plum A+		13897820
Intraoral x-ray unit	250-001450	Nova Institution	Dental suite	Gendex/Visualix		1514B-3945ABG
LCM control module	250-001440	Nova Institution	Dental suite	Sinclair Dental	LUX115V	
Nebulizer	250-001559	Nova Institution	Ax room	Sunmark	4075	045143
Oximeter		Nova Institution	Ax room	black		
Oximeter	250-001117	Nova Institution	Emerg kit		NPB40	G03860573
Oximeter	250-001225	Nova Institution	Secure Unit	Rad-5V		538887
Panoramic x-ray unit	250-001438	Nova Institution	Dental suite	Gendex	8500DDE	35-1858141DP
Physio bed	250-001460	Nova Institution	Tx room	Triton	TRT-600	2217
Pulse rate monitor	250-001462	Nova Institution	Ax room	Dynamap	DPC400N	2019205-001
Pulse rate monitor	250-001265	Nova Institution	Secure Unit	Dynamap	DCP400N	SDT10110608SP
Stretcher	250-001457	Nova Institution	Tx room	Hill Rom	P8000	K043AQ4641
TENS machine	250-001575	Nova Institution	Tx room	Physiotens		5330588
TENS machine	250-001574	Nova Institution	Tx room	Physiotens		447523
Therapeutic ultrasound	250-001460	Nova Institution	Tx room	Intelect	2762CC	4948
Traction unit	250-001460	Nova Institution	Tx room	Chatanooga	2841	3317
Ultraclave	250-001455	Nova Institution	Dental suite	Ritter	Midmark M9-022	V460343
Ultrasonic cleaner (large)	250-001454	Nova Institution	Dental suite	Sulton	Prosonic 1000MT	QDG020710129E

Ultrasonic cleaner (small)	250-001444	Nova Institution	Dental suite	Densply	Cavitron GEN-131	9788
Vibrator		Nova Institution	Dental suite	Buffalo	No. 1A	
X-ray arm	250-001453	Nova Institution	Dental suite	Gendex	Expert DC	12-1852698DP
Exam table	25000-001381	Nova Institution	Ax room	Ritter	Midmark 222	V1168926
Wheelchair (wide)	250-001459	Nova Institution	Healthcare Unit	AMG	700-631	100050-0512-0026
Examination lamp	250-001551	Nova Institution	Ax Room	Welch Allyn	GS300	
Thermometer	250-001553	Nova Institution	Tx Room	Braun	Thermoscan Plus	10105478202
Thermometer	250-001554	Nova Institution	Ax Room	Braun	Thermoscan 6022	11110850664
Handheld BP Monitor	250-001555	Nova Institution	Emerg kit- stockroom	Welch Allyn	3400	100000050811
Scale - Adult	250-001556	Nova Institution	Inmate Bathroom	Detecto		
Integrated Diagnostic System	250-001558	Nova Institution	Ax Room	Welch Allyn		
Thermometer	250-001563	Nova Institution	Emerg kit - stockroom	Braun	6022	11110890735
Nebulizer	250-001565	Nova Institution	Secure Unit	Sunmark	4075	040471
Thermometer	250-001567	Nova Institution	Secure Unit emerg kit	Braun	Thermoscan	11110846300
Scale - infant	250-001573	Nova Institution	Tx Room	Detecto		
Integrated Diagnostic System	250-001580	Nova Institution	Tx Room	Welch Allyn		
Wheelchair (regular)	250-001570	Nova Institution	Tx Room	Medline		
Speculum light with charger		Nova Institution	Ax Room	Welch Allyn	739 series	20100720/20100714
TENS machine		Nova Institution	Tx Room	PhysioSD	TS-1211	11-049683
TENS machine		Nova Institution	Tx Room	PhysioSD	TS-1211	11-049685
Telehealth stethoscope		Nova Institution	Tx Room	Dr. Camscope	AMD-3700	21150
Telehealth control unit		Nova Institution	Tx Room	Dr. Camscope	DCSM(XGA)	DCSM-10G0169
Telehealth opthal/autoscope		Nova Institution	Tx Room	Dr. Camscope	DCS-104T(NTSC)	DCSN-10E-285
Telehealth ultrasound wand	20100-004578	Nova Institution	Tx Room	Dr. Camscope	DCS-P101	DCSN-08B-341

CPAP machine		Nova Institution	Storeroom	ResMed	H5i/S9	22112011401/22111819
Desk viewer, X-ray backlight		Nova Institution	Dental suite	Dux		
Thermometer		Nova Institution	Secure Unit	Braun	6013	10100054947
Oxygen tank regulator		Nova Institution	Emerg kit - stockroom	Medline		
Oxygen tank regulator		Nova Institution	H/C Unit emerg kit	Medline		
Sphygmomanometer		Nova Institution	Tx Room	TYCOS		90121134041
Sphygmomanometer		Nova Institution	Ax Room	TYCOS		9508090077
TENS Machine		Nova Institution		PhysioTens		5330585
Bladder Scanner		Nova Institution	Treatment Room	Verathon	Phantom BVI-6100	SS03792
Hydrocollator		Nova Institution	Treatment Room	Chattanooga	2402	T7783C
Blanket warmer		Nova Institution	Treatment Room	Olympic	Warmette	13507
Sphygmomanometer		Nova Institution	Emerg kit- stockroom	Welch Allyn	analog	
Sphygmomanometer		Nova Institution	Emerg kit - secure unit	Welch Allyn	analog	
Optometry Recline Chair	25000-003999	Nova Institution	Dr office	Axis Medical	CS-418	CP-PZY-0418AJ
Visual Acuity Kit	25000-003998	Nova Institution	Dr office	Huvitz	HDC9000N	9HC00N13I0009
Ophthalmic Stand	25000-003997	Nova Institution	Dr office	Marco	Deluxe II	13I0499
Slit Lamp Microscope	25000-003996	Nova Institution	Dr office	Huvitz	HS-5000	5HS00013I0002
Applanation Tonometer	25000-003995	Nova Institution	Dr office	Axis Medical	R-Type	05117
Cylinders Phoropter	25000-003994	Nova Institution	Dr office	Reichert	Minus	582150613
Half Shadow Mirror		Nova Institution	Dr office	Axis Medical	Plain white	none
Pneumatic Exam Stool		Nova Institution	Dr office	Brewer	none	61609
Wheelchair (regular)		Nova Institution	Stockroom	Airgo	Procare IC	
ECG		Nova Institution	Tx Room	Welch Allyn	CP 100	10021867
Natus Neurology Elite - Doppler		Nova Institution		Natus	200	ECO1776

Appendix I – Equipment List Regional Pharmacy

Item	Equipment Code / Code d'équipement	Institution / Établissement	Location / Localisation	Manufacturer / Nom du fabricant	Model / Modèle	Serial Number / Numéro de série
BIISTER PACKAGER - Automed Heat Sealer	201-0045554	Regional Pharmacy	Packaging Room	Healthmark	Autogen	AG-009C
BIISTER PACKAGER - Automed Heat Sealer	201-004550	Regional Pharmacy	Packaging Room	Healthmark	Autogen	AG-003C
BIISTER PACKAGER - Automed Heat Sealer	201-004552	Regional Pharmacy	Packaging Room	Healthmark	MTS -350	350-TC-004
Blister Packager - Manual HEAT SEALER	210-005023	Regional Pharmacy	Shipping / Receiving	Hunt Corporation	Bienfang	MH-00804
COUNTING MACHINE	201-004553	Regional Pharmacy	Dispensary	McKesson	Parata Accucount	NTEP CC No. 96-116A3P
REFRIGERATOR	201-005024	Regional Pharmacy	Methadone Room	True Food Services	TS-49 - True	4978832
REFRIGERATOR	201-005025	Regional Pharmacy	Dispensary	True Food Services	TS-49 - True	4947312
REFRIGERATOR	201-004595	Regional Pharmacy	Basement	ESBE Scientific	SAN-MPR1411	10020036
REPEATER PUMP	201-004551	Regional Pharmacy	Methadone Room	Baxa	Baxa	2952
REPEATER PUMP	201-004555	Regional Pharmacy	Methadone Room	Baxa	Baxa	2039

Appendix J – Equipment List Regional Headquarters

Item	Equipment Code Code d'équipement	Institution Établissement	Location Localisation	Manufacturer Nom du fabricant	Model Modèle	Serial Number Numéro de série
General Video Exam Scope	20100-004579	Regional Headquarters		Sometech	456384	DCSM-11C-015
IntraOral Hand Piece DCS-P	20100-004578	Regional Headquarters		Sometech	456383	DCSN-08B-341
Refrigerator (SAN-MPR1411)	20100-004597	Regional Headquarters	Currently in storage (basement @ RHQ)	ESBE Scientific	458353	40202019
Telephonic Stethoscope AMD	20100-004547	Regional Headquarters		AMD	448659	SN 21151
Telephonic Stethoscope AMD	20100-004548	Regional Headquarters		AMD	448660	SN 22150