

Questions and Answers:

Annex B 1.3.1 – how many 75 litre nominal capacity clear bags are to be picked up weekly?

- The exact amount is not set and will vary. A weekly estimate provided for the purpose of the tender (see Annex C- Financial Proposal, 1.3) is one (1)- 75 litre bag of glass bottles plus one (1)- 75 litre bag of metal cans plus one (1)- 75 litre bag of type 1&2 plastics. Therefore three (3) 75 litre clear plastic bags per week is the estimate to be used. Each bag will have an estimated weight of 4.536 kg.

Is the collection driver required to collect the bags throughout the facility or collect all bags from the recycling room?

- The Collection driver is required to collect all bags from the recycling room only. They do not have to collect throughout the facility.

Annex B 1.4.1 – how many 75 litre nominal capacity clear bags are to be picked up weekly?

- The exact amount is not set and will vary. A weekly estimate provided for the purpose of the tender (see Annex C- Financial Proposal, 1.4.1) is one (1)- 75 litre bag of shredded material. Therefore one (1)- 75 litre clear plastic bag per week is the estimate to be used.

Annex B 1.5.1 – how many 75 litre nominal capacity clear bags are to be picked up weekly?

- The exact amount is not set and will vary. A weekly estimate provided for the purpose of the tender (see Annex C- Financial Proposal, 1.5.1) is one (1)- 75 litre bag of white Styrofoam material. Therefore one (1)- 75 litre clear plastic bag per week is the estimate to be used.