





Adjustable Vacuum Switch
Adjustable Vacuum Switch
Outlet
Outlet
Vacuum Alarm
Vacuum Meter (0; -1 bar)




Vacuum Pump
Vacuum Pump



Test Valve 1/2"



Connection Points	
Connection	Size
Inlet	DN40
Outlet	DN40

<div>Gertsen & Olufsen A/S</div> <div>Tlf: +45 4817 4055 - Solvang 22 - 3450 - DK</div>						General Tolerances: DS/ISO 2768-vl		General welding: DS/ISO 5817:2007 -D		
						Scale	1 : 12		Projection	
						Format	A3		Model No.	
<div>Title</div> <div>Assembly Vacuum Shelf</div> <div>Customer</div>						Weight	N/A		Date	17-02-2017
						Mail material:	ST37			
						Drawing No.				1504V001
Design	BJR	Check		Rev. No.	0					
Confidential: Property of Gertsen & Olufsen A/S. Bioreactor. Not to be handed over to, copied or used by third party. Reproduction of contents to be authorized by G&O Maritimegroup.										

Latest Update: 03-03-2017