

Appendix F - Equipment List - Atlantic Institution

Item	Equipment Code-Code d'équipement	Institution-Établissement	Location-Localisation	Manufacturer-Nom du fabricant	Model-Modèle	Serial Number-Numéro de série
AMALGAMATOR-DENTAL	231-24502	Atlantic Institution	D308	Caulk		
ASPIRATOR	231-004505	Atlantic Institution	D303	Advantage	L190-GR	NONE
ASPIRATOR	231-004510	Atlantic Institution	Unit 3		L190-GR	NONE
BATH-WATER	231-308	Atlantic Institution		Eagle		
BED-ELECTRIC	231-003076	Atlantic Institution	D2017 - Dry Cell	Hill Room	P8400D3867	E086AL6139
Braun Thermoscan		Atlantic Institution	Unit 5	Cardinal Health		
CAMERA-CONTROLLER-LIGHT-SOURCE		Atlantic Institution	D112	AMD Global Telemedicine		
CAMERA-CONTROLLER-LIGHT-SOURCE	BMD-7999	Atlantic Institution	D112	Sometech Inc.	DCSM (XGA)	DCSM-11C-030
CAMERA-VIDEO	BMD-8008	Atlantic Institution	D112	AMD Global Telemedicine	AMD-2500	
CAMERA-VIDEO	BMD-8004	Atlantic Institution	D112	Sometech Inc.	DCS-P101	DCSN-11C-118
CAMERA-VIDEO	BMD-8003	Atlantic Institution	D112	Sometech Inc	DCS-104 (NTSC)	DCSM-118-059
Cane - J-Hook 250LB 30-30IN Push Button		Atlantic Institution	Health Care	Cardinal Health		
CENTRIFUGE		Atlantic Institution	Pharmacy	Unico	Power Spln FX	250-520-02
CHAIR-DENTAL	231-003832	Atlantic Institution	D308	DENTECH	ULTIMA2	CM2315-3-K1
CHAIR-OPHTHALMIC	231-004576	Atlantic Institution	D104	TOPCON		
CLEANER-ULTRASONIC	231-004519	Atlantic Institution	D308	L & R	T14B	Q87122194T
CONSOLE-XRAY-UNIT	BMD-8011	Atlantic Institution	D304	GE	2259976	
CRAP Mask, Full face Simplus Medium		Atlantic Institution	Unit 5	Quality Respiratory Care		
CRAP, Aisense 10 Auto		Atlantic Institution	Unit 5	Quality Respiratory Care		
Crutch 300LB 62-70 Adult		Atlantic Institution	Health Care	Cardinal Health		
Crutch 300LB 62-70 Adult		Atlantic Institution	Health Care	Cardinal Health		
Crutch Axil 300LB 70-80		Atlantic Institution	Health Care	Cardinal Health		
DENTAL-UNIT-ULTRASONIC	BMD-7997	Atlantic Institution	Room D308	South East- Henry Schein		
DOPPLER-FETAL-HR	BMD 1456	Atlantic Institution			4070	
ELECTROCARDIOGRAPH	231-003476	Atlantic Institution	D303	Welch Allyn Canada Ltd.	CD200	20004343
ELECTROCARDIOGRAPH & TROLLY		Atlantic Institution	Treatment Room	Cardinal Health	MAC 2000	
IMAGER-COMPUTED-RADIOGRAPHY	23100-004522	Atlantic Institution	D304	Agfa	CR 30-X	
INCUBATOR-BIOLOGICAL	231-3335	Atlantic Institution	D306	3M Attest	305807	208585
KEROMETER	231-284	Atlantic Institution	D104	Bausch & Lomb	B & L 71-21-35	18256
LATHE-DENTAL	231-000312	Atlantic Institution	D308	BDMI		816-22-1287
LENSOMETER	231-000281	Atlantic Institution	D104		322056	322056
Light with mobile stand		Atlantic Institution	Unit 5	Cardinal Health		
LIGHT-EXAMINING	231-0001697	Atlantic Institution	D104	Hein		
LIGHT-EXAMINING	231-000273	Atlantic Institution	D303	Heraeus	55002286	100005
LIGHT-ULTRAVIOLET	BMD-8002	Atlantic Institution		DSP	RADII PLUS	
MICROSCOPE	231-000284	Atlantic Institution	D104	Amtek	CP 40	105072
MONITOR-NIBP-SPO2-TEMP	231-003761	Atlantic Institution		General Electric Mesystems Limited	2009284001	AAW080502335
MONITOR-NIBP-SPO2-TEMP	231-004512	Atlantic Institution	Unit 3	Dinamap	V100-PBXA-AXAA-DX	SH6111400455A
MONITOR-NIBP-SPO2-TEMP	23100-004170	Atlantic Institution	Unit 5	Welch Allyn Canada Ltd.		20140201184
MONITOR-NIBP-SPO2-TEMP	23100-004172	Atlantic Institution	Unit 5	Welch Allyn Canada Ltd.		20140201183
MONITOR-SPO2	BMD-5207	Atlantic Institution		Nonin	8500	500257365
MONITOR-SPO2	231-004527	Atlantic Institution	Unit 3	Nonin	8500	500714515
MONITOR-SPO2	23100 004169	Atlantic Institution	Unit 5	NONIN	SPECTRO2	4047289
MONITOR-SPO2		Atlantic Institution	Unit 5	Nonin	8000AA	500716022

Appendix F - Equipment List - Atlantic Institution

MONITOR-SPO2	231-004513	Atlantic Institution	Unit 3		8500	500327436
MONITOR-TEMP-TYMPANIC	BMD-7993	Atlantic Institution		Welch Allyn	Braun 6021	
MONITOR-TEMP-TYMPANIC	BMD-7993	Atlantic Institution			Braun 6021	
NX PC Platform	23100-004230	Atlantic Institution	D-304-X-Ray Room	Hewlett Packard	NX PC Platform	
OPHTHALMIC-INSTRUMENT-STAND	231-004504	Atlantic Institution	D104	Topcon (Koenigkramer)	IS2000	351541
OPHTHALMOSCOPE	BMD-5210	Atlantic Institution			767	
OPHTHALMOSCOPE	23100-004171	Atlantic Institution	Unit 5	Welch Allyn Canada Ltd.		05259-M
Oximeter Digital Finger		Atlantic Institution	Unit 5	Cardinal Health		
Oximeter Digital Finger		Atlantic Institution	Unit 3	Cardinal Health		
Oximeter Digital Finger		Atlantic Institution	Health Care	Cardinal Health		
Oxygen cylinder holder for stretcher		Atlantic Institution	Unit 5	FERNO		
Oxygen cylinder holder for stretcher		Atlantic Institution	Health Care	FERNO		
PHOROPTER	231-001698	Atlantic Institution	D104	Reichert	11625 AR	10806-2
PRINTER-XRAY	23100-004521	Atlantic Institution	D304	AGFA	Drystar 5300	
PUMP-COMPRESSOR	BMD-11771	Atlantic Institution			31DB0000	
PUMP-COMPRESSOR	BMD-1449	Atlantic Institution			3410	
PUMP-COMPRESSOR	BMD-7880	Atlantic Institution	D303	Nedquip	RT500	170682
PUMP-INFUSION	231-001972	Atlantic Institution	D303	Abbot Laboratories	XL	11203144
PUMP-INFUSION	231-004501	Atlantic Institution	D303	Hospira Worldwide Inc.	123910403	13898447
REFRIGERATION-UNIT	200-004118	Atlantic Institution	Pharmacy		428730	EW3317486
REFRIGERATION-UNIT	CSC-004229	Atlantic Institution		Danby	2306	
REFRIGERATOR-BIOLOGICAL	231-004141	Atlantic Institution	Pharmacy	Panasonic	SR-L6111W-PA	120700476
REFRIGERATOR-BIOLOGICAL	231-003427	Atlantic Institution	Pharmacy	Follett	128934	889235-35605
REFRIGERATOR-BIOLOGICAL	23100-004141	Atlantic Institution	Pharmacy	Panasonic	SR-L6111W-PA	120700476
Scale Digital Eye Level		Atlantic Institution	Unit 5	Cardinal Health		
Scale Digital Eye Level		Atlantic Institution	Unit 5	Cardinal Health		
Scale Digital Eye Level 500/LB/227KG Capacity		Atlantic Institution	Health Care	Cardinal Health		
SCALE-PATIENT	BMD-7952	Atlantic Institution				
SCALE-PATIENT	BMD-7897	Atlantic Institution			SECA	
SPHYGMOMANOMETER-ANALOG		Atlantic Institution	Clinic 5	Welch-Allyn	DS66	
SPHYGMOMANOMETER-ANALOG		Atlantic Institution	Clinic 1	Welch-Allyn	DS66	
SPHYGMOMANOMETER-ANALOG		Atlantic Institution	Clinic 1 (spare)	Welch-Allyn	DS66	
Stand Mobile Spot/LX1 Basket		Atlantic Institution	Unit 5	Cardinal Health		
Standard Quality Wide Touch Monitor	23100-004231	Atlantic Institution	D-304-X-Ray Room	Hewlett Packard (Canada) CO	HP Monitor	
STERILIZER-STEAM-TABLE-TOP	23100-004140	Atlantic Institution	D308	Statim	Statim 2000 G4	101212J00026
STERILIZER-STEAM-TABLE-TOP	231-001988	Atlantic Institution	D308	STATim	1102	228426A
STETHESCOPE-ELECTRONIC	BMD-8009	Atlantic Institution	D112	AMD Global Medecine		21151
STRETCHER	231-3270	Atlantic Institution	Hall	FERNO	287857	785076
STRETCHER	231- 003579	Atlantic Institution	D303	Stryker	1007	
STRETCHER	23001-004142	Atlantic Institution	Halfway	FERNO		238-1493-018
STRETCHER		Atlantic Institution	Basement-For use in Unit 5	FERNO	287784	888113616
STRETCHER	23100-004142	Atlantic Institution	Halfway	Ferno	35X proflex	12N232511
STRETCHER		Atlantic Institution	Unit 5	Stryker	747	1502033580
TABEL-EXAMINING	23100-004036	Atlantic Institution	D112	Global Medco Ltd.	9105-30	
TABEL-EXAMINING	23100-004035	Atlantic Institution	Unit 3	Global Medco Ltd.	9105-30	
TABLE-EXAMINING	231-002886	Atlantic Institution	NOT IN USE - In Basement(Tunnel)		056-605	2702002001
TABLE-PHYSIO-TREATMENT	231-004508	Atlantic Institution	D304	Centry Medical Dist	TRT600	2245

Appendix F - Equipment List - Atlantic Institution

TABLE-XRAY	23100-04516	Atlantic Institution	D304	GE		
Thermometer Caretemp Touchfree		Atlantic Institution	Unit 5	Cardinal Health		
Thermometer Caretemp Touchfree		Atlantic Institution	Unit 3	Cardinal Health		
Thermometer Caretemp Touchfree		Atlantic Institution	Health Care	Cardinal Health		
TONOMETER	231-000285	Atlantic Institution	D104	Topcon (Koenigkramer)	SL4E	643280
TRACTION-UNIT	231-004509	Atlantic Institution	D304	Triton-Chattanooga	MP1	3358
ULTRASOUND-THERAPEUTIC	231-004506	Atlantic Institution	D304	Intelect Advanced	2762CC	4879
WARMER-BLANKET-SOLUTION	231-3582	Atlantic Institution	D306	Pedigo	398442	545085-00
WHEELCHAIR	BMD-7879	Atlantic Institution			700-631	
WHEELCHAIR	BMD-7878	Atlantic Institution			924XL	
XRAY-UNIT	231-289	Atlantic Institution	D304		260-076-R	1462
XRAY-UNIT	23100-004515	Atlantic Institution	D304	Ge	AL01F	02246
XRAY-UNIT	BMD-8016	Atlantic Institution	D304	GE	2269646	
XRAY-UNIT	231-003766	Atlantic Institution	D308	Sinclair Dental	BELRAY 096	EX08K0769

Appendix G - Equipment List - Dorchester Institution (Medium)

Item	Equipment Code-Code d'équipement	Institution-Établissement	Location-Localisation	Manufacturer-Nom du fabricant	Model-Modèle	Serial Number-Numéro de série
AMALGAMATOR-DENTAL	220-004041	Dorchester Institution	Dentist's Office	The LD Caulk Comp.	VM-C	54154
ASPIRATOR	220-009018	Dorchester Institution	Treatment room 1	Contemporary Products	6260	AS21745
ASPIRATOR		Dorchester Institution	Health Services		L190-GR	0813
ASPIRATOR		Dorchester Institution	Clean Utility Room	BMD-10561	ALLUED L190-G4	714
ASPIRATOR	BMD-8064	Dorchester Institution	Trauma Room	Gomco by Allied	405	20080331005
Assembly cam-14		Dorchester Institution	Treatment room 1	GE Medical Systems	CAM-14	RNC11015386PA
Auto refractometer Cradle		Dorchester Institution	Optometry Room	NIDEK	AR20ST1	130493
Auto refractometer Hand Held		Dorchester Institution	Optometry Room	NIDEK	AR20	230454
BATH-WATER		Dorchester Institution	Dentist's Office	Whip Mix		07110195
BED-ELECTRIC (BARIATRIC)HOSPITAL	BMD-16319	Dorchester Institution	REGIONAL HOSPITAL		153008V	
BED-ELECTRIC (HOSPITAL BED)	BMD-10489	Dorchester Institution	REGIONAL HOSPITAL	BMD-10489	HILL-ROM P870	G017AJ
BED-ELECTRIC (HOSPITAL BED)	BMD-10490	Dorchester Institution	REGIONAL HOSPITAL	BMD-10490	HILL-ROM P870	G349AJ
BED-ELECTRIC (HOSPITAL BED)	BMD-10491	Dorchester Institution	REGIONAL HOSPITAL	BMD-10491	HILL-ROM P870	H219AJ
BED-ELECTRIC (HOSPITAL BED)	BMD-10492	Dorchester Institution	REGIONAL HOSPITAL	BMD-10492	HILL-ROM P870	H219AJ
CAMERA-CONTROLLER-LIGHT-SOURCE	220-009022	Dorchester Institution	Nursing supervisor's office	Sometec	DCSM(X6A)	DCSM-10G-167
Camscope wand (angled)	220-009022	Dorchester Institution	Nursing supervisor's office	Sometec	DCS-104T (NTSC)	DCSN-10E-217
Camscope wand (straigh)	220-009022	Dorchester Institution	Nursing supervisor's office	Sometec	P101	DCSN-11C-116
CART-DENTAL	220-009024	Dorchester Institution	Treatment room 1	-	447192	C912329
CART-DENTAL		Dorchester Institution	Dentist's office	Belmont		12057301
CART-EMERGENCY	220-009042	Dorchester Institution	Storage room 3rd floor	Harloff Co.		46269
CART-VIDEO-MONITOR	BMD-11663	Dorchester Institution	IDN Office	Polycom	PPCE-1208P	N091215032
CENTRIFUGE	BMD-10562	Dorchester Institution	Treatment room 1	Thermo IEC	CENTRA CL-2	42613405
CHAIR RECLINE TILT	BMD-10563	Dorchester Institution	CLEAN EQUIPMENT STORAGE ROOM	BMD-10563	LUMEX FR574G 454-PRO-LAB REV	7170760 / HSP321
CHAIR RECLINE TILT	BMD-10552	Dorchester Institution	Trauma Room	BMD-10552	WINCO 574	574100
CHAIR RECLINE TILT	BMD-10553	Dorchester Institution	IDN Office	BMD-10553	WINCO 570	570810
CHAIR-DENTAL	220-009023	Dorchester Institution	Dentist's office	Takara Belmont	447193	012057301
CHAIR-OPHTHALMIC	220-07672	Dorchester Institution			FX920	93010707007
CHAIR-OPHTHALMIC	220-009031	Dorchester Institution	X-ray room	Power Tilt	FX920L	93010707007
CLEANER-ULTRASONIC	220-000201	Dorchester Institution	Dentist's lab	L & R Manufacturing	T-14B	150
CLEANER-UITRASONIC	220-07243	Dorchester Institution	Dentist's Office	Densply	cAVITRON SPS	119-27469
COLLIMATOR-XRAY	220-009038	Dorchester Institution	X-ray room	GE Hualun Medical System	2379827	53339H18
COLLIMATOR-XRAY	220-009027	Dorchester Institution	Dentist's office	Belmont	097-H	EH10J0510
COMPRESSION-UNIT-SEQUENTIAL		Dorchester Institution	Dentist's office	Dominic Hunter	DHDRD10	3701300004
DISPENCER-ICE		Dorchester Institution	B-6 End of hallway	Whirlpool	G15NDXTB1	EW1501128
DRILL-CONSOLE-SURGICAL	220-009024	Dorchester Institution	Dentist's office	-	-	--
ELECTROCARDIOGRAPH	220-009002	Dorchester Institution	Treatment room 1	GE Medical Systems	MAC5500	11-41583
Frame Heater		Dorchester Institution	Optometry Room	GFC	A600	682
HYFRECATOR	220-000209	Dorchester Institution	Treatment Room 1	Birtcher		2157354
IMAGER-COMPUTED-RADIOGRAPHT	220-009034 + 220-07878	Dorchester Institution	X-ray room	AGFA	CR30-X	2676
INTERFERENTIAL-ULTRASOUND	220-009015	Dorchester Institution	Treatment room 2	Chattanooga Canada	Intellect Adv	4978
KERATOMETER	220-006033	Dorchester Institution	Optometry Room	Marco		21226
LAMP-SLIT	221-006503	Dorchester Institution	Dr room	Top Con PS-20A	SL-2E	6211049
LAMP-SLIT	220-009040	Dorchester Institution	Optometry Room	HAAG Streit	BQ900	10654
LATHE-DENTAL	BMD-8061	Dorchester Institution	Dentist's Office	Red Wing	16B	243/10/11/044
LENSOMETER	220-009037	Dorchester Institution	X-ray room	Marco	201	M2474

Appendix G - Equipment List - Dorchester Institution (Medium)

LENSOMETER		Dorchester Institution	Optometry Room			
LIFT-PATIENT-FLOOR	BMD-8067	Dorchester Institution	Range	Steady Aid	PWR TSA3510X	7159
LIFT-PATIENT-FLOOR		Dorchester Institution	Range	BHM Medical	Ergolift	ERL2746
LIFT-PATIENT-FLOOR	BMD-7407	Dorchester Institution	Range		0102-8851	96511054100
LIFT-PATIENT-FLOOR	7989	Dorchester Institution	Range	Steady Aid	NA016D	HARN
LIFT-PATIENT-FLOOR	BMD-10539	Dorchester Institution	RANGE	BMD-10539	INVACARE RPL450-1	14KM01
LIFT-PATIENT-FLOOR-BATH	BMD-8066	Dorchester Institution	Range	TR EQUIPMENT	9650-0542020	1210-20929
LIGHT-EXAMINING	220-009025	Dorchester Institution	Dentist's office	Belmont	ALS01T	T040461
LIGHT-EXAMINING	220-009014	Dorchester Institution	Treatment room 2	Welch Allyn Canada Ltd	44200	--
LIGHT-EXAMINING		Dorchester Institution	DOCTOR'S CLINIC	BMD-10550	LUXO KFM-17213	E24877
LIGHT-HEAD	220-009033	Dorchester Institution	X-ray room	Heine	M Pack	2157354
LIGHT-HEAD	220-009033	Dorchester Institution	X-ray room	Heine	--	1130545
LIGHT-ULTRAVIOLET		Dorchester Institution	Dentist's office	Coltene/Whaledent	Coltoflux 75	060128666
LIGHT-ULTRAVIOLET	BMD-10383	Dorchester Institution	Dentist's Office	Elipar	S10	939112017246
LIGHT-ULTRAVIOLET	BMD-11674	Dorchester Institution	Treatment room 1	Luxo	0316040	89705
MONITOR		Dorchester Institution	Nursing supervisor's office	Samsung	LN26B460B2D	
MONITOR-NIBP		Dorchester Institution	NURSES CLINIC	BMD-10551	OMRON HEM-907XL	2014030
MONITOR-NIBP		Dorchester Institution		BMD-10422		
MONITOR-NIBP-SPO2-TEMP	220-009011	Dorchester Institution	Treatment room 2	Welch Allyn	Spot Vital Signs Ltd	20110201548
MONITOR-NIBP-SPO2-TEMP	220-009181	Dorchester Institution	Storage room 3rd floor	Critikon	Dynamap Compact BP monitor	ref. 117207
MONITOR-NIBP-SPO2-TEMP	220-009004	Dorchester Institution	Treatment room 1	GE Medical Systems	Dynamap DPC400N-CE	2019205-001
MONITOR-SPO2	BMD-10423	Dorchester Institution	To be confirmed			
MONITOR-SPO2	BMD-7413	Dorchester Institution	To be confirmed			
MONITOR-SPO2	BMD-10420	Dorchester Institution	To be confirmed			
MONITOR-SPO2		Dorchester Institution	EMERGENCY BAG-HC	BMD-10418	MASIMO	526205
MONITOR-SPO2		Dorchester Institution	EMERGENCY BAG-FRONT ENTRANCE	BMD-10413	WELCH ALLYN WW1020EN	4025973
MONITOR-SPO2	220-009003	Dorchester Institution	Treatment room 2	Masimo Corporation	RADS	532646
MONITOR-SPO2	220-009008	Dorchester Institution	Treatment room 1	Masimo Corporation	RADS	530779
MONITOR-TEMP		Dorchester Institution	Clean Utility Room	BMD-10419		
MONITOR-TEMP		Dorchester Institution	Emergency Bag HC	BMD-10421	690	
MONITOR-TEMP	BMD-7414	Dorchester Institution	IDN Office	Welch Allyn	690	08272340
MONITOR-TEMP	BMD-11676	Dorchester Institution		Welch Allyn	690	11110835
MONITOR-TEMP	BMD-8058	Dorchester Institution	To be confirmed	Welch Allyn		
MONITOR-TEMP		Dorchester Institution	Clean Utility Room	Welch Allyn BMD-10544	690	
MONITOR-TEMP	BMD-8059	Dorchester Institution	Trauma Room	Welch Allyn	690	12512584
MONITOR-TEMP	220-009010	Dorchester Institution		Welch Allyn	690 Thermometer	08047847
MONITOR-TEMP	BMD-11675	Dorchester Institution	To be confirmed	Welch Allyn		
MONITOR-TEMP		Dorchester Institution	EMERGENCY BAG-FRONT ENTRANCE	BMD-10412	WELCH ALLYN 690	124804
MONITOR-VITALS		Dorchester Institution	TBD			
MONITOR-VITALS		Dorchester Institution	TBD			
MONITOR-VITALS		Dorchester Institution	TBD			
MONITOR-VITALS		Dorchester Institution	TBD			
OPHTHALMIC-INSTRUMENT	220-002017	Dorchester Institution	Optometry Room	Masimo		500255
OPHTHALMIC-INSTRUMENT-STAND	220-009032	Dorchester Institution	X-ray room	Refiance Medical Products	7900	79010702014
OPHTHALMIC-INSTRUMENT-STAND	220-07671	Dorchester Institution	Optometry Room		7900	79010702014
OPHTHALMOSCOPE	220-009021	Dorchester Institution	Treatment Room 2	Welch Allyn		
OPHTHALMOSCOPE	220-009006	Dorchester Institution	Treatment Room 1	Welch Allyn		

Appendix G - Equipment List - Dorchester Institution (Medium)

Ophthalmoscope		Dorchester Institution	Trauma Room	Welch Allyn	GS 777 WA 7791-MX	
Optometry equipment	220-09170	Dorchester Institution	Optometry Room	NIDEK		130493
PHOROPTER	220-0004032	Dorchester Institution	X-ray room	Marco	RT-I	10975
PROJECTOR-OPHTHALMIC	220-009030 + 220-006025 =	Dorchester Institution	X-ray room	Marco	--	22143
PROJECTOR-OPHTHALMIC	220-009030	Dorchester Institution	X-ray room	Marco	50	M22143
PUMP-COMPRESSOR		Dorchester Institution	Dentist's office	Belmont	4000	20213
PUMP-COMPRESSOR	220-002007	Dorchester Institution	Storage 1st Floor	Schuco	5711102	
PUMP-COMPRESSOR	BMD-11655	Dorchester Institution	Trauma Room	AirLife	RT500	169542
PUMP-COMPRESSOR		Dorchester Institution		BMD-11656	CAREFUSION (PULMONETICS) RTS	169548
PUMP-INFUSION	BMD-8060	Dorchester Institution	Trauma Room	Hospira	1239-104-07	17446358
PUMP-INFUSION	220-009001	Dorchester Institution	Treatment room 1	Hospira	1239-104-03	13898441
PUMP-INFUSION	220-009028	Dorchester Institution	Medical supply room 1st floor	Hospira	1239-104-03	13898345
PUMP-INFUSION	220-006718	Dorchester Institution	Medical supply room 1st floor	Abbott Laboratories	PLUMXL 11555-04-01	11204399
PUMP-VACUUM		Dorchester Institution	Hall outside Dental Office	VacStar	V520	208990
REFRIGERATOR-BIOLOGICAL	220-07993	Dorchester Institution	IDN Office		SAN-SRL6111W	120200257
REFRIGERATOR-BIOLOGICAL	220-07992	Dorchester Institution	Secondary Pharmacy		SAN-MPRS14	12010029
REFRIGERATOR-BIOLOGICAL		Dorchester Institution		BMD-11657		
REFRIGERATOR-BIOLOGICAL		Dorchester Institution		BMD-11883	LG ELECTRONICS CANADA	GC-151RCW
REFRIGERATOR-BIOLOGICAL		Dorchester Institution		BMD-11881	HAIER	NSA36L
SAW-CAST	220-009016	Dorchester Institution	Treatment room 2	Stryker	8127	E34984
SCALE-PATIENT	220-009017	Dorchester Institution	Treatment Room 1	Pelstar	500KL	5000005579
SCALE-PATIENT	BMD-16324	Dorchester Institution	Trauma Room		SECA 703 1321993	5703160
STERILIZER-STEAM-TABLE-TOP	220-006702	Dorchester Institution	Dentist's lab	SciCan	1102	228756A
STERILIZER-STEAM-TABLE-TOP		Dorchester Institution	Dentist's Office	Midmark	M9 UltraClave M9-022	V546354
STERILIZER-STEAM-TABLE-TOP	BMD-8070	Dorchester Institution	Dentist's Office	SciScan	Statim 2000S - 01-121105	101212J00019
STRETCHER	220-006439	Dorchester Institution	Treatment room 1	Midmark	520	
STRETCHER	220-009087	Dorchester Institution		Stryker	0747-000	1105034705
STRETCHER	220-009088	Dorchester Institution		Stryker	0747-000	1105034706
STRETCHER		Dorchester Institution		BMD-11882	MEDITECH 500	TMPR13
STRETCHER		Dorchester Institution		BMD-10540	FERNO CANADA 35X	10N-14
TABLE-EXAMINING		Dorchester Institution		BMD-10556	MIDMARK 204-001	V166734
TABLE-PHYSIO-TREATMENT	220-009013	Dorchester Institution	Treatment room 2	Chattanooga Canada	Triton DTS TRT-600	2214
TABLE-XRAY	220-009038	Dorchester Institution	X-ray room	GE Hualun Medical System	2259988	S0508H19
TONOMETER		Dorchester Institution	X-ray room	Clement Clarke	Perkins tonometer MK2	1999T2073
TRACTION-UNIT		Dorchester Institution	Treatment room 2	Chattanooga Canada	2841/3336	
Trolley		Dorchester Institution	Treatment room 1	GE Medical Systems	MAC 5000	SCD11051932PA
Tub - Hospital		Dorchester Institution	REGIONAL HOSPITAL			
ULTRASOUND-DIAGNOSTIC-BLADDER	220-009012	Dorchester Institution	Treatment room 2	Bladder Scan	0570-0188	PA011179
ULTRASOUND-THERAPEUTIC		Dorchester Institution	Treatment room 2	Chattanooga Canada	27335	37018
VENTILATOR-BIPAP-CRAP		Dorchester Institution	Range	Phillips		C04772616
VENTILATOR-BIPAP-CRAP-MARLOW		Dorchester Institution	RANGE	BMD-10541	FISHER PAYKEL HC2346H0	1005140
VENTILATOR-BIPAP-CRAP-MARLOW		Dorchester Institution	REGIONAL HOSPITAL	Resmed Corporation	33039	20091
WARMER-BLANKET-SOLUTION		Dorchester Institution	Clean Utility Room	Pedigo	P-20105	545088-000
WASHER-DISINFECTOR		Dorchester Institution	Hall	FRANKE		
WHEELCHAIR - ELECTRIC		Dorchester Institution	REGIONAL HOSPITAL	Quantum Edge	WC Q6EDGE HD	JC703
WHEELCHAIR - ELECTRIC		Dorchester Institution	REGIONAL HOSPITAL	Chignecto		06-222
X-RAY - RAD	220-009038	Dorchester Institution	X-ray room	GE Hualun Medical System	2259981	95296-7U

Appendix G - Equipment List - Dorchester Institution (Medium)

XRAY-BUCKY	220-009035	Dorchester Institution	X-ray room	GE Hualun Medical System	600-0301	50997HL6
XRAY-CONTROLLER	220-009027	Dorchester Institution	Dentist's office	Belmont	097-R	ER10/0510
XRAY-PRINTER	220-009039	Dorchester Institution	X-ray room	Sharp AGFA	Drystar 5300	
XRAY-PROCESSOR-FILM	220-009026	Dorchester Institution	Dentist's office	Deil	Studio Hybrid 140g	18HV8M1
XRAY-UNIT	220-009036	Dorchester Institution	X-ray room	GE Medical System	2320993	07115877
XRAY-UNIT	220-009038	Dorchester Institution	X-ray room	GE Hualun Medical System	2259976	50529HL7
XRAY-UNIT	220-009027	Dorchester Institution	Dentist's office	Belmont	097-CM	EX10/0141
XRAY-UNIT	220-009038 + 220-07879	Dorchester Institution	X-ray room	GE Hualun Medical System	2259973-2	50567HL7 + 5795
XRAY-UNIT		Dorchester Institution	XRAY ROOM	BMD-11669	GE MEDSYSTEMS LTD 2269646	101445

Appendix H - Equipment List - Dorchester Institution (Minimum)

Item	Equipment Code-Code d'équipement	Institution-Établissement	Location-Localisation	Manufacturer-Nom du fabricant	Model-Modèle	Serial Number-Numéro de série
AMALGAMATOR-DENTAL	BMD-10361	Dorchester Minimum	Dentists rm			G-10-1-4-V
AMALGAMATOR-DENTAL		Dorchester Minimum	DENTAL OFFICE	BMD-10373	DENTSPLY CANADA LTD	VM-D (VARL-MIX 11
ASPIRATOR	221-006303	Dorchester Minimum	Dr room	Gomco	300	L10395
ASPIRATOR	220-09121	Dorchester Minimum	NURSES' CLINIC	BMD-10560	LIFESUPPORT PRODUCT	L190-GR
Camscope angled wand	221-002023	Dorchester Minimum	Dr room	Sometech	DCS-104T (NTSC)	DSCN-11B-061
Camscope receiver	221-002023	Dorchester Minimum	Dr room	Sometech	DCSM (XGA)	DCSM-11C-029
Camscope straight wand	221-002023	Dorchester Minimum	Dr room	Sometech	P101	DCSM-11C-115
CART-DENTAL	221-002027	Dorchester Minimum	Dentists rm	Aidec	PACI	C912329
CART-DENTAL	220-09169	Dorchester Minimum	DENTAL ROOM	DEXIS	DEXIS LLC-PLUS102C	13463
CART-VIDEO-MONITOR	221-001623	Dorchester Minimum	Dr room	Polycorn	PPCE-1208P	L100614525
Cauk Dentsply	221-0005000	Dorchester Minimum	Dentists rm	Dentsply International	VM-D	71970
CENTRIFUGE	221-002014	Dorchester Minimum	Nurses Office	Thermo Scientific	CL-2	42619248
CHAIR-DENTAL	BMD-7422	Dorchester Minimum	Dentist rm	Belmont	093G	36911
CHAIR-DENTAL	BMD-7425	Dorchester Minimum	Dentist rm	Belmont	092Q	36910
CHAIR-DENTAL	221-009078	Dorchester Minimum	Dentists rm	Belmont	X-Caliber	AN09A0073
CHAIR-DENTAL	BMD-10354	Dorchester Minimum	Dentists rm	Belmont	BDS0050	13344
CHAIR-DENTAL		Dorchester Minimum	DENTAL ROOM	BMD-1532		SDP1C
CHAIR-EXAMINATION	221-002015	Dorchester Minimum	Nurses Office	Winco	574	574100103
COMPRESSOR-UNIT-SEQUENTIAL	221-009082	Dorchester Minimum	storage rm	Midmark	2770CE50-660A	
Dental X-ray Arm	221-009079	Dorchester Minimum	Dentist rm	Belmont	096-A	EA08K0779
Dental X-Ray Controler	221-009079	Dorchester Minimum	Dentists rm	Belmont	96	EX08K0779
DENTAL-UNIT-ULTRASONIC	221-000499	Dorchester Minimum	Dentists rm	Dental Equipment	2002	103-05362
DOPPLER-FETAL-HR		Dorchester Minimum	DOCTOR'S OFFICE	BMD-10565	CARDIO BEAT	CB0879
ELECTROCARDIOGRAPH		Dorchester Minimum	Doctor's Clinic	Welch Allyn	(01)00732094237740	(21) 100032813917
ELECTROCARDIOGRAPH Cart		Dorchester Minimum	Doctor's Clinic	Welch Allyn	105341 / Lot 17269	
LASER-ULTRASOUND	221-002007	Dorchester Minimum	Dr room	Dynatronics Research Corp.	D850PLUS	860726
LATHE-DENTAL	221-000535	Dorchester Minimum	Dentists rm	Red-wing	26A	816-227 84
LIGHT MAGNIFYING		Dorchester Minimum	NP'S OFFICE	BMD-10564	LUXO	316040
LIGHT-EXAMINING	BMD-10375	Dorchester Minimum	Dentist rm	Belmont	AL501TSCR3	AV09A0040
LIGHT-EXAMINING	221-002009	Dorchester Minimum	Dr room	Welch Allyn Canada Ltd.	LS200	ref # 74227
LIGHT-EXAMINING	221-002008	Dorchester Minimum	Dr room	Welch Allyn Canada Ltd.	48100	
LIGHT-EXAMINING	BMD-10555	Dorchester Minimum	Nurse practitioner's office	BMD-10555	WELSH ALLYN G5600	150710
LIGHT-ULTRAVIOLET	BMD-10369	Dorchester Minimum	Dentists rm	Coitene/Whaledent INC.	C7970100115	110507017
LIGHT-ULTRAVIOLET		Dorchester Minimum	Dentists clinic	3M ESPE	Elipar 510	
LIGHT-ULTRAVIOLET	BMD-8034	Dorchester Minimum	Dentist rm	3M ESPE	Elipar 510	939112017014
MONITOR-NIBP	221-002013	Dorchester Minimum	Dr room	Critikon / Johnson & Johnson	100V2	043M0299009
MONITOR-NIBP	221-002026	Dorchester Minimum	Nurses office	Critikon/Johnson & Johnson	DCP400N-CE	2019205-001
MONITOR-NIBP		Dorchester Minimum	DOCTOR'S CLINIC	BMD-10557	OMRONHEM-907XL	201508
MONITOR-SPO2	221-002017	Dorchester Minimum	Nurses Office	Masimo Corporation	RAD5	500255
MONITOR-SPO2	BMD-7533	Dorchester Minimum	Dr room	Pulse Oximeter	MS0E09CS0	9060429
MONITOR-SPO2		Dorchester Minimum	EMERGENCY BAG	BMD-5202		
MONITOR-SPO2		Dorchester Minimum		BMD-5201	MASIMO CORPORATION	RAD5
MONITOR-SPO2		Dorchester Minimum	Doctor's clinic	Smiths Medical	SpectrOW 10 -	2015-
MONITOR-TEMP	221-002018	Dorchester Minimum	Nurses office	Welch Allyn Canada Ltd.	690	4106441

Appendix H - Equipment List - Dorchester Institution (Minimum)

MONITOR-TEMP	221-002010	Dorchester Minimum	Dr room	Welch Allyn Canada Ltd.	690	09361171
MONITOR-TEMP	220-002019	Dorchester Minimum	EMERGENCY BAG	Welch Allyn BMD-1571	690	14310025
MONITOR-TEMP		Dorchester Minimum		BMD-10493	WELCH ALLYN CANADA	690 SURETEMPPLUS
OPHTHALMOSCOPE	221-002018	Dorchester Minimum	Nurses office	Welch Allyn Canada Ltd.	74710	92192
OPHTHALMOSCOPE	BMD-1409	Dorchester Minimum	Dr room	Welch Allyn Canada Ltd.	767	5290
OPHTHALMOSCOPE	BMD-10554	Dorchester Minimum	Nurse practitioner's office	BMD-10554	WELSH ALLYN 777	
PUMP COMPRESSOR		Dorchester Minimum	Storage rm	AirLife	RT500	169547
PUMP COMPRESSOR	BMD-8051	Dorchester Minimum	Storage rm	Invacare	IRC1196	99H328019
PUMP COMPRESSOR	BMD-8049	Dorchester Minimum	Storage rm	AirLife	RT500	169547
PUMP-COMPRESSOR	221-006303	Dorchester Minimum	Dr. room	Schuco	5711102	188532
PUMP-COMPRESSOR		Dorchester Minimum		BMD-11672	SCHUCO	571102 SCHUCO - M
PUMP-COMPRESSOR	BMD-8050	Dorchester Minimum		BMD-8050	INVACARE CANADA	IRC1196-PASSPORT
PUMP-COMPRESSOR		Dorchester Minimum		BMD-8051	INVACARE CANADA	IRC1196-PASSPORT
PUMP-COMPRESSOR		Dorchester Minimum		BMD-10558	AIR TECHNIQUES	VS40VL (VACSTAR 40)
PUMP-INFUSION	221-009028	Dorchester Minimum	Nurses office	Plum A+	84095270002	13898345
REFRIGERATOR-BIOLOGICAL	BMD-8035	Dorchester Minimum	Inf Dis Office	Panasonic	SR-L6111W-PA	121101390
SAW-CAST	221-002021	Dorchester Minimum	Dr room	Stryker Canada	840	92111273
SCALE-PATIENT	BMD-16327	Dorchester Minimum	Nurses Clinic	SECA	8761321004	
SCALE-PATIENT	BMD-16327	Dorchester Minimum	Nurses Clinic	Health O Meter	500KL	
SCALE-PATIENT		Dorchester Minimum	HALL	BMD-10559	SECA 7031321993	570316
STERILIZER-STEAM-TABLE-TOP	221-002024	Dorchester Minimum	Dentists rm	L&R	T-148	285 83120
STERILIZER-STEAM-TABLE-TOP	221-002020	Dorchester Minimum	Nurses Office	Midmark	KS20	TAV003996
STERILIZER-STEAM-TABLE-TOP		Dorchester Minimum	DENTIST'S OFFICE	BMD-8038	SCICAN	STATIM 2000
TABLE-EXAMINING		Dorchester Minimum	Nurse practitioner's office	BMD-10556	RITTER / MIDMARK 201-	V16673
TABLE-PHYSIO-TREATMENT	221-002005	Dorchester Minimum	Dr room	Cardon Rehabilitation & medical equipm	R28535	039-08011
VENTILATOR BIPAP CRAP		Dorchester Minimum		BMD-10494	RESPIRONICS INC	551P
VENTILATOR BIPAP CRAP		Dorchester Minimum		BMD-10498	RESMED CORP.	S9 ELITE
VIEWER XRAY FILM		Dorchester Minimum		BMD-X045	OLIVE CRAIG	XRAY VIEWER
WARMER-BLANKET-SOLUTION	221-002028	Dorchester Minimum	Medication Room	Pedigo	P-20105	545084-000
XRAY-UNIT	221-009079	Dorchester Minimum	Dentists rm	Belmont	096-H	EH08K0779

Appendix I - Equipment List Regional Headquarters

Item	Equipment Code-Code d'équipement	Institution-Etablissement	Location-Localisation	Manufacturer-Nom du fabricant	Model-Modèle	Serial Number-Numéro de série
CAMERA-CONTROLLER		Regional Headquarters				
CAMERA-CONTROLLER-LIGHT-SOURCE	20100-004579	Regional Headquarters		Sometech	456384	DCSM-11C-015
CAMERA-VIDEO	20100-004578	Regional Headquarters		Sometech	456383	DCSN-088-341
REFRIGERATOR-BIOLOGICAL	20100-004597	Regional Headquarters	Currently in storage (basement @ RHQ)	ESBE Scientific	458353	40Z07019
STETHOSCOPE-ELECTRONIC		Regional Headquarters			AMD-3700	
STETHOSCOPE-ELECTRONIC		Regional Headquarters			AMD-3700	
STETHOSCOPE-ELECTRONIC	20100-004547	Regional Headquarters		AMD	448659	SN 21151
STETHOSCOPE-ELECTRONIC	20100-004548	Regional Headquarters		AMD	448650	SN 22150

Appendix J - Equipment List - Regional Pharmacy

Item	Equipment Code-Code d'équipement	Institution-Établissement	Location-Localisation	Manufacturer-Nom du fabricant	Model-Modèle	Serial Number-Numéro de série
COUNTING MACHINE	201-004553	Regional Pharmacy	Dispensary	McKesson	Parata Accucount	NTEP CC No. 96-116A3P
PACKAGE-SEALER-PHARMACY	201-004554	Regional Pharmacy	Packaging Room	Healthmark	Autogen	AG-009C
PACKAGE-SEALER-PHARMACY	201-004550	Regional Pharmacy	Packaging Room	Healthmark	Autogen	AG-003C
PACKAGE-SEALER-PHARMACY	201-004552	Regional Pharmacy	Packaging Room	Healthmark	MTS -350	350-TC-004
PACKAGE-SEALER-PHARMACY	210-005023	Regional Pharmacy	Shipping / Receiving	Hunt Corporation	Bienfang	MH-00804
PUMP-DOSING-PHARMACY	201-004551	Regional Pharmacy	Methadone Room	Baxa	Baxa	
PUMP-DOSING-PHARMACY	201-004555	Regional Pharmacy	Methadone Room	Baxa	Baxa	
PUMP-DOSING-PHARMACY	201-004775	Regional Pharmacy	Methadone Room	Baxa	Baxa	
PUMP-DOSING-PHARMACY	201-004810	Regional Pharmacy	Methadone Room	Baxa	Baxa	
REFRIGERATOR-BIOLOGICAL	201-005024	Regional Pharmacy	Methadone Room	True Food Services	TS-49 - True	4978832
REFRIGERATOR-BIOLOGICAL	201-005025	Regional Pharmacy	Dispensary	True Food Services	TS-49 - True	4947312
REFRIGERATOR-BIOLOGICAL	201-004595	Regional Pharmacy	Basement	ESBE Scientific	SAN-MPR1411	10020036
SCALE		Regional Pharmacy		PTA	Accu Count	

Appendix K - Equipment List - Shepody Healing Center

Item	Equipment Code-Code d'équipement	Institution-Établissement	Location-Localisation	Manufacturer-Nom du fabricant	Model-Modèle	Serial Number-Numéro de série
BED-ELECTRIC		Shepody Healing Centre	UNIT 5	BMD-10486	ROLOTEC BED	
BED RESTRAINT		Shepody Healing Centre	UNIT 5	Kinesic (Max-Secure)	EN 6000-12RB	
DYNAMOMETER-HAND		Shepody Healing Centre	Gisele's office 4th floor	BMD-10485	JAMAR	
DYNAMOMETER-HAND		Shepody Healing Centre	PSR Unit			
MONITOR NIBP SPO2 TEMP		Shepody Healing Centre	Acute Unit	Welch Allyn Canada Ltd	42NTB	
MONITOR NIBP SPO2 TEMP		Shepody Healing Centre	Acute Unit	Welch Allyn Canada Ltd	42 MTB	201419692
MONITOR NIBP SPO2 TEMP		Shepody Healing Centre	PSR Unit	Welch Allyn Canada Ltd	42NTB	
MONITOR NIBP SPO2 TEMP		Shepody Healing Centre	PSR Unit	Welch Allyn Canada Ltd	42NMTB	201419223
MONITOR SPO2		Shepody Healing Centre	Acute Unit	O-Two Systems International Inc	01MN1000	
MONITOR SPO2		Shepody Healing Centre	Acute Unit	O-Two Systems International Inc	01MN1000	
MONITOR SPO2		Shepody Healing Centre	PSR Unit	O-Two Systems International Inc	01MN1000	
MONITOR TEMP		Shepody Healing Centre	Acute Unit	Welch Allyn Canada Ltd	690	
MONITOR TEMP		Shepody Healing Centre	Acute Unit	Welch Allyn Canada Ltd	690	
MONITOR TEMP		Shepody Healing Centre	PSR Unit	Welch Allyn Canada Ltd	690	
MONITOR TEMP TYMPANIC		Shepody Healing Centre	Acute Unit	Welch Allyn Canada Ltd	PRO4000	
MONITOR-TEMP		Shepody Healing Centre	UNIT 5		WELCH ALLYN	690
MONITOR-TEMP		Shepody Healing Centre	UNIT 6		WELCH ALLYN	690
REFRIGERATOR-BIOLOGICAL		Shepody Healing Centre	Acute Unit	Haier (LG Electronics)		NSA36LACG
REFRIGERATOR-BIOLOGICAL		Shepody Healing Centre	PSR Unit	LG	GC151RCW	
SCALE PATIENT		Shepody Healing Centre	PSR Unit	Health-O-Meter	402KL	4020079413
SCALE-PATIENT		Shepody Healing Centre	UNIT 5		SECA	8761321004
SCALE-PATIENT		Shepody Healing Centre	UNIT 5		HEALTH-O-METER	160KL
SPHYGMOMANOMETER ANALOG		Shepody Healing Centre	Acute Unit	Heine	GAMMA G5	
SPHYGMOMANOMETER-ANALOG		Shepody Healing Centre	UNIT 6	BMD-10542	WELCH ALLYN TYCOS	
SPHYGMOMANOMETER-ANALOG		Shepody Healing Centre	UNIT 6	BMD-10543	WELCH ALLYN	
STRETCHER		Shepody Healing Centre	PSR Unit	Techlem	500	TMR1303029
STRETCHER		Shepody Healing Centre			FERNO PRO FLEX	35X

Appendix L - Equipment List Jamieson CCC

Item	Equipment Code-Code d'équipement	Institution-Établissement	Location-Localisation	Manufacturer-Nom du fabricant	Model-Modèle	Serial Number-Numéro de série
LIGHT-EXAMINING		Jamieson CCC				
MONITOR-TEMP		Jamieson CCC				
OPHTHALMOSCOPE		Jamieson CCC				
SPHYGMOMANOMETER-ANALOG		Jamieson CCC				
TABLE-EXAMINING		Jamieson CCC				

Appendix M - Equipment List - Nova Institution

Item	Equipment Code-Code d'équipement	Institution-Établissement	Location-Localisation	Manufacturer-Nom du fabricant	Model-Modèle	Serial Number-Numéro de série
AMALGAMATOR-DENTAL	250-001441	Nova Institution	Dental suite (pump room)	Henry Schein	Solmetex/Vacstar/Internat	25669/89824/5858009
AMALGAMATOR-DENTAL	BMD-11632	Nova Institution	Dental suite	Silamat	55	
AMALGAMATOR-DENTAL		Nova Institution	Dental suite	Buffalo	No. 1A	
ASPIRATOR	250-001576	Nova Institution	Tx room	Gomco	400S	20071218030
BREAST-PUMP	25000 001752	Nova Institution	HEALTH CARE UNIT		MEDELA	MEDELA
CAMERA-VIDEO-TELEHEALTH		Nova Institution	Tx Room	Dr. Camscope	DCS-104T(NTSC)	DCSN-10E-285
CAMERA-VIDEO-WAND-TELEHEALTH	20100-004578	Nova Institution	Tx Room	Dr. Camscope	DCS-P101	DCSN-08B-341
CART-VIDEO-MONITOR		Nova Institution	Tx Room	Dr. Camscope	DCSM(XGA)	DCSM-10G0169
CHAIR-DENTAL	250-001449	Nova Institution	Dental suite	Xcalibur		94492407
CHAIR-DENTAL		Nova Institution	Dr office	Brewer	none	61609
CHAIR-OPHTHALMIC	25000-003999	Nova Institution	Dr office	Axis Medical	CS-418	CP-PZY-0418AJ
DENTAL-UNIT-ULTRASONIC	250-001454	Nova Institution	Dental suite	Suifon	Prosonic 1000MT	QDGO20710129E
DENTAL-UNIT-ULTRASONIC	250-001444	Nova Institution	Dental suite	Densphy	Cavitron GEN-131	9788
DOPPLER-FETAL-HR	250-001560	Nova Institution	Ax room	Medisonics		FB06949
DOPPLER-FETAL-HR		Nova Institution		Natus	200	ECO1776
DRILL-CONSOLE-SURGICAL	250-001452	Nova Institution	Dental suite	Adec Pac	3420	C888511
ELECTROCARDIOGRAPH	250-001458	Nova Institution	Tx room	Welch Allyn	CP2AS	20004345
ELECTROCARDIOGRAPH		Nova Institution	Tx Room	Welch Allyn	CP 100	10021867
ELECTROCAUTERY-THERMAL	250-001443	Nova Institution	Dental suite	System 8		
Half Shadow Mirror		Nova Institution	Dr office	Axis Medical	Plain white	none
HYDROCOLLATOR-HOT-PACK	250-001118	Nova Institution	Tx room		E-1	137706
HYDROCOLLATOR-HOT-PACK		Nova Institution	Treatment Room	Chattanooga	2402 BMD-9379	T7783C
Integrated Diagnostic System	250-001558	Nova Institution	Ax Room	Welch Allyn		
Integrated Diagnostic System	250-001580	Nova Institution	Tx Room	Welch Allyn		
LAMP-SLIT	25000-003996	Nova Institution	Dr office	Huvitz	HS-5000	5HS0001310002
LCM control module	250-001440	Nova Institution	Dental suite	Sinclair Dental	LUX115V	
LIGHT-EXAMINING	250-001446	Nova Institution	Dental suite	ALT	508T12N	CN07L0139
LIGHT-EXAMINING	250-001578	Nova Institution	Tx room	Luxo		
LIGHT-EXAMINING	250-001566	Nova Institution	Secure Unit		LS-135	
LIGHT-EXAMINING	250-001551	Nova Institution	Ax Room	Welch Allyn	GS300	
LIGHT-SPECULA		Nova Institution	Ax Room	Welch Allyn	739 series	20100720/20100714
LIGHT-ULTRAVIOLET	250-001445	Nova Institution	Dental suite	Demetron	A.2	
MONITOR-BIOFEEDBACK	250-001442	Nova Institution	Dental suite	Sinclair Dental	Root2XII	
MONITOR-NIBP	250-001555	Nova Institution	Emerg kit- stockroom	Welch Allyn	3400	100000050811
MONITOR-SPO2	25000 001652	Nova Institution	Treatment Room	Dynamap	2038172-001-219568	SH 615045495A
MONITOR-SPO2		Nova Institution	Minimum House Emergency kit	black	BMD-7594	
MONITOR-SPO2	250-001117	Nova Institution	Emerg kit		NPB40	G03860573
MONITOR-SPO2	250-001225	Nova Institution	Secure Unit	Rad-5V		538887
MONITOR-SPO2	250-001462	Nova Institution	Ax room	Dynamap	DPC400N	2019205-001
MONITOR-SPO2	250-001265	Nova Institution	Secure Unit	Dynamap	DCP400N	SDT101106085P
MONITOR-SPO2		Nova Institution	SECURE UNIT CART		VANTAGE RED	756/NN9590
MONITOR-TEMP	250-001553	Nova Institution	Tx Room	Braun	Thermoscan Plus	10105478202
MONITOR-TEMP	250-001554	Nova Institution	Ax Room	Braun	Thermoscan 6022	11110850664
MONITOR-TEMP	250-001563	Nova Institution	Emerg kit - stockroom	Braun	6022	11110890735
MONITOR-TEMP	250-001567	Nova Institution	Secure Unit emerg kit	Braun	Thermoscan	11110846300

Appendix M - Equipment List - Nova Institution

NEBULIZER-ULTRASONIC	250-001559	Nova Institution	Ax room	Sunmark	4075	045143
NEBULIZER-ULTRASONIC	250-001565	Nova Institution	Secure Unit	Sunmark	4075	040471
OPHTHALMIC-INSTRUMENT-STAND	25000-003997	Nova Institution	Dr office	Marco	DeLuxe II	1310499
OPHTHALMOSCOPE		Nova Institution	MINIMUM HOUSE		WELCH ALLYN	911/97150-M
PHOROPTER	25000-003994	Nova Institution	Dr office	Reichert	Minus	582150613
PUMP-INFUSION	250-001577	Nova Institution	Tx room	Plum A+		13897820
REFRIGERATOR-BIOLOGICAL	25000 002873	Nova Institution	Healthcare Medication room		PANASONIC	SANSR-L611W
REGULATOR-FLOWMETER-O2		Nova Institution	Emerg kit - stockroom	Medline		BMD-11634
REGULATOR-FLOWMETER-O2		Nova Institution	H/C Unit emerg kit	Medline		
SCALE-ELECTRONIC	250-001572	Nova Institution	Tx room	Health-o-Meter	2842KL	7E62045000069
SCALE-INFANT	250-001573	Nova Institution	Tx Room	Detecto		
SCALE-PATIENT	250-001556	Nova Institution	Inmate Bathroom	Detecto		
SPYGMOMANOMETER-ANALOG		Nova Institution	Tx Room	TYCOS	BMD-7554	90121134041
SPYGMOMANOMETER-ANALOG		Nova Institution	Ax Room	TYCOS	BMD-7557	9508090077
SPYGMOMANOMETER-ANALOG		Nova Institution	Emerg kit- stockroom	Welch Allyn	analog	
SPYGMOMANOMETER-ANALOG		Nova Institution	Emerg kit - secure unit	Welch Allyn	analog	
STERILIZER-STEAM-TABLE-TOP	250-001455	Nova Institution	Dental suite	Ritter	Midmark M9-022	V460343
STETHOSCOPE-ELECTRONIC		Nova Institution	Tx Room	Dr. Camscope	AMD-3700	21150
STIMULATOR-TENS	250-001575	Nova Institution	Tx room	Physiotens		5330588
STIMULATOR-TENS	250-001574	Nova Institution	Tx room	Physiotens		447523
STIMULATOR-TENS		Nova Institution	Tx Room	PhysioSD	TS-1211	11-049683
STIMULATOR-TENS		Nova Institution	Tx Room	PhysioSD	TS-1211	11-049685
STIMULATOR-TENS		Nova Institution		PhysioTens		5330585
STRETCHER	250-001457	Nova Institution	Tx room	Hill Rom	P8000	K043AQ4641
STRETCHER	BMD-11769	Nova Institution	Secure Unit Closet	Techlem	Premier Series 5000	
TABLE-EXAMINING	250-001568	Nova Institution	Secure Unit	Midmark	Ritter 104	
TABLE-EXAMINING	25000-001381	Nova Institution	Ax room	Ritter	Midmark 222	V1168926
TABLE-PHYSIO-TREATMENT	250-001460	Nova Institution	Tx room	Triton	TRT-600	2217
TONOMETER	25000-003995	Nova Institution	Dr office	Axis Medical	R-Type	05117
TRACTION-UNIT	250-001460	Nova Institution	Tx room	Chatanooga	2841	3317
ULTRASOUND-THERAPEUTIC	250-001460	Nova Institution	Tx room	Intelect	2762CC	4948
VENTILATOR-BIPAP-CPAP	2500-001722	Nova Institution	Storeroom	ResMed	HSI/S9	22112011401/22111819
VIEWER-XRAY-FILM		Nova Institution	Dental suite	Dux		
VISUAL-RESPONSE-ASSESSMENT	25000-003998	Nova Institution	Dr office	Huvitz	HDC9000N	9HC00N1310009
WHEELCHAIR	250-001570	Nova Institution	Tx Room	Medline		
WHEELCHAIR		Nova Institution	Stockroom	Airgo	Procare IC	
WHEELCHAIR	250-001459	Nova Institution	Healthcare Unit	AMG	700-631	100050-0512-0026
WHEELCHAIR	BMD-10517	Nova Institution	Physio Asst Room	AMG	700-628	
XRAY-UNIT	250-001450	Nova Institution	Dental suite	GenDex/Visualix		15148-3945ABG
XRAY-UNIT	250-001438	Nova Institution	Dental suite	GenDex	8500DDE	35-1858141DP
XRAY-UNIT PART	250-001453	Nova Institution	Dental suite	GenDex	Expert DC	12-1852698DP

Appendix N - Equipment List - Springhill Institution

Item	Equipment Code-Code d'équipement	Institution-Établissement	Location-Localisation	Manufacturer-Nom du fabricant	Model-Modèle	Serial Number-Numéro de série
AMALGAMATOR-DENTAL		Springhill Institution	dentist	Ultra Mat 2		71903
ANALYSER-GAS	004382	Springhill Institution	Fit Testing Office	Portacount Pro+	8038	80389433
ASPIRATOR-DENTAL	210-005009	Springhill Institution	Treatment	Gomco	4005	20071218026
ASPIRATOR-DENTAL	210-005004	Springhill Institution	Treatment	Devil Bless	7305	PD235889
BATH-PARAFFIN		Springhill Institution		Chattanooga Canada	Para Care	24753
CAMERA-CONTROLLER-LIGHT-SOURCE	BMD-11646	Springhill Institution	Back Office		DCSM(XGA)	DCSN-11B-002
CAMERA-VIDEO	BMD-11647	Springhill Institution	Back office	Sometech Inc	DCSN	DCSN-11C-117
CAMERA-VIDEO	MBD-11648	Springhill Institution	Back Office	SomeTech	DCS-104(NTSC)	DCSN-11B-058
CART-DENTAL	BMD-7855	Springhill Institution	Dentist	ADEC		
CHAIR-DENTAL		Springhill Institution	dentist	Belmont	X-Calibur	73203603
CHAIR-OPHTHALMIC	210-003191	Springhill Institution	Optometrist	Topcon	238717	222466
DENTAL - DRILL HANDPIECE		Springhill Institution	Dentist	Kavo	LUX E675L	
DENTAL - DRILL HANDPIECE		Springhill Institution	Dentist	Kavo	LUX E675L	
DENTAL-UNIT-ULTRASONIC	BMD-7850	Springhill Institution	Dentist	L&R Medical	Quantrex140	196231
DENTAL-UNIT-ULTRASONIC		Springhill Institution	dentist	DENSPLY	Cavitrion	11546484
ELECTROCARDIOGRAPH	210-004364	Springhill Institution	Treatment	GE	MAC5500	SCD10087696PA
GENERATOR-PARTICLE	004382	Springhill Institution	Fit Testing Office	Particle Generator	8026	80263697
HYDROCOLLATOR-HOT-PACK	210-005014	Springhill Institution		Chattanooga Canada	Hydrocollator -2102C	2328C
KERATOMETER	210-003190	Springhill Institution	Optometrist	Topcon	238714	01074G
LENSOMETER	BMD-7543	Springhill Institution		Bausch & Lomb	712655	512045
LENSOMETER		Springhill Institution	Health Services	Right Med	Model 1 S4OPTIK	
LIGHT-EXAMINING	BMD-7863	Springhill Institution	Treatment Room	Welch Allyn	44300/LS135	
LIGHT-EXAMINING-CEILING		Springhill Institution	dentist	Belmont 501T	105000 Dental Light	H71903
LIGHT-HEAD	210-003429	Springhill Institution	Optometrist	Heine	268299	0/11848
LIGHT-ULTRAVIOLET		Springhill Institution	dentist	Bobcat	506 QHL75	506-30012
LIGHT-ULTRAVIOLET		Springhill Institution	dentist	Kerr	Demi Plus	
MONITOR-NIBP-SPO2-TEMP	210-005006	Springhill Institution	Treatment	GE	2009284001	AAW080502158A
MONITOR-NIBP-SPO2-TEMP	BMD-7978	Springhill Institution		Welch Allyn	42MTB Spot Vital	201103450
MONITOR-NIBP-SPO2-TEMP	BMD-8084	Springhill Institution	RRC - Nursing Office	Welch Allyn	420 Series	201216120
MONITOR-NIBP-SPO2-TEMP	BMD-10430	Springhill Institution		Welch Allen	42NTB	
MONITOR-NIBP-SPO2-TEMP	BMD-8083	Springhill Institution	Dr. Office	Welch Allyn	420 Series	201216780
MONITOR-NIBP-SPO2-TEMP		Springhill Institution	Treatment Room/Port	Welch Allyn	42NTB	201403441
MONITOR-NIBP-SPO2-TEMP		Springhill Institution	Treatment Room	Welch Allyn	Vital Signs LXI	20150100638
MONITOR-NIBP-SPO2-TEMP		Springhill Institution	A&D	WELCH ALLYN		297 201403328
MONITOR-NIBP-SPO2-TEMP		Springhill Institution	Health Services	Welch Allyn	Connex VSM 6400	
MONITOR-SPO2	BMD-5337	Springhill Institution		Neilcor Puritan Bennett	NPB40	G00842991
MONITOR-SPO2	BMD-8080	Springhill Institution	Dr. Office	Neilcor	N-65P	P12827521
MONITOR-SPO2	210-003983	Springhill Institution			365367	G00842987
MONITOR-SPO2		Springhill Institution	TREATMENT ROOM	GE	TUFFSAT	FCB081100
MONITOR-SPO2		Springhill Institution	TREATMENT ROOM	GE	TUFFSAT	FCB081101
MONITOR-TEMP	BMD-5343	Springhill Institution		Welch Allyn	690 Suretemp	743382
MONITOR-TEMP	BMD-5345	Springhill Institution		Welch Allyn	679 Suretemp	3137764
MONITOR-TEMP		Springhill Institution	TREATMENT ROOM	BRAUN		6022 111111402
MONITOR-TEMP-TYMPANIC	BMD-10426	Springhill Institution		Welch Allen	6021	

Appendix N - Equipment List - Springhill Institution

MONITOR-TEMP-TYMPANIC	BMD-10425	Springhill Institution		Welch Allan	6022	
OPHTHALMIC-INSTRUMENT-STAND	210-003192	Springhill Institution	Optometrist	Topcon	238721	350236
OPHTHALMOSCOPE		Springhill Institution	Reception Nurse	Welch Allyn	GS777	
OPHTHALMOSCOPE		Springhill Institution	NP Office	Welch Allyn	GS777	
OPHTHALMOSCOPE	BMD-5346	Springhill Institution		Welch Allyn	767	
OPHTHALMOSCOPE	BMD-7871	Springhill Institution		Welch Allyn	767	
OPHTHALMOSCOPE	BMD-11645	Springhill Institution	Dr Office	Welch Allyn Canada	767	
PHOROPTOR	210-003189	Springhill Institution	Optometrist	Topcon	238712	3617864
PROJECTOR-OPHTHALMIC	210-003193	Springhill Institution	Optometrist	Topcon	238723	262933G
PUMP-COMPRESSOR	BMD-7862	Springhill Institution	Treatment Room	Devilbiss	5650C Pulmoaide	C170892
PUMP-COMPRESSOR	BMD-7865	Springhill Institution	Treatment	Invacare Corp	IRC1196 - Passport	59400000890
PUMP-COMPRESSOR	BMD-7866	Springhill Institution	Treatment	Invacare Corp	IRC1196 - PASSPORT	59400000908
PUMP-COMPRESSOR	BMD-7865	Springhill Institution	Treatment	Invacare Corp	IRC1196 - PASSPORT	59400000897
PUMP-COMPRESSOR	BMD-8081	Springhill Institution		Innovocorp	IRC1196	
PUMP-INFUSION	210-005008	Springhill Institution	Treatment	Hospira PlumA+	12910405	13898186
PUMP-VACCUUM		Springhill Institution	dentist	Henry Schien		07344-563850
REFRIGERATOR-BIOLOGICAL	BMD-10351	Springhill Institution	Medication Room	Sanyo	MPR721	11110815
REFRIGERATOR-BIOLOGICAL		Springhill Institution	Treatment Room	Panasonic		121001043
REFRIGERATOR-BIOLOGICAL	BMD-8079	Springhill Institution	Nursing Office - RRC	Panasonic	SR-L6111W	120700510
SAW-CAST		Springhill Institution	Treatment Room	BSN Medical Inc.	0295-200	28330
SCALE-PATIENT	BMD-10427	Springhill Institution		DET	6437	
SCALE-PATIENT		Springhill Institution	hallway	Health O Meter	2500KL	2500L00415
STERILIZER-STEAM-TABLE-TOP		Springhill Institution	dentist	Midmark	M7011 Speed Clave	MH003986
STRETCHER		Springhill Institution	NP Office	Stryker	CA LP 748	16
STRETCHER		Springhill Institution	Reception Nurse	Stryker	CA LP 748	16
STRETCHER	210-003427	Springhill Institution	Treatment	Hill Rom	268261	F093AW7612
STRETCHER	BMD-10416	Springhill Institution	Treatment Room	Pedigo	PP7500-W-Spec	0199025
TABLE-XRAY	210-00436	Springhill Institution	X-ray	GE	X171137 Picker	B25863
TONOMETER	210-001322	Springhill Institution	Optometrist	TONO-PENXL		T5932
TONOMETER	BMD-10433	Springhill Institution		CLC	MK2	
TOPCON COUNTERBALANCED MIDD	210-003194	Springhill Institution	Optometrist	Topcon	238725	2031000178
ULTRASOUND-THERAPEUTIC	210-005015	Springhill Institution		Chattanooga Canada	2762CC Intellect Advanced	4295
WARMER-BLANKET-SOLUTION	210-005008	Springhill Institution	Treatment	Pedigo	P-2010S	545087
WHEELCHAIR		Springhill Institution		AMG	700-625	
WHEELCHAIR		Springhill Institution		AMG	700-627	
WHEELCHAIR		Springhill Institution		AMG	700-631	
WHEELCHAIR		Springhill Institution		AMG	700-632	
XRAY-UNIT		Springhill Institution	X-Ray room	AGFA	CF30X	1380
XRAY-UNIT		Springhill Institution	X-ray room	AGFA	Drystar 5300	
XRAY-UNIT	210-004213	Springhill Institution	X-ray	GENDEX	Orthorifix 8500	35-1858143DP
XRAY-UNIT	210-004318	Springhill Institution	dentist	Belmont	428384	EH070201
XRAY-UNIT	210-000434	Springhill Institution	X-ray	GE	5892919	4034
XRAY-UNIT	210-000435	Springhill Institution	X-ray	GE	600-030	46335HL6