

[illegible]

NOTES:

- 1- QUANTITIES GIVEN IS FOR 1 (ONE) STRAINER , 2 (TWO) ARE REQUIRED.
- 2- ALL DIMENSIONS ARE IN MILLIMETERS EXCEPT AS NOTED.
- 3- STRAINER TO BE GALVANIZED AFTER FABRICATION.
- 4- BASKET RATIO : 6:4:1
- 5- WEIGHT : 780 KG
- 6- PRESSURE TEST 1000 KPA



3		AS FITTED FOR C-451	R.B.	07-05-
	S-D	AS FITTED UP TO 450	DP	06/01/84
	S-D	DIM. 172 WAS 766		
4	S-D	DIM. 732 WAS 734		
	S-C	DIM. 33 WAS 31	D.B.	06/17/84
3		CONNECTION FOR VENT MODIFIED	MN	06-11-84
		NOTE G ADDED	D.B.	06-18-84
		MATERIAL FOR ITEMS 3,12 AND 27 MODIFIED	MN	23/12/84
2				
1		COVER ATTACHMENT AND MATERIAL FOR THE BASKET MODIFIED	MN	5/12/84
NO	ZN	CODE	DESCRIPTION	PAR BY DATE

REVISIONS

APPROBATIONS/APPROVALS							
	SOU/SUB	APP	REV 1	REV 1	REV	REV	REV
CLASSN	83-10-02	83-12-81					
GCC/CGG	83-10-02	83-12-16					
PROP/OWN	83-10-02	83-12-02	83-12-12	83-12-15			



MARINE INDUSTRIE LIMITEE
SOREL P.Q. CANADA

TYPE 1100 NAVAIDS TENDER/LIGHT ICEBREAKER

SEA SUCTION STRAINER

PAR/BY <u>DENIS B.</u>	APPR <i>Johnson</i>	CONT NO <u>450</u> 451	DESSIN DWG	REV 5
VER/CK	ECN/SC <u>AS SHOWN</u> DATE 19-9-83	- 465-21		F/SH 1