



RETURN BIDS TO:

RETOURNER LES SOUMISSIONS À:

Bid Receiving - PWGSC / Réception des soumissions -
TPSGC

11 Laurier St./ 11, rue Laurier

Place du Portage, Phase III

Core 0B2 / Noyau 0B2

Gatineau, Québec K1A 0S5

FAX pour soumissions: (819) 997-9776

Revision to a Request for a Standing Offer

Révision à une demande d'offre à commandes

Regional Individual Standing Offer (RISO)

Offre à commandes individuelle régionale (OCIR)

The referenced document is hereby revised; unless
otherwise indicated, all other terms and conditions of
the Offer remain the same.

Ce document est par la présente révisé; sauf
indication contraire, les modalités de l'offre demeurent
les mêmes.

Comments - Commentaires

Vendor/Firm Name and Address

Raison sociale et adresse du
fournisseur/de l'entrepreneur

Issuing Office - Bureau de distribution

Clothing and Textiles Division / Division des
vêtements et des textiles
L'Esplanade Laurier,
East Tower 7th Floor
Tour est 7e étage
140 O'Connor, rue O'Connor,
Ottawa
Ontario
K1A 0R5

Title - Sujet Rib knit and spandex fabric Tricot côtelé et du tissu Lycra	
Solicitation No. - N° de l'invitation 21C31-202805/A	Date 2021-08-12
Client Reference No. - N° de référence du client 3332805	Amendment No. - N° modif. 002
File No. - N° de dossier pr759.21C31-202805	CCC No./N° CCC - FMS No./N° VME
GETS Reference No. - N° de référence de SEAG PW-\$\$PR-759-80255	
Date of Original Request for Standing Offer Date de la demande de l'offre à commandes originale 2021-07-29	
Solicitation Closes - L'invitation prend fin at - à 02:00 PM Eastern Daylight Saving Time EDT on - le 2021-10-05 Heure Avancée de l'Est HAE	
Address Enquiries to: - Adresser toutes questions à: Lafleur, Mario	Buyer Id - Id de l'acheteur pr759
Telephone No. - N° de téléphone (873) 354-0072 ()	FAX No. - N° de FAX () -
Delivery Required - Livraison exigée	
Destination - of Goods, Services, and Construction: Destination - des biens, services et construction:	
Security - Sécurité This revision does not change the security requirements of the Offer. Cette révision ne change pas les besoins en matière de sécurité de la présente offre.	

Instructions: See Herein

Instructions: Voir aux présentes

Acknowledgement copy required Accusé de réception requis	Yes - Oui <input type="checkbox"/>	No - Non <input type="checkbox"/>
The Offeror hereby acknowledges this revision to its Offer. Le proposant constate, par la présente, cette révision à son offre.		
Signature	Date	
Name and title of person authorized to sign on behalf of offeror. (type or print) Nom et titre de la personne autorisée à signer au nom du proposant. (taper ou écrire en caractères d'imprimerie)		
For the Minister - Pour le Ministre		

This amendment 002 is raised to answer questions from the industry.

QUESTION 1: About, Item 3 – Bias cut Spandex strips, 1-13/16 in. – White, in the specification it is noted that the mass must be 215.g/m² ± 3%; but the table shows a different required Unit mass of fabric. Could you please clarify which one is correct?

ANSWER 1: Yes it has to be 215.g/m² +3% . The annex B- Technical specifications is corrected as follows:

Item 3 – Bias cut Spandex strips, 1-13/16 in. - White

White 1/1 spandex strips, 48% polyester / 48% cotton / 4% spandex ± 5%;

The width of the strip must be 1-13/16 inch, bias-cut ± ¼ in.;

The mass must be 215.g/m² ± 3%;

The spandex must be packed in 10-inch diameter rolls.

The spandex must be the same white as the tubular rib knit and the bias rib knit.

- The fabric must meet all the requirements specified here under within the minimum and maximum acceptable range.
- The fabric must be tested in accordance with the test methods specified in the table below. **Only the test methods listed in the table below will be accepted.**

Property	Evaluation and acceptance procedure	Required specifications	Minimum acceptable	Maximum acceptable
Quantitative analysis of fibre mixtures	ISO / TR 11827 :2012 ASTM D276	48% polyester 48% cotton 4% spandex	- 5% (45.6%) - 5% (45.6%) -5% (3.8%)	+ 5% (52.50%) + 5% (52.50%) +5% (4.2%)
Unit mass of fabrics	CAN/CGSB 4.2 No. 5.1 – M90 (R2013) ISO 3801:1977 Method 5	215 g/m ²	-3% : 208g/m ²)	+3% : 222g/m ²)
Fabric width		1-13/16 in.	- ¼ in. (1-9/16 in)	+ ¼ in. (2-1/16 in.)

Solicitation No. - N° de l'invitation
21C31-202805/A
Client Ref. No. - N° de réf. du client
21C31--202805

Amd. No. - N° de la modif.
002
File No. - N° du dossier
pr759. 21C31--202805

Buyer ID - Id de l'acheteur
pr759
CCC No./N° CCC - FMS No./N° VME

QUESTION 2: Is that possible to send us some pictures of these items for our reference?

ANSWER 2: Yes, see below:

Item 1 - Tubular rib knit – White

White tubular 1/1 rib knit, 50% polyester / 50% brushed cotton $\pm 5\%$;
The mass must be of 175 g/m² - 3%;
The tube widths must be of 31 inches $\pm \frac{1}{2}$ inch;
The tubular rib knit must be the same white as the bias rib knit and the spandex.



Item 2 – Bias cut Rib knit 1 ½ in. strips – White

White 1/1 rib knit strips, 50% polyester / 50% brushed cotton $\pm 5\%$;
The width of the strip must be 1½-inch, bias-cut $\pm \frac{1}{4}$ in.;
The mass must be 175 g/m² $\pm 3\%$;
The bias rib knit must be packed in 10-inch diameter rolls.
The bias rib knit must be the same white as the tubular rib knit and the spandex.

Solicitation No. - N° de l'invitation
21C31-202805/A
Client Ref. No. - N° de réf. du client
21C31--202805

Amd. No. - N° de la modif.
002
File No. - N° du dossier
pr759. 21C31--202805

Buyer ID - Id de l'acheteur
pr759
CCC No./N° CCC - FMS No./N° VME



Item 3 – Bias cut Spandex strips, 1-13/16 in. - White

White 1/1 spandex strips, 48% polyester / 48% cotton / 4% spandex \pm 5%;
The width of the strip must be 1-13/16 inch, bias-cut \pm 1/4 in.;
The mass must be 215.g/m² \pm 3%;
The spandex must be packed in 10-inch diameter rolls.
The spandex must be the same white as the tubular rib knit and the bias rib knit.



ALL OTHER TERMS AND CONDITIONS REMAIN THE SAME.