



**RETURN BIDS TO:**

**RETOURNER LES SOUMISSIONS À:**

Bid Receiving Public Works and Government  
Services Canada/Réception des soumissions/Travaux  
publics et Services gouvernementaux Canada  
See herein for bid submission  
instructions/

Voir la présente pour les  
instructions sur la présentation  
d'une soumission

NA  
Ontario

**SOLICITATION AMENDMENT  
MODIFICATION DE L'INVITATION**

The referenced document is hereby revised; unless otherwise  
indicated, all other terms and conditions of the Solicitation  
remain the same.

Ce document est par la présente révisé; sauf indication contraire,  
les modalités de l'invitation demeurent les mêmes.

**Comments - Commentaires**

**Vendor/Firm Name and Address**  
**Raison sociale et adresse du**  
**fournisseur/de l'entrepreneur**

**Issuing Office - Bureau de distribution**  
Public Works and Government Services Canada  
Ontario Region  
10th Floor, 4900 Yonge Street  
Toronto  
Ontario  
M2N 6A6

<b>Title - Sujet</b> De l'équipement musical	
<b>Solicitation No. - N° de l'invitation</b> W3027-222298/A	<b>Amendment No. - N° modif.</b> 001
<b>Client Reference No. - N° de référence du client</b> W3027-222298	<b>Date</b> 2022-01-10
<b>GETS Reference No. - N° de référence de SEAG</b> PW-\$TOR-302-8153	
<b>File No. - N° de dossier</b> TOR-1-44103 (302)	<b>CCC No./N° CCC - FMS No./N° VME</b>
<b>Solicitation Closes - L'invitation prend fin</b> <b>at - à 02:00 PM</b> Eastern Standard Time EST <b>on - le 2022-01-19</b> Heure Normale du l'Est HNE	
<b>F.O.B. - F.A.B.</b> <b>Plant-Usine:</b> <input type="checkbox"/> <b>Destination:</b> <input checked="" type="checkbox"/> <b>Other-Autre:</b> <input type="checkbox"/>	
<b>Address Enquiries to: - Adresser toutes questions à:</b> Berends, Robert	<b>Buyer Id - Id de l'acheteur</b> tor302
<b>Telephone No. - N° de téléphone</b> (647) 308-5031 ( )	<b>FAX No. - N° de FAX</b> ( ) -
<b>Destination - of Goods, Services, and Construction:</b> <b>Destination - des biens, services et construction:</b>	

**Instructions: See Herein**

**Instructions: Voir aux présentes**

<b>Delivery Required - Livraison exigée</b>	<b>Delivery Offered - Livraison proposée</b>
<b>Vendor/Firm Name and Address</b> <b>Raison sociale et adresse du fournisseur/de l'entrepreneur</b>	
<b>Telephone No. - N° de téléphone</b> <b>Facsimile No. - N° de télécopieur</b>	
<b>Name and title of person authorized to sign on behalf of Vendor/Firm</b> <b>(type or print)</b> <b>Nom et titre de la personne autorisée à signer au nom du fournisseur/</b> <b>de l'entrepreneur (taper ou écrire en caractères d'imprimerie)</b>	
<b>Signature</b>	<b>Date</b>

La présente modification vise à reporter la date de clôture de l'invitation à soumissionner, à modifier les annexes A et B, et à répondre aux questions reçues jusqu'à présent.

Sur la page couverture, à la date de clôture de l'invitation à soumissionner :

Supprimer : 14 janvier 2022

Ajouter :

19 janvier 2022

À l'annexe A, Besoin :

Supprimer : en entier.

Ajouter :

Annexe A révisée, Besoin

## 1.0 SCOPE

### 1.1 Purpose

The Department of National Defence (DND) 31 Canadian Brigade Group has a requirement for the supply and delivery of music equipment in support of multiple professional military bands in the Southwestern Ontario area.

### 1.2 Background

Due to the professional nature of these bands and their importance to unit history and heritage all instruments must be of professional caliber and suitable for parade, concert and small ensemble work. The instruments supplied must be brand new and come with a manufacturer warranty.

## 2.0 GENERAL REQUIREMENTS

### 2.1 Scope of Requirement

The Contractor must be able to supply and deliver all music equipment prior to 31 March 2022.

### 2.2 Deliverables

#### 2.2.1 Brass/Reed

2.2.1.1 EUPHONIUM, Key of Bb- Besson 967 or Equivalent. Equivalent is defined as:

- Key of Bb
- 4 Valve (3 top, 1 side)
- .590" bore
- 12" bell
- Silver Plated
- Lyre post
- Quantity of 2

2.2.1.2 FRENCH HORN, Double Horn (Bb/F) - Hoyer G10L2 or Equivalent. Equivalent defined as:

- Double Horn (Bb/F)
- 4 Rotary Valves
- .468" bore
- 310mm Bell
- String Mechanism
- Quantity of 2

2.2.1.3 Bb CLARINET, Key of Bb- Buffet R13 Nickel keys or Equivalent. Equivalent is defined as:

- Key of Bb
- .577" bore
- 66mm barrel
- Premium quality granadilla wood
- 17 Key
- 6 ring nickel plated keys
- Quantity of 1

2.2.1.4 OBOE – Fox 335 or Equivalent. Equivalent is defined as:

- Grenadilla wood body
- Silver plated keys
- Quantity of 1

2.2.1.5 TROMBONE, Key of Bb- Bach 42BO or Equivalent. Equivalent is defined as:

- Key of Bb
- .547" bore
- 8.5" bell
- Lacquered Brass
- Open Wrap F attachment with standard rotor valve
- Quantity of 3

2.2.1.6 FLUTE, Key of C- Muramatsu or Equivalent. Equivalent is defined as:

- Key of C
- B Foot
- Offset G
- Solid Silver Head Joint
- Solid Silver Body
- Quantity of 3

2.2.1.7 ALTO SAX, Key of Eb- Saxophone Alto Yamaha YAS62III or Equivalent. Equivalent is defined as:

- Key of Eb
- High F# Key
- Rocking Table Mechanism
- Fully ribbed construction
- Pro Leather pads
- Yellow Brass
- Quantity of 1

2.2.1.8 TRUMPET, Key of Bb- Yamaha YTR8335SII Xeno Series or Equivalent. Equivalent is defined as:

- Key of Bb
- 11.65mm bore
- 123mm bell
- Heavy Weight
- Silver Plated
- Quantity of 2

2.2.1.9 TENOR SAX, Key of Bb- Yamaha YTS62III or Equivalent. Equivalent is defined as:

- Key of Bb
- Front F and High F# Key
- Leather Pads
- Yellow Brass
- Quantity of 1

2.2.1.10 PICCOLO, Key of C- Yamaha YPC62R or Equivalent. Equivalent is defined as:

- Key of C
- Split E mechanism
- Undercut tone holes
- Conical bore
- Grenadilla Wood
- Quantity of 3

2.2.1.11 TUBA, Key of Bb- Yamaha YBB632S or Equivalent. Equivalent is defined as:

- Key of Bb
- 4 Valves (3 Upright, 1 side)
- 500mm bell
- Silver Plated
- Quantity of 1

2.2.1.12 TUBA, Key of Bb- Yamaha YBB641E or Equivalent. Equivalent is defined as:

- Key of Bb
- 4 Rotary Valves
- 4/4 Size
- 16-1/2" bell
- .811" bore
- Brass Lacquered
- Comes with Mouthpiece and case
- Quantity of 1

2.2.1.13 MARIMBA with Soft Bag Case, 3.3 Octaves- 3.3 Academy AMPD3RR or Equivalent. Equivalent is defined as:

- 3.3 Octaves
- Alpha Tuning System
- Aluminum Beam
- Integrated Resonator
- Rail System
- 3 Position Adjustable X Stand
- Soft Bag Case
- Quantity of 1

2.2.1.14 FRENCH HORN with Case, Key of F/Bb, Yamaha YHR-671D or Equivalent. Equivalent is defined as:

- Key of F/Bb
- M Bell Size
- Detachable Bell
- 12mm (0.472") bore
- 4 Rotary Valves
- String Lever Action
- Yellow Brass Clear Lacquer
- Nickel Silver Valve Casings
- Case
- Quantity of 1

2.2.1.15 EUPHONIUM with Hard Case, Key of Bb/F, Yamaha YEP642 or Equivalent. Equivalent is defined as:

- Key of Bb/F
- 15-16.8mm (0.591-0.661) Dual Bore
- 300mm (11-4/5") Bell
- 3 monel valve top action/ 1 monel valve side action
- Auto-compensating system
- Large Shank Mouthpiece
- Yellow Brass with Clear Epoxy Lacquer
- Hard Case
- Quantity of 1

## 2.2.2 PIPE AND DRUMS

- 2.2.2.1 15" Tenor Drum, ANDANTE Pro Series Model 361.
- Professional Class Pipe Band Tenor Drum
  - 15" heads
  - 12" shell depth
  - Red Woodgrain Finish
  - Silver Metalwork
  - With Hard Carrying Case
  - Quantity of 4
- 2.2.2.2 16" Tenor Drums, ANDANTE Pro Series Model 362 .
- Professional Class Pipe Band Tenor Drum
  - 16" heads
  - 12" shell depth
  - Red Woodgrain Finish
  - Silver Metalwork
  - With Hard Carrying Case
  - Quantity of 4
- 2.2.2.3 18" Tenor Drums, ANDANTE Pro Series Model 363 .
- Professional Class Pipe Band Tenor Drum
  - 18" Heads
  - 12" Shell Depth
  - Red Woodgrain Finish
  - Silver Metalwork
  - With Hard Carrying Case
  - Quantity of 2
- 2.2.2.4 Side Drums, ANDANTE Next Generation Reactor.
- Professional Class Pipe Band Side Drum
  - High-Tension
  - Top and Bottom Snares
  - 14" Heads
  - 12" Shell Depth
  - Red Woodgrain Finish
  - Silver Metalwork
  - With Hard Carrying Case
  - Quantity of 8
- 2.2.2.5 28" Bass Drum, ANDANTE Pro Series Model 260.
- Professional Class Pipe Band Bass Drum
  - 28" Heads
  - 14" Shell Depth
  - Red Woodgrain Finish
  - Silver Metalwork
  - With Hard Carrying Case
  - Quantity of 1

- 2.2.2.6 14" Drum Heads, ANDANTE CORE-TEC Elite White Drum Head.
- 14" White Drum Head
  - Quantity of 6
- 2.2.2.7 Bb Infantry Pattern Bugles, Viscount Bb Bugle or Equivalent. Equivalent is defined as:
- Tunable with 5C Silver Plate Mouthpiece
  - Polished Brass
  - Carry Case
  - Quantity of 28
- 2.2.2.8 Great Highland Bagpipes, Naill PN1 Model Fully Assembled or Equivalent. Equivalent is defined as:
- Polypenco fully combed/beaded
  - Nickel ferrules and ring caps
  - Imitation ivory projecting mounts
  - Quantity of 6
- 2.2.2.9 Bagpipes, Piper's Choice Full Imitation Bagpipes Fully Assembled or Equivalent. Equivalent is defined as:
- Plastic Bagpipes
  - Bag - Bannatyne Zip Medium
  - Bag Cover - Green – Velveteen
  - Trim - Silver – Metallic fringe
  - Cords - Silver – Silk
  - Chanter Reed – Medium
  - Chanter - Poly Chanter
  - Drone Reeds - Ezeedrone
  - Quantity of 1
- 2.2.2.10 Great Highland Bagpipes, Naill DN5 Model Fully Assembled or Equivalent. Equivalent is defined as:
- African Blackwood
  - Hand engraved ferrules, slides, ring caps and mouthpiece tube
  - Imitation ivory projecting mounts
  - 0.925 Sterling
  - Quantity of 3
- 2.2.2.11 Great Highland Bagpipes, McCallum AB4 Model Fully Assembled or Equivalent. Equivalent is defined as:
- African Blackwood
  - Engraved Nickel Ferrules, Slides and ring caps
  - Imitation Ivory projecting mounts
  - Quantity of 6

- 2.2.2.12 Great Highland Bagpipe Bb Orchestral Chanter, Shepherd or Equivalent. Equivalent is defined as:
- Bb orchestral Bagpipe Chanter
  - Quantity of 39
- 2.2.2.13 Great Highland Bagpipe Standard Chanter, McCallum or Equivalent. Equivalent is defined as:
- Standard Chanter
  - Polypenco without engraving MKIII
  - Quantity of 45
- 2.2.2.14 Electronic Bagpipe, Blair or Equivalent. Equivalent is defined as:
- Plain sole on chanter without engraving
  - With carrying case
  - Quantity of 18
- 2.2.2.15 Great Highland Bagpipe Bag (Extended Small), Bannatyne Hybrid Bag or Equivalent. Equivalent is defined as:
- 9.5" Deep by 26.5" Long
  - Hybrid Hide Synthetic Bag
  - With zipper for easy access to inside of bag
  - Bottle Water trap to be included
  - Quantity of 2
- 2.2.2.16 Great Highland Bagpipe Bag (Medium), Bannatyne Hybrid Bag or Equivalent. Equivalent is defined as:
- 10.5" Deep by 26.5" Long
  - Hybrid Hide Synthetic Bag
  - With zipper for easy access to inside of bag
  - Bottle Water trap to be included
  - Quantity of 6
- 2.2.2.17 Great Highland Bagpipe Bag Cover (Medium), Bag Cover Deluxe or Equivalent. Equivalent is defined as:
- MacGregor Tartan as associated with the Essex and Kent Regiment
  - Size Medium
  - With zipper for easy access
  - Velcro close at back
  - Quantity of 12
- 2.2.2.18 Great Highland Bagpipe Bag Cover (Extended Small), Bag Cover Deluxe or Equivalent. Equivalent is defined as:
- MacGregor Tartan as associated with the Essex and Kent Regiment
  - Size Small
  - With zipper for easy access
  - Velcro close at back
  - Quantity of 3



### 3.0 TECHNICAL AUTHORITIES

31 Canadian Brigade Group G1 Bands

Name: Warrant Officer Ben Jaeckle

Phone: 905-870-3695

Email: [argylldm@hotmail.com](mailto:argylldm@hotmail.com) or [Bernhard.Jaeckle@forces.gc.ca](mailto:Bernhard.Jaeckle@forces.gc.ca)

31 Canadian Brigade Group G4 MM 2

Name: Master Warrant Officer Michael York

Phone: 519-660-5275 ext 5315

Email: [Michael.York@forces.gc.ca](mailto:Michael.York@forces.gc.ca)

À l'annexe B, Base de paiement :  
 Supprimer : en entier.  
 Ajouter :

Annexe B révisée, Base de paiement

Prices are firm, all-inclusive in Canadian dollars, FOB Destination. Transportation charges, Travel and material cost, Customs duties and Excise taxes are included and Applicable Taxes are extra.

INSTRUMENT TYPE	INSTRUMENT DESCRIPTION	TOTAL QUANTITY (COLUMN A)	FIRM PRICE PER UNIT (COLUMN B)	EXTENDED PRICE COLUMN C=(AxB)
BRASS AND REED INSTRUMENTS 2.2.1	Euphonium Key of Bb- Besson 967 or Equivalent	2		
	French Horn Double Horn Bb/F- Hoyer G10L2 or Equivalent	2		
	Clarinet Bb- Buffet R13 or Equivalent	1		
	Oboe Grenadilla Wood- Fox 335 or Equivalent	1		
	Trombone Bb- Bach 42BO or Equivalent	3		
	Flute Key of C- Muramatsu or Equivalent	3		
	Alto Sax Eb- Saxophone Alto Yamaha YAS62III or Equivalent	1		
	Trumpet Bb- Yamaha YTR8335SII or Equivalent	2		
	Tenor Sax Bb- Yamaha YTS62III or Equivalent	1		
	Piccolo key of C- Yamaha YPC62R or Equivalent	3		
	Tuba Bb Type 1- Yamaha YBB632S or Equivalent	1		
	Tuba Bb Type 2- Yamaha YBB641E or Equivalent	1		
	Marimba 3.3 Octaves w/ case- 3.3 Academy AMPD3RR or Equivalent	1		
	French Horn with Case Key of F/Bb- Yamaha YHR671D or Equivalent	1		
	Euphonium key of Bb/F w/ Hard Case- Yamaha YEP642 or Equivalent	1		

INSTRUMENT TYPE	INSTRUMENT DESCRIPTION	TOTAL QUANTITY (COLUMN A)	FIRM PRICE PER UNIT (COLUMN B)	EXTENDED PRICE COLUMN C=(AxB)
PIPE AND DRUM INSTRUMENTS 2.2.2	Tenor Drum (15")- ANDANTE Pro Series Model 361	4		
	Tenor Drum (16")- ANDANTE Pro Series Model 362	4		
	Tenor Drum (18")- ANDANTE Pro Series Model 363	2		
	Side Drums- ANDANTE Next Generation Reactor or Equivalent	8		
	Bass Drum (28")- ANDANTE Pro Series Model 260	1		
	Drum Heads (14")- ANDANTE Core-Tec Elite White	6		
	Infantry Pattern Bugles Bb- Viscount Bb Bugle or Equivalent	28		
	Great Highland Bagpipes Polypenco- Naill PN1 Model Fully Assembled or Equivalent	6		
	Plastic Bagpipes- Piper's choice full imitation Bagpipes Fully Assembled, Ready to Play or Equivalent	1		
	Great Highland Bagpipes African Blackwood- Naill DN5 Model Fully Assembled or Equivalent	3		
	Great Highland Bagpipes African Blackwood- McCallum AB4 Model Fully Assembled or Equivalent	6		
	Bagpipe Bb Orchestral Chanter- Shepherd brand or Equivalent	39		
	Bagpipe Standard Chanter- McCallum or Equivalent	45		
	Electronic Bagpipe- Blair or Equivalent	18		
	Bagpipe Hybrid Bag Extended Small- Bannatyne Hybrid Bag or Equivalent	2		
	Bagpipe Hybrid Bag Medium- Bannatyne Hybrid Bag or Equivalent	6		
	Bagpipe Bag Cover Medium MacGregor Tartan- as associated with the Essex and Kent Regiment - Bannatyne Medium or Equivalent	12		

	Bagpipe Bag Cover Small MacGregor Tartan- as associated with the Essex and Kent Regiment -Bannatyne Extended Small or Equivalent	3		
--	--	---	--	--

Questions et réponses

Q1. – Y aura-t-il une préférence pour la marque Andante?

R1. – La sélection de chaque article sera fondée sur le prix le plus bas pour l'article demandé ou un substitut acceptable lorsqu'il est indiqué « ou l'équivalent ».

Q2. – Concernant le point 2.2.2.17, de quel tartan MacGregor s'agit-il?

R2. – Le tartan MacGregor associé au régiment Essex et Kent (Essex and Kent Regiment).

Q3. – Concernant le point 2.2.2.18, veulent-ils une taille « petite » (small) pour les housses de sac (covers) ou veulent-ils dire « petite allongée » (extended small)?

R3. – Voir les annexes révisées.

Q4. – Concernant les points 2.2.2.8, 2.2.2.9, 2.2.2.10, et 2.2.2.11, on peut lire « ensemble complet avec poche d'air (sac), housse, cordons et anches et prêt à être utilisé (complete with pipe bag, cover, cords and reeds and ready to play) ». Est-ce exact?

R4. – Oui, un ensemble complet prêt à être utilisé est requis.

Q5. – Les soumissions partielles sont-elles acceptables?

R5. – Oui.

Q6. – Le coût de l'expédition à partir d'ici (île de Vancouver) jusqu'à London (Ontario) doit-il être séparé ou inclus?

R7. – Voir l'annexe B.

Q7. – Concernant le point 2.2.2.11, la description concerne une cornemuse McCallum AB4 (AB4 McCallum Pipe). Le modèle AB3 Deluxe n'est livré qu'avec des glissières et des plaques gravées (engraved slides and plates), mais pas avec des viroles et des bagues (ferrules and ring caps) comme demandé. Je crois que vous voulez une AB4.

R7. – Voir les annexes révisées.

Q8. – Les points 2.2.2.8, 2.2.2.10 et 2.2.2.11 ne précisent pas quel groupe : 1) tuyaux (sticks) seulement, 2) tuyaux et anche (sticks and chanter) 3) entièrement assemblée?

R8. – Entièrement assemblée.

Q9. – Concernant le point 2.2.2.10, Naill DN5 Pipes : la description doit-elle correspondre au modèle DN5?

R9. – Voir les annexes révisées.