

**Nature of Requirements:** Millipore Milli-Q Reference Ultrapure Water Purification System, 2022-12707, Jessica Chartrand, Jessica.Chartrand2@forces.gc.ca

THIS REQUEST IS RESERVED FOR HOLDERS OF SUPPLY ARRANGEMENT # E60PV-19EQUI ONLY, FOR THE PURCHASE OF LABORATORY AND SCIENTIFIC EQUIPMENT, PARTS AND ACCESSORIES, SERVICES AND SUPPLIES.

Only Suppliers currently pre-qualified on Supply Arrangement E60PV-19EQUI have been invited to bid.

As a requirement of the Supply Arrangement, this notice is published on Buy and Sell Canada for a period of 40 calendar days. The closing date published on this notice identifies how long the notice will be published. For the closing date of any solicitation under the supply arrangement, invited suppliers should refer to the solicitation documents.

Suppliers that do not have a Supply Arrangement for the supply of Laboratory and Scientific Equipment, Parts and Accessories, Services and Supplies with Public Works and Government Services Canada, cannot submit a bid. Any bids received from suppliers not pre-qualified on the Supply Arrangement will not be evaluated.

Suppliers may qualify under Supply Arrangement E60PV-19EQUI for Laboratory and Scientific Equipment, Parts and Accessories, Services and Supplies, at any time. Interested suppliers should download solicitation document E60PV-19EQUI /C from Buyandsell.gc.ca and submit a response as per the requirements of the Request for Supply Arrangement.

The Department of National Defence has a requirement for:

Items:

1. Quantity: 2 - Z00QSV0US Millipore Milli-Q Reference Ultrapure Water Purification System
2. Quantity: 10 - QTUM0TEX1 Millipore Quantum TEX Cartridge PK1
3. Quantity: 10 - CDUFBI001 Millipore biopak
4. Quantity: 2 - WMBSMT002 Millipore wall mounting bracket -system
5. Quantity: 2 - ZSPRVBG01 millipore NPTF Pressure regulator with guage
6. Quantity: 2 - ZMQSP0D02 Millipore Q-POD 2 remote dispenser pk/1 (included with water system)
7. Quantity: 10 - QGARDT1x1 Millipore Q-GARD T1 pack pk/1
8. Quantity: 3 - ZMQUVLP01 UV lamp for Millipore milli-q reference
9. Quantity: 8 - 04141 Thermo (Barnstead) UV lamp
10. Quantity: 8 - 25001-DB Barnstead Tank vent filter
11. Quantity: 2 - Ventgard air filter, Barnstead H3120 Protects stored water against airborne particulates, organics and CO2
12. Quantity: 1 - Barnstead water seal H3130 Replacement water seal

This is a no substitution solicitation due to infrastructure requirements. The responsive bidder already holding a Supply Arrangement for Laboratory and Scientific Equipment, Parts and Accessories, Services and Supplies from Public Works and Government Services Canada who meets the mandatory technical criteria with the lowest aggregate price will be recommended for contract award.