



National Defence

Défense nationale

National Defence
Headquarters
Ottawa, Ontario
K1A 0K2

Quartier général de la Défense nationale
Ottawa (Ontario)
K1A 0K2

**REQUEST FOR STANDING OFFER
DEMANDE D'OFFRE À
COMMANDES**

**AMENDMENT NO.
N° DE LA MODIFICATION
001**

**RETURN BIDS TO /
RETOURNER LES SOUMISSIONS À :**

Attn: Adam Cheung, DLP 8-1-1
Adam.Cheung@forces.gc.ca

Proposal To: National Defence Canada

We hereby offer to sell to His Majesty the King in right of Canada, in accordance with the terms and conditions set out herein, referred to herein or attached hereto, the goods and services listed herein and on any attached sheets at the price(s) set out therefore.

Proposition à : Défense nationale Canada

Nous offrons par la présente de vendre à Sa Majesté la Reine du chef du Canada, aux conditions énoncées ou incluses par référence dans la présente et aux annexes ci-jointes, les biens et services énumérés ici et sur toute feuille ci-annexée, au(x) prix indique(s).

Title/Titre Forensic Sampling Equipment		Solicitation No – N° de l'invitation W6399-24LP79/A	
Date of Solicitation – Date de l'invitation 31 May 2024			
Address Enquiries to – Adresser toutes questions à Adam.Cheung@forces.gc.ca			
Telephone No. – N° de téléphone 613-945-2730		FAX No – N° de fax	
Destination Specified herein			

Instructions:

Municipal taxes are not applicable. Unless otherwise specified herein all prices quoted must include all applicable Canadian customs duties, GST/HST, excise taxes and are to be delivered Delivery Duty Paid including all delivery charges to destination(s) as indicated. The amount of the Goods and Services Tax/Harmonized Sales Tax is to be shown as a separate item.

Instructions: Les taxes municipales ne s'appliquent pas. Sauf indication contraire, les prix indiqués doivent comprendre les droits de douane canadiens, la TPS/TVH et la taxe d'accise. Les biens doivent être livrés « rendu droits acquittés », tous frais de livraison compris, à la ou aux destinations indiquées. Le montant de la taxe sur les produits et services/taxe de vente harmonisée doit être indiqué séparément.

Solicitation Closes – L'invitation prend fin
At – à : 1400hr EDT
On - le : 17 July 2024

Delivery required - Livraison exigée See herein/ Précisé dans les présentes	Delivery offered - Livraison proposée
Vendor Name and Address - Raison sociale et adresse du fournisseur	
Name and title of person authorized to sign on behalf of vendor (type or print) - Nom et titre de la personne autorisée à signer au nom du fournisseur (caractère d'imprimerie)	
Name/Nom _____	Title/Titre _____
Signature _____	Date _____

QUESTIONS AND ANSWERS

Question 1

Line 1 UV Lights (Battlefield). The title lists UV light but the Description mentions Infrared light. Can you confirm it is the UV light that is required? UV and IR are at opposite ends of the light spectrum and are quite different.

Answer 1

DND requires the UV light.

Question 2

D100 Light. This spec outlines the white light only version, but we have in the past supplied your troops with the 265 nm versions used to fluoresce Sciluminate. I want to ensure the spec here is also correct. The light specified will not fluoresce the Sciluminate powder that is used in your operations. If the use is for operation with Sciluminate the correct light is a D100 with 265 nm and 450 nm wavelengths in one D100 light.

Answer 2

DND requires the dual wavelength version (265&450 nm).

Question 3

Our company is a supplier of high end digital photographic imaging technology to institutional and Gov. applications. In reading through the solicitation I understand we are not best suited to help in the entire request... however I believe we might have a rather unique technology in capturing images for the UV, IR and white light spectrum. (this tech. is currently in use for CDN FIS law enforcement and various cultural heritage, Museums, globally)

So my question is this... can we participate in just this portion of the tender request with possible hardware solution?

Answer 3

As stated in the RFSO, DND is seeking a single Offeror that provides the best overall offer for the entirety of the requirement. Offers that do not meet all mandatory technical evaluation criteria will be declared non-responsive and not considered any further.

Question 4

In reference to 3.0 MANDATORY REQUIREMENTS, section 1.6(a) v "Quantity one (1) package of twenty-five (25) item C12: Fingerprint Lift Tape Pad" and section 1.6(b) iii "Quantity one (1) package of twenty-five (25) item C12: Fingerprint Lift Tape Pad",

Q1: Can you please confirm that this refers to one (1) package of twenty-five (25) pre-cut strips of clear, air bubble free fingerprint tape attached to a still poly pad?

Answer 4

Yes, quantity of one package containing one pad of 25 strips for each.

Note: To address the questions 1 & 2 above, Annex C & Annex D are being revised as follows:

Annex C - TECHNICAL REQUIREMENTS FOR FORENSIC SAMPLING EQUIPMENT

Ref: 1.4 Forensic Light

Delete: Item (a) Item L1: Ultra-Violet (UV) Forensic Light in its entirety and **insert** the following.

- (a) Item L1: Ultra-Violet (UV) Forensic Light – The ultra-violet forensic light must be the Lynn Peavey Company BattleLite UV model #56038 (or equivalent) as follows:
- i. 3-Watt forensic light;
 - ii. LED UV light head (395 nm wavelength);
 - iii. Hi-Lo power switch;
 - iv. Single thumb-operated on/off switch;
 - v. Shock, dust and moisture resistant; and
 - vi. Quantity: Package of one (1);

Delete: Item (c) Item L3: Forensic Photography Light in its entirety and **insert** the following.

- (c) Item L3: Forensic Photography Light – The forensic photography light must be the FoxFury D100 Photography Scene LED Light System model #700-DH-341 (or equivalent) as follows:
- i. Output: 2300 lumen, 15-degree beam, Dual Wavelength (265nm and 450 nm);
 - ii. Rechargeable lithium ion batteries (battery charger included);
 - iii. Life: Up to 70 hours on a single charge;
 - iv. Hand-held with included handle;
 - v. Permit mounting to a standard camera tripod; and
 - vi. Quantity: Package of one (1).

Annex D - EVALUATION PLAN FOR FORENSIC SAMPLING EQUIPMENT

Delete: Line Item #1 (Item L1: Ultra-Violet (UV) Forensic Light) and Line Item #3 (Item L3: Forensic Photography Light) in its entirety and **insert** the following:

Item #	Annex B Ref.	Requirement	Instructions to Offerors	Offer Reference
1	1.4(a)	<p><u>Item L1: Ultra-Violet (UV) Forensic Light</u> – The ultra-violet forensic light must be the Lynn Peavey Company BattleLite UV model #56038 (or equivalent) as follows:</p> <ul style="list-style-type: none"> i. 3-Watt forensic light; ii. LED UV light head (395 nm wavelength); iii. Hi-Lo power switch; iv. Single thumb-operated on/off switch; v. Shock, dust and moisture resistant; and vi. Quantity: Package of one (1); 	<p>The Offeror must provide:</p> <p>(1) Written confirmation they will be providing the Lynn Peavey Company BattleLite UV model #56038;</p> <p>OR</p> <p>(2) Equivalent Only - If equivalent components are being offered, the Offeror must provide a Certificate of Compliance that confirms the Ultra-Violet (UV) Forensic Light being offered is equivalent or better in every way to the Lynn Peavey Company BattleLite UV model #56038;</p>	
3	1.4(c)	<p><u>Item L3: Forensic Photography Light</u> – The forensic photography light must be the FoxFury D100 Photography Scene LED Light System model #700-DH-341 (or equivalent) as follows:</p>	<p>The Offeror must provide:</p> <p>(1) Written confirmation they will be providing the FoxFury</p>	

		<ul style="list-style-type: none">i. Output: 2300 lumen, 15-degree beam, Dual Wavelength (265 nm and 450 nm);ii. Rechargeable lithium ion batteries (battery charger included);iii. Life: Up to 70 hours on a single charge;iv. Hand-held with included handle;v. Permit mounting to a standard camera tripod; andvi. Quantity: Package of one (1).	D100 Photography Scene LED Light System model #700-DH-341; OR (2) Equivalent Only - If equivalent components are being offered, the Offeror must provide a Certificate of Compliance that confirms the Forensic Photography Light being offered is equivalent or better in every way to the FoxFury D100 Photography Scene LED Light System model #700-DH-341;	
--	--	--	--	--

Note: Please refer to RFP AMD001 for the changes noted above.

All other 'Terms and Conditions' remain the same.