

Annex F - Equipment List - Dorchester Institution (Medium)/ Annexe F - Liste d'équipements - Établissement Dorchester (Médium)

| | A | B | C | D | E | F | G |
|----|---|---------------------------------------|--------------------------------|----------------------------|------------------------------------|---------------------------------|------------------------------------|
| 1 | Item | ~~Equipment Code~~Code d'équipement~~ | ~~Institution~~Établissement~~ | ~~Location~~Localisation~~ | ~~Manufacturer~~Nom du fabricant~~ | ~~Model~~Modèle~~ | ~~Serial Number~~Numéro de série~~ |
| 2 | Amalgamator - Dental | BMD-8055 | Dorchester Medium Institution | Dental Office | The L.D Caulk Comp. | VM-C | 54154 |
| 3 | Aspirator | BMD-10561 | Dorchester Medium Institution | Clean Utility Room | Allied | L190-G4 | 0714 |
| 4 | Aspirator | BMD-16033 | Dorchester Medium Institution | Health Services | Allied | L190-GR | 0813 |
| 5 | Aspirator | BMD-8064 | Dorchester Medium Institution | Trauma Room | Gomco by Allied | 405 | 20080331005 |
| 6 | Aspirator | BMD-5218 | Dorchester Medium Institution | Treatment room 1 | Contemporary Products | 6260 | AS21745 |
| 7 | Assembly cam-14 | N/A | Dorchester Medium Institution | Treatment room 1 | GE Medical Systems | CAM-14 | RNC11015386PA |
| 8 | Auto Refractometer Cradle | N/A | Dorchester Medium Institution | Optometry Room | NIDEK | AR20ST1 | 130493 |
| 9 | Auto refractometer Hand Held | N/A | Dorchester Medium Institution | Optometry Room | NIDEK | AR20 | 230454 |
| 10 | Back board - Scoop Yellow | N/A | Dorchester Medium Institution | Back Hall | Ferno | 06-010-8002-65EXL | 088752 |
| 11 | Bath - Water | BMD-8071 | Dorchester Medium Institution | Dental Office | Whip Mix | 5350 | 07110195 |
| 12 | Camera - Controller - Light - Source - Camscope box | BMD-10377 | Dorchester Medium Institution | IDN office | Sometec | DCSM(X6A) | DCSM-10G-167 |
| 13 | Camera - Controller - Light - source (Camscope angled wand) | BMD-10378 | Dorchester Medium Institution | IDN office | Sometec | DCS-104T (NTSC) | DCSN-10E-217 |
| 14 | Camera - controller - light - source (Camscope straight wand) | BMD-10379 | Dorchester Medium Institution | IDN office | Sometec | P101 | DCSN-11C-116 |
| 15 | Cart - Dental (Chair Instruments) | | Dorchester Medium Institution | Dental Office | Belmont | | 12057301 |
| 16 | Cart - Emergency | BMD-11665 | Dorchester Medium Institution | Storage room 3rd floor | Harloff Co. | N/A | 46269 |
| 17 | Cart - Medication | BMD-11964 | Dorchester Medium Institution | ? | Health Care Logistics | | |
| 18 | Cart - Pandemic | BMD-12003 | Dorchester Medium Institution | Trauma Room | Harloff | Harloff Classic Line Red - 6510 | N/A |
| 19 | Cart - Video - Monitor | BMD-11663 | Dorchester Medium Institution | IDN office | Polycom | PPCE-120BP | N091215032 |
| 20 | Centrifuge | BMD-10562 | Dorchester Medium Institution | Trauma Room | Eppendorf | 5702 | 5702EN339263 |
| 21 | Chair - Blood work - Hydraulic | BMD-11965 | Dorchester Medium Institution | Trauma Room | Med Care | | 125924 |
| 22 | Chair - Dental | BMD-7418 | Dorchester Medium Institution | Dental Office | Takara Belmont | 447193 | 012057301 |
| 23 | Chair - Ophthalmic | BMD-1483 | Dorchester Medium Institution | | | FX920 | 93010707007 |

Annex F - Equipment List - Dorchester Institution (Medium)/ Annexe F - Liste d'équipements - Établissement Dorchester (Médium)

| | A | B | C | D | E | F | G |
|----|--|---------------------------------------|--------------------------------|-----------------------------------|------------------------------------|---------------------------|------------------------------------|
| 1 | Item | ~~Equipment Code~~Code d'équipement~~ | ~~Institution~~Établissement~~ | ~~Location~~Localisation~~ | ~~Manufacturer~~Nom du fabricant~~ | ~~Model~~Modèle~~ | ~~Serial Number~~Numéro de série~~ |
| 24 | Chair - Ophthaomalic | BMD-5225 | Dorchester Medium Institution | X-ray room | Power Tilt | FX920L | 93010707007 |
| 25 | Chair - Recline - Tilt (Blood Drawing Chair) | BMD-10552 | Dorchester Medium Institution | Trauma Room | Winco | 574 | 574100108 |
| 26 | Chair - Recline - Tilt (Blood Drawing Chair) | BMD-10553 | Dorchester Medium Institution | IDN Office | Winco | 570 | 570B102846 |
| 27 | Chair Recline Tilt | BMD-10563 | Dorchester Medium Institution | Clean Equipment Storage Room | Lumex | FR574G 454-PRO-Lab RevA06 | 717076030168 / HSP321SS0234 |
| 28 | Cleaner - Ultrasonic | BMD-11652 | Dorchester Medium Institution | Dental Office | L & R Manufacturing | T-14B | 150 |
| 29 | Cleaner - Ultrasonic (Dental) | BMD-8056 | Dorchester Medium Institution | Dental Office | Densply | cAVITRON SPS | 119-27469 |
| 30 | Collimator - XRAY (Eclipse Proteus Collimator) | BMD-11668 | Dorchester Medium Institution | X-ray room | GE Hualun Medical System | 2379827 | 53339HL8 |
| 31 | Compressed Air | | Dorchester Medium Institution | Dental Office | Belmont | 4000 | 20213 |
| 32 | Compressor | | Dorchester Medium Institution | Dental Office | Dominic Hunter | DHDRD10 | 3701300004 |
| 33 | Controls for X-ray machine | | Dorchester Medium Institution | X-ray room | GE Medical System | 2320993 | 07115877 |
| 34 | COVID 19 SWAB TESTING MACHINE | N/A | Dorchester Medium Institution | Infectious Disease Nurse's office | Abbott | Abbott ID Now | 14CDD81C |
| 35 | COVID-19 Swab Testing Machine | N/A | Dorchester Medium Institution | IDN Office | Abbot | Abbott ID NOW | 14CDD81C |
| 36 | Doppler | BMD-15928 | Dorchester Medium Institution | | Life Dop | | 2018060091 |
| 37 | Drill - Console - Surgical | BMD-0266 | Dorchester Medium Institution | Dental Office | Buffalo Manufacturing Comp. | X35Micromotor | N/A |
| 38 | Electrocardiograph | BMD-10382 | Dorchester Medium Institution | Treatment room 1 | GE Medical Systems | MAC5500 | 11-41583 |
| 39 | Frame Heater | N/A | Dorchester Medium Institution | Optometry Room | GFC | A600 | 682 |
| 40 | Hyfrecator | BMD-0262 | Dorchester Medium Institution | Treatment Room 1 | Birtcher | | 2157354 |
| 41 | Ice maker | N/A | Dorchester Medium Institution | B-6 End of hallway | Whirlpool | GI15NDXTB1 | EW1501128 |
| 42 | ID NOW Instrument Printer | N/A | Dorchester Medium Institution | Infectious Disease Nurse's office | ABBOTT | ID NOW | 2009A722604 |
| 43 | Imager - Computed - Radiography (Digitizer) | | Dorchester Medium Institution | X-ray room | AGFA | CR30-X | 2676 |
| 44 | Interferential - Ultrasound - Laser (Electro therapy unit) | BMD-11673 | Dorchester Medium Institution | Treatment room 2 | Chattanooga Canada | 2762CC Intellect Adv | 4978 |
| 45 | Lamp - Slit | BMD-11661 | Dorchester Medium Institution | Optometry Room | HAAG Streit | BQ900 | 10654 |

Annex F - Equipment List - Dorchester Institution (Medium)/ Annexe F - Liste d'équipements - Établissement Dorchester (Médium)

| | A | B | C | D | E | F | G |
|----|---|---------------------------------------|--------------------------------|----------------------------|------------------------------------|-------------------|------------------------------------|
| 1 | Item | ~~Equipment Code~~Code d'équipement~~ | ~~Institution~~Établissement~~ | ~~Location~~Localisation~~ | ~~Manufacturer~~Nom du fabricant~~ | ~~Model~~Modèle~~ | ~~Serial Number~~Numéro de série~~ |
| 46 | Lamp - Slit | BMD-1568 | Dorchester Medium Institution | Dr room | Top Con PS-20A | SL-2E | 6211049 |
| 47 | Lathe - Dental (Grinder / Polisher) | BMD-8061 | Dorchester Medium Institution | Dental Office | Red Wing | 16B | 243/10/11/044 |
| 48 | Lensmeter | BMD-11650 | Dorchester Medium Institution | X-ray room | Marco | 201 | M2474 |
| 49 | Lensometer | DOR-A026 | Dorchester Medium Institution | Optometry Room | | | |
| 50 | Light - Examining | BMD-10550 | Dorchester Medium Institution | Doctor's Clinic | Luxo | KFM-17213 | E24877 |
| 51 | Light - Examining | BMD-1542 | Dorchester Medium Institution | Treatment room 2 | Welch Allyn Canada Ltd | 44200 | -- |
| 52 | Light - Examining (Dental) | | Dorchester Medium Institution | Dental Office | Belmont | AL501T | T040461 |
| 53 | Light - Examining (UV inspection lamp) | BMD-11674 | Dorchester Medium Institution | Treatment room 1 | Luxo | 0316040 | 89705 |
| 54 | Light - Head | BMD-7410 | Dorchester Medium Institution | X-ray room | Heine | -- | 1130545 |
| 55 | Light - Head (Battery pack for head lamp) | BMD-7410 | Dorchester Medium Institution | X-ray room | Heine | M Pack | 2157354 |
| 56 | Light - Ultraviolet | BMD-10383 | Dorchester Medium Institution | Dental Office | Elipar | S10 | 939112017246 |
| 57 | Light - Ultraviolet (Curing) | | Dorchester Medium Institution | Dental Office | Coltène/Whaledent | Coltolux 75 | 060128666 |
| 58 | Monitor | N/A | Dorchester Medium Institution | IDN office | Samsung | LN26B460B2D | Z1HC3CXZ105915A |
| 59 | Monitor - Bloodpressure (with stand) | BMD-11178 | Dorchester Medium Institution | TBD | Welch Allyn | WA-71XX-B | 21100009864117 |
| 60 | Monitor - Bloodpressure (with stand) | BMD-11179 | Dorchester Medium Institution | TBD | Welch Allyn | WA-71XX-B | 21100009494117 |
| 61 | Monitor - bloodpressure (with stand) | BMD-11180 | Dorchester Medium Institution | TBD | Welch Allyn | WA-71XX-B | 100009694117 |
| 62 | Monitor - Bloodpressure (with stand) | BMD-11181 | Dorchester Medium Institution | TBD | Welch Allyn | WA-71XX-B | 10009924117 |
| 63 | Monitor - BP | BMD-15930 | Dorchester Medium Institution | Emergency Bag? | ? | BP-742CANN | 20190500113LG |
| 64 | Monitor - BP | | Dorchester Medium Institution | Emergency Bag? | ? | BP-742CANN | 20190500120LG |
| 65 | Monitor - Finger SPO2 | BMD-15894 | Dorchester Medium Institution | | Biolight | M70 | M010E523903 |
| 66 | Monitor - Finger SPO2 | BMD-15945 | Dorchester Medium Institution | | Biolight | M70 | M010E523929 |
| 67 | Monitor - Finger SPO2 | BMD-15946 | Dorchester Medium Institution | | Biolight | M70 | M010E523928 |

Annex F - Equipment List - Dorchester Institution (Medium)/ Annexe F - Liste d'équipements - Établissement Dorchester (Médium)

| | A | B | C | D | E | F | G |
|----|--------------------------------|---------------------------------------|--------------------------------|--------------------------------|------------------------------------|----------------------------|------------------------------------|
| 1 | Item | ~~Equipment Code~~Code d'équipement~~ | ~~Institution~~Établissement~~ | ~~Location~~Localisation~~ | ~~Manufacturer~~Nom du fabricant~~ | ~~Model~~Modèle~~ | ~~Serial Number~~Numéro de série~~ |
| 68 | Monitor - NIBP | BMD-10422 | Dorchester Medium Institution | Unknown - checking | | | |
| 69 | Monitor - NIBP | BMD-10551 | Dorchester Medium Institution | Nurses Clinic | Omron | HEM-907XL | 20140301019AF |
| 70 | Monitor - NIBP | BMD-15890 | Dorchester Medium Institution | B7 Clinic Room | Welsh Allyn | | 160928222325 |
| 71 | Monitor - NIBP - SPO2 - Temp | BMD-10381 | Dorchester Medium Institution | Treatment room 1 | GE Medical Systems | Dynamap DPC400N-CE | 2019205-001 |
| 72 | MONITOR - NIBP - SPO2 - TEMP | BMD-7963 | Dorchester Medium Institution | Health Care | Welch Allyn Canada Ltd | 42NTB | |
| 73 | Monitor - NIBP-SPO2-Temp | BMD-0251 | Dorchester Medium Institution | Storage room 3rd floor | Critikon | Dynamap Compact BP monitor | ref. 117207 |
| 74 | Monitor - Oximeter (hand held) | N/A | Dorchester Medium Institution | Chief's office | Smiths Medical ASD, Inc. | REF # WW1020EN | 4108067 |
| 75 | Monitor - Oximeter (hand held) | N/A | Dorchester Medium Institution | Chief's office | Smiths Medical ASD, INC. | REF # WW1020EN | 4108068 |
| 76 | Monitor - SPO2 | BMD-10413 | Dorchester Medium Institution | Emergency Bag - Front Entrance | Welch Allyn | WW1020EN | 4025971 |
| 77 | Monitor - SPO2 | BMD-10418 | Dorchester Medium Institution | Emergency Bag - HC | Masimo | | 526205 |
| 78 | Monitor - SPO2 | BMD-10423 | Dorchester Medium Institution | Unknown - checking | | | |
| 79 | Monitor - SPO2 | BMD-16004 | Dorchester Medium Institution | | Spectro 2 | Model 10 | 4088792 |
| 80 | Monitor - SPO2 | BMD-5215 | Dorchester Medium Institution | Treatment room 2 | Masimo Corporation | RAD5 | 532646 |
| 81 | Monitor - SPO2 | BMD-5220 | Dorchester Medium Institution | Treatment room 1 | Masimo Corporation | RAD5 | 530779 |
| 82 | Monitor - SPO2 | | Dorchester Medium Institution | | Spectro 2 | Model 10 | |
| 83 | Monitor - Temp | BMD-10412 | Dorchester Medium Institution | Emergency Bag - Front Entrance | Welch Allyn | 690 | 12480415 |
| 84 | Monitor - Temp | BMD-10419 | Dorchester Medium Institution | Clean Utility Room | Welsh Allyn | | |
| 85 | Monitor - Temp | BMD-10421 | Dorchester Medium Institution | Emergency Bag HC | Welch Allyn | 690 | 08011188 |
| 86 | Monitor - Temp | BMD-10544 | Dorchester Medium Institution | Clean Utility Room | Welsh Allyn | 690 | 14460118 |
| 87 | Monitor - Temp | BMD-11675 | Dorchester Medium Institution | Unknown - checking | Welch Allyn | | |
| 88 | Monitor - Temp | BMD-11676 | Dorchester Medium Institution | | Welch Allyn | 690 | 11110835 |
| 89 | Monitor - Temp | BMD-15892 | Dorchester Medium Institution | B7 - Clinic Room | Welch Allyn | 1690 | (21)16454159 |

Annex F - Equipment List - Dorchester Institution (Medium)/ Annexe F - Liste d'équipements - Établissement Dorchester (Médium)

| | A | B | C | D | E | F | G |
|-----|---|---------------------------------------|--------------------------------|----------------------------|------------------------------------|---------------------|------------------------------------|
| 1 | Item | ~~Equipment Code~~Code d'équipement~~ | ~~Institution~~Établissement~~ | ~~Location~~Localisation~~ | ~~Manufacturer~~Nom du fabricant~~ | ~~Model~~Modèle~~ | ~~Serial Number~~Numéro de série~~ |
| 90 | Monitor - Temp | BMD-16099 | Dorchester Medium Institution | | Welsh Allyn | SureTemp | (21)18491361 |
| 91 | Monitor - Temp | BMD-5219 | Dorchester Medium Institution | | Welch Allyn | 690 Thermometer | 08047847 |
| 92 | Monitor - Temp | BMD-8059 | Dorchester Medium Institution | Trauma Room | Welch Allyn | 690 | 12512584 |
| 93 | Monitor - Temp - Infrared | BMD-15932 | Dorchester Medium Institution | | Berrcom | JXB-178 | M22200475202 |
| 94 | Monitor - Temp - Infrared | BMD-15933 | Dorchester Medium Institution | | Berrcom | JXB-178 | M22200475212 |
| 95 | Monitor - Temp - Tympanic | BMD-15929 | Dorchester Medium Institution | | Welch Allyn | | (21) 18473111 |
| 96 | Monitor - Temp - Tympanic | BMD-15974 | Dorchester Medium Institution | IDN office | Covidien | Genius 3 | N18792640 |
| 97 | Monitor - vitals | BMD-15893 | Dorchester Medium Institution | TBD | Smiths Medical | SpectrO2 | 4081837 |
| 98 | Monitor - Vitals | BMD-15894 | Dorchester Medium Institution | TBD | Smiths Medical | Spectro2 - WW1000EN | 4082070 |
| 99 | Monitor - Vitals | BMD-15917 | Dorchester Medium Institution | ? | Welch Allyn | 42MTB | 202001103 |
| 100 | Monitor - Vitals | BMD-15945 | Dorchester Medium Institution | TBD | Smiths Medical | SpectrO2 - WW1000EN | 4082304 |
| 101 | Monitor - Vitals | BMD-15946 | Dorchester Medium Institution | TBD | Smiths Medical | SpectrO2 - WW1000EN | 4082303 |
| 102 | Ophthalmic - Instrument - stand | BMD-1482 | Dorchester Medium Institution | Optometry Room | | 7900 | 79010702014 |
| 103 | Ophthalmoscope | BMD-10424 | Dorchester Medium Institution | Treatment Room 1 | Welch Allyn | | |
| 104 | Ophthalmoscope | BMD-8075 | Dorchester Medium Institution | Treatment Room 2 | Welch Allyn | | |
| 105 | Ophthalmoscope (Integrated Diagnostic System) | BMD-8057 | Dorchester Medium Institution | Trauma Room | Welch Allyn | GS 777 WA 7791-MX | |
| 106 | Ophthalmoscope | BMD-15891 | Dorchester Medium Institution | B7 Rm 109A - Clinic Room | Welsh Allyn | G5777 | 408282 |
| 107 | Optometry equipment | | Dorchester Medium Institution | Optometry Room | NIDEK | | 130493 |
| 108 | Phoropter - Refractor | BMD-11662 | Dorchester Medium Institution | X-ray room | Marco | RT-I | 10975 |
| 109 | Projector - Ophthalmic | BMD-5223 | Dorchester Medium Institution | X-ray room | Marco | 50 | M22143 |
| 110 | Projector - Ophthalmic | BMD-5223 / DOR-A025 | Dorchester Medium Institution | X-ray room / Optometrist | Marco | -- | 22143 |
| 111 | Projector - Ophthalmic (Keratometer) | DOR-A025 | Dorchester Medium Institution | Optometry Room | Marco | | 21226 |

Annex F - Equipment List - Dorchester Institution (Medium)/ Annexe F - Liste d'équipements - Établissement Dorchester (Médium)

| | A | B | C | D | E | F | G |
|-----|-------------------------------------|---------------------------------------|--------------------------------|----------------------------|------------------------------------|-------------------------------|------------------------------------|
| 1 | Item | ~~Equipment Code~~Code d'équipement~~ | ~~Institution~~Établissement~~ | ~~Location~~Localisation~~ | ~~Manufacturer~~Nom du fabricant~~ | ~~Model~~Modèle~~ | ~~Serial Number~~Numéro de série~~ |
| 112 | Pump - Compressor | BMD-11655 | Dorchester Medium Institution | Trauma Room | AirLife | RT500 | 169542 |
| 113 | Pump - Compressor | BMD-11656 | Dorchester Medium Institution | | Carefusion (Pulmonetics) | RT500 | 169548 |
| 114 | Pump - Compressor | BMD-11672 | Dorchester Medium Institution | Storage 1st Floor | Schuco | 5711102 Schuco - Mist | 112586 |
| 115 | Pump - Compressor - O2 | BMD-15948 | Dorchester Medium Institution | 3rd Floor storage | NewLife | Ref: AS099-101, M766AS0991010 | BUB0116040440 |
| 116 | Pump - Compressor - O2 Concentrator | BMD-15949 | Dorchester Medium Institution | 3rd Floor Storage | | NewLife elite | 305603 |
| 117 | Pump - Compressor - O2 Concentrator | BMD-15950 | Dorchester Medium Institution | 3rd Floor Storage | Respironic | EverFlo Q OPI, Ref: 0120014 | 0091429 |
| 118 | Pump - Nebulizer | BMD-15952 | Dorchester Medium Institution | 3rd Floor Storage | DeVilbiss | PulmoAide, 5650C | C131112 |
| 119 | Pump - Vacume - Dental | BMD-10414 | Dorchester Medium Institution | Hall outside Dental Office | VacStar | VS20 | 208990 |
| 120 | PUMP-INFUSION | BMD-12018 | Dorchester Medium Institution | | Abbott | Plum 360 Large volume | 43457377 |
| 121 | PUMP-INFUSION | BMD-12019 | Dorchester Medium Institution | | Abbott | Plum 360 Large volume | 43468648 |
| 122 | Refrigerator - Biological | BMD-11657 | Dorchester Medium Institution | Unknown - checking | | | |
| 123 | Refrigerator - Biological | BMD-11881 | Dorchester Medium Institution | | Haier | None | NSA36LACG |
| 124 | Refrigerator - Biological | BMD-11883 | Dorchester Medium Institution | | LG Electronics Canada | GC-151RCW | Unknown |
| 125 | Refrigerator - Biological | BMD-8062 | Dorchester Medium Institution | Secondary Pharmacy | | SAN-MPR514 | 12010029 |
| 126 | Refrigerator - Biological | BMD-8063 | Dorchester Medium Institution | IDN Office | | SAN-SRL6111W | 120200257 |
| 127 | Retina Camera - Portable | | Dorchester Medium Institution | Optomitrist's Clinic | | VISUSCOUT 100 | 4519521151126 |
| 128 | Saw - Cast | BMD-1411 | Dorchester Medium Institution | Treatment room 2 | Stryker | 8127 | E34984 |
| 129 | Scale - Patient | BMD-16324 | Dorchester Medium Institution | Trauma Room | Seca | 703 1321993 | 5703160154691 |
| 130 | Scale - Patient | BMD-8076 | Dorchester Medium Institution | Treatment Room 1 | Pelstar | 500KL | 5000005579 |
| 131 | Screeener - Spot Vision | | Dorchester Medium Institution | Optomitrist Clinic | | VS100 | (21) 1000025994420 |
| 132 | SPHYGMOMANOMETER-ANALOG | BMD-12049 | Dorchester Medium Institution | | Welch Allyn Canada Ltd. | TYCOS | |
| 133 | Stand - optometrist equipment | BMD-5224 | Dorchester Medium Institution | X-ray room | Reliance Medical Products | 7900 | 79010702014 |

Annex F - Equipment List - Dorchester Institution (Medium)/ Annexe F - Liste d'équipements - Établissement Dorchester (Médium)

| | A | B | C | D | E | F | G |
|-----|--|---------------------------------------|--------------------------------|--------------------------------|------------------------------------|--------------------------|------------------------------------|
| 1 | Item | ~~Equipment Code~~Code d'équipement~~ | ~~Institution~~Établissement~~ | ~~Location~~Localisation~~ | ~~Manufacturer~~Nom du fabricant~~ | ~~Model~~Modèle~~ | ~~Serial Number~~Numéro de série~~ |
| 134 | Sterilizer - Steam - Table - top | BMD-11824 | Dorchester Medium Institution | Dental Office | SciCan | Statim 5000, G4-201102 | 511017I00018 |
| 135 | Sterilizer - Steam - Table - Top | BMD-8070 | Dorchester Medium Institution | Dental Office | SciScan | Statim 2000S - 01-121105 | 101212J00019 |
| 136 | Sterilizer - Steam - Table - top (Cassette autoclave) | BMD-11653 | Dorchester Medium Institution | Dental Office | SciCan | 1102 | 228756A |
| 137 | Stretcher | BMD-10540 | Dorchester Medium Institution | | Ferno Canada | 35X | 10N-149328 |
| 138 | Stretcher | BMD-11882 | Dorchester Medium Institution | Unknown - checking | Meditech | 5000 | TMPR1303029 |
| 139 | Stretcher (Stryker) | BMD-8068 | Dorchester Medium Institution | Trauma Room | Stryker | 0747-000 | 1105034705 |
| 140 | Stretcher (Stryker) | BMD-8073 | Dorchester Medium Institution | Doctor's Clinic | Stryker | 0747-000 | 1105034706 |
| 141 | Table - examining | BMD-15899 | Dorchester Medium Institution | B7 | | UMF5245 | 134817-0001 |
| 142 | Table - Physio treatment | BMD-22926 | Dorchester Medium Institution | Treatment Room 2 | Cardon | CTT | 1130-32321 |
| 143 | Table - XRAY (Proteus elevating table) | BMD-11666 | Dorchester Medium Institution | X-ray room | GE Hualun Medical System | 2259988 | 50508HL9 |
| 144 | Tonometer | BMD-11664 | Dorchester Medium Institution | X-ray room | Clement Clarke | Perkins tonometer MK2 | 1999T2073 |
| 145 | Tonometer | BMD-15900 | Dorchester Medium Institution | Optometry Room | icare | | 1839AJ024 |
| 146 | Traction unit - disposed | BMD-10384 | Dorchester Medium Institution | Storage Room | Chattanooga Canada | 2841/3336 | |
| 147 | Trolley | BMD-10382 | Dorchester Medium Institution | Trauma Room | GE Medical Systems | MAC 5000 | SCD11051932PA |
| 148 | Ultrasound - Diagnostic - Aorta - Bladder (Bladder scan) | BMD-11651 | Dorchester Medium Institution | Trauma room | Bladder Scan | 0570-0188 | PA011179 |
| 149 | Ultrasound - Therapeutic | | Dorchester Medium Institution | Treatment room 2 | Chattanooga Canada | 27335 | 37018 |
| 150 | Warmer - Blanket - Solution | BMD-7419 | Dorchester Medium Institution | Clean Utility Room | Pedigo | P-20105 | 545088-000 |
| 151 | Wheelchair | BMD-11818 | Dorchester Medium Institution | Equipment Storage Closet | Drive | Silver Sport II | 4S1512211210 |
| 152 | Wheelchair | BMD-11820 | Dorchester Medium Institution | Clean Equipment Storage Closet | Drive | | 4S1608142168 |
| 153 | Wheelchair | BMD-1187 | Dorchester Medium Institution | Equipment Storage Closet | Drive | Silver Sport II | 5D141100734 |
| 154 | Wheelchair | BMD-11962 | Dorchester Medium Institution | | | PCP Mobility | JL19121037 |
| 155 | Wheelchair | BMD-15939 | Dorchester Medium Institution | Clean Equipment Closet | Cardinal Health | | 8400001069-SH |

Annex F - Equipment List - Dorchester Institution (Medium)/ Annexe F - Liste d'équipements - Établissement Dorchester (Médium)

| | A | B | C | D | E | F | G |
|-----|---|--|-----------------------------------|-------------------------------|---------------------------------------|------------------------------|---------------------------------------|
| 1 | Item | ~~~Equipment Code~~~Code d'équipement~~~ | ~~~Institution~~~Établissement~~~ | ~~~Location~~~Localisation~~~ | ~~~Manufacturer~~~Nom du fabricant~~~ | ~~~Model~~~Modèle~~~ | ~~~Serial Number~~~Numéro de série~~~ |
| 156 | Wheelchair | BMD-15953 | Dorchester Medium Institution | 3rd Floor Storage | Tango Medical | Solara | 08LE000269 |
| 157 | Wheelchair | BMD-16006 | Dorchester Medium Institution | | Drive | Silver Sport II, SSP218FA-SF | 4S1805021446 |
| 158 | Wheelchair | BMD-16007 | Dorchester Medium Institution | | Drive | Silver Sport II, SSP218FA-SF | 4S1805021445 |
| 159 | Wheelchair | BMD-16008 | Dorchester Medium Institution | | Drive | Silver Sport II, SSP218FA-SF | 4S1805021442 |
| 160 | Wheelchair | | Dorchester Medium Institution | | Drive | Silver Sport II, SSP218FA-SF | 4S1805021437 |
| 161 | Wheelchair | | Dorchester Medium Institution | | Drive | Silver Sport II, SSP218FA-SF | 4S1805021441 |
| 162 | Wheelchair | | Dorchester Medium Institution | | Drive | Silver Sport II, SSP218FA-SF | 4S1805021443 |
| 163 | XRay - Bucky (Diagnostic X-ray wallstand) | BMD-11667 | Dorchester Medium Institution | X-ray room | GE Hualun Medical System | 600-0301 | 50997HL6 |
| 164 | XRay - RAD-14 tube | BMD-11699 | Dorchester Medium Institution | X-ray room | GE Hualun Medical System | 2259981 | 95296-7U |
| 165 | XRay - Unit - Suspension | BMD-11669 | Dorchester Medium Institution | Xray Room | GE Medsystems Ltd | 2269646 | 1014457WK4 |
| 166 | XRAY - Unit (Console) | | Dorchester Medium Institution | X-ray room | GE Hualun Medical System | 2259976 | 50529HL7 |
| 167 | XRAY - Unit (Proteus) | | Dorchester Medium Institution | X-ray room | GE Hualun Medical System | 2259973-2 | 50567HL7 + 5795 |
| 168 | X-ray Controller | | Dorchester Medium Institution | Dental Office | Belmont | 097-CM | EX10J0141 |
| 169 | X-ray Head | | Dorchester Medium Institution | Dental Office | Belmont | 097-H | EH10J0510 |
| 170 | X-ray printer | | Dorchester Medium Institution | X-ray room | Sharp AGFA | Drystar 5300 | |
| 171 | X-Ray processor - Dental | | Dorchester Medium Institution | Dental Office | Dell | Studio Hybrid 140g | 18HVBM1 |
| 172 | X-ray Regular Cone | | Dorchester Medium Institution | Dental Office | Belmont | 097-R | ER10J0510 |